AS SEEN IN

LAB MANAGEMENT TODAY

JANUARY · 2014

BUSINESS STRATEGIES FOR DENTAL LABORATORY DECISION MAKERS

in LAB TECHNOLOGY TODAY

Lab Technology Today

How to Fabricate Occlusally Screw-Retained Prettau® Bridges Using the Hybrid Technique

Zirkonzahn's Georg Walcher, MDT, offers his step-by-step technique for rehabilitating an edentulous patient with implants and fixed Zirconia Prettau Bridges.



AFTER: The happy patient with Prettau Zirconia bridges



clasps and wanted a more sented a major challenge.

The dental team—including Drs. Fernando Rojas-Vizcaya and Homavoun Zadeh, and Georg Walcher, MDT—decided to restore her smile with eight implants, four upper and four lower, and fixed, metal-free suprastructures. Due to the limited bone available, the tear chose occlusally screw-retained Prettau Bridges—fully ana-tomical restorations made with ranslucent Prettau Zirconia tha are ideal for implant-supported restorations, limited-space with a gingival flange.



at sites #4, #6, #11 and #13. #23, #26 and #28 and, since there were only four implant and an extensive distal eder tulous space, the team de-cided to use a titanium bar to quarantee bridge stability. that the esthetics and handling benefits of zirconia oupled with titanium's high flexural strength would resul in strong, esthetic and bio



the articulator and scanned using the fully automated, optical structured light scan ner S600 ARTI. The software ises the articulator scan to and align the casts in perfect occlusion, allowing the identification and elimination of static and dynamic




the screw-access channels The software then analyzed and adapted it to the gingival

THE DENTAL TEAM



Georg Walcher, MDT, has been a zirconia Fernando Rojas-Vizcaya, DDS, MS, Adjunct Assistant and CAD/CAM expert at Dental Laboratory Steger in Bruneck, Italy since 2007. His education includes diplomas in Dental Technology and Dental Technical Assistant from Homayoun H. Zadeh, DDS, Ph.D, Associate the Istituto Professionale per l'Industria e L'Artigianato in Bolzano, Italy. He earned his Ostrow School of Dentistry, Los Angeles, CA, and master's certificate in Baden, Vienna in 2000.



Next, to determine the optimal position of the titanium support bar and define the emergence profile, the mandibular digital waxup was displayed on screen (left photo). The corresponding implant types in the software library were used to plan the bar construction that would be directly screw retained on the implants (right photo). The software stores all files as open .stl files, which can be easily exported. The bar and connectors were then milled in titanium using the M5, a 5+1-axis simultaneous milling unit, and the bar was polished. Before designing the mandibular bridge, the fit of the milled titanium bar was examined on the master cast.



After some minor adjustments, the restorations were milled in Prettau Zirconia with the M5 milling unit using the simultaneous 5-axis processing in combination with a 0.3 C milling bur that milled the extra-fine details.



In preparation for the anterior porcelain veneer, the zirconia bridge was manually cut back and characterized. The incisal edge of the maxillary bridge retained its full anatomic contour in the functional loading area, protecting the edge and preventing chipping.



The bar was scanned and the mandibular bridge designed on screen. To ensure adequate fit, function, esthetics and phonetics, temporary bridges were milled from Temp Basic resin (left photo) and veneered with gingival composites (right photo). The titanium bar and temporaries were sent to the dentist for the try-in.



For a natural color, the restorations were brush-painted with acid-free Colour Liquids Prettau Aquarell; each stroke was manually applied with a thin brush and then the restorations were placed under a drying lamp. For visual aid, the liquids contain special biopigments that help in realistic shade distribution and color grading. Individual custom coloring of finer details such as mamelons, cervical or interproximal areas is considerably easier and can be performed with pinpoint precision.



After sintering in the Zirkonofen 700 oven, the Prettau Zirconia bridges exhibited natural esthetics and high translucency.



The fit of the titanium bar was tested in the maxillary bridge. Once the fit was confirmed. the gingival area and the anteriors were veneered with ICE Zirconia ceramic. In the anterior region, only vestibular ICE Zirconia dentines, effects and incisals were used; the incisal edge remained unveneered to avoid chipping. After stain and glaze firings, the restoration was complete.



Since the titanium bar gave the restoration a high value and grav hue, the bar and attached titanium bases were anodized with Titanium Spectral Colouring Anodizer to produce a biocompatible, gold-colored oxide layer with a warm tone.



After steam cleaning, the titanium bases in the maxillary bridge (shown here) and the titanium bar in the mandibular bridge were bonded to the zirconia structures with selfadhesive luting composite.



The finished Prettau Bridges were tried on the master cast and sent to the dentist for delivery.

LMT • January 2014 • www.LMTmag.com 103

AS SEEN IN

LAB MANAGEMENT TODAY

JANUARY · 2014

BUSINESS STRATEGIES FOR DENTAL LABORATORY DECISION MAKERS

in LAB DAY PRODUCT CATALOG on page 24

MILLING UNIT COMPACT LINE M1

Exhibitor: Zirkonzahn USA, Inc.

ZIRKONZAHN's newly developed Milling Unit Compact Line M1 makes CAD/CAM technology accessible to every laboratory. At the top of the line is the five-axis Milling Unit M1 Wet Heavy that can process zirconia, titanium, chrome/cobalt, PMMA, wax, prefabricated titanium abutments as well as glass ceramics. The intelligent five-axis orbit technology and a milling spindle with special high torque offer extra stability, a necessity for milling hard metals.

800-989-8931, 678-441-9419 infousa@zirkonzahn.com www.zirkonzahn.com

