



PRETTAU® 3 DISPERSIVE®

The Zirkonzahn Culture



PRETTAU® 3 DISPERSIVE® ZIRCONIA STRUCTURE ON ANODISED TITANIUM BAR

Available data:

Photos, 3D facial scan data, digital oral situation, digitised master models

Planned restoration:

Full arch bridge for the mandible made of Prettau® 3 Dispersive® zirconia on anodised titanium bar

Realisation:

- Digital mounting of the oral situation with Face Hunter 3D facial scanner and PlaneSystem®
- Initial virtual tooth set-up with Zirkonzahn. Modifier, individualisation of tooth shapes with selection from Heroes Collection virtual tooth library
- Titanium bar design in the Zirkonzahn. Modellier software using the planned restoration as a situ scan; milling and digitisation of the titanium bar, adaptation of the planned lower jaw zirconia structure
- Positioning of the structure in Zirkonzahn.Nesting: the colour gradient visualisation of the Prettau® 3 Dispersive® Gradual-Triplex-Technology provides an optimal alignment of the incisal aspect in the highly translucent blank area and the cervical area in the high strength region
- Milling of the Prettau® Bridge in the M2 Dual Teleskoper milling unit and sintering at 1500°C
- Glazing with 3D Base Glaze and characterisation with ICE Stains 3D by Enrico Steger
- Minimal veneering of the gingiva (0.4–0.5 mm) with Fresco Gingiva ceramic pastes; firing at 790 °C in the ceramic furnace; slight polishing
- Anodisation and cementation of the titanium bar in the zirconia structure; final colour check with the Zirkonzahn Individual Shade Guide Tooth 11









100% MONOLITHIC DESIGN, VENEERED ONLY IN THE GINGIVAL AREA

DT Alexander Lichtmannegger - Zirkonzahn Education Center Brunico, South Tyrol, Italy














NEW! PRETTAU® 3 DISPERSIVE® WITH **GRADUAL-TRIPLEX-TECHNOLOGY**

- Already during the manufacture of the material are worked into colour, translucency and flexural strength:
 - 1. Cervically increasing flexural strength; extremely high flexural strength at the neck of the tooth
 - 2. Incisally increasing translucency; highly translucent incisal edge
 - 3. Natural colour gradient from dentine to enamel
- Indicated for reduced or monolithic single crowns, inlays, onlays, veneers and bar-supported multi-unit bridges; especially suitable for monolithic design
- No ceramic chipping (thanks to monolithic design); fast sintering of single crowns possible
- Can be characterised individually for each patient with Colour Liquid Prettau®Aquarell Intensive, ICE Ceramics, Fresco Ceramics and ICE 3D Stains by Enrico Steger



HUMAN ZIRCONIUM TECHNOLOGY

Zirkonzahn Worldwide - Tel +39 0474 066 680 - info@zirkonzahn.com - www.zirkonzahn.com

NEW! 1 ORBIT – 4 DIAMETERS 125, 106, 98 AND 95 mm

M2 MILLING UNIT COMFORT LINE, WITH EXTRA LARGE TELESKOPER ORBIT. FULLY AUTOMATIC, FLEXIBLE, VIBRATION-FREE







