



PRETTAU[®] 2 DISPERSIVE[®]

The Zirkonzahn Culture

Zirkonzahn

PRETTAU® 2 DISPERSIVE® ZIRCONIA RESTORATION ON ANODISED TITANIUM BAR

Initial situation:
Patient request:
Patient data:
Planned restoration:
Realisation

Elderly patient with maxillary edentulism; four implants
A zirconia bridge to replace the old restoration (aesthetically and functionally inadequate)
Photos, master models with implants, situ scans of the old restoration
Monolithic Prettau[®] 2 Dispersive[®] zirconia bridge on a titanium bar, veneered only in the gingiva
Initial tooth set-up in Zirkonzahn.Modifier, based on the digitised models and situ scans

- of the presented restoration; individual adaptation of the tooth forms from the Heroes Collection rooted tooth library
- Titanium bar design in Zirkonzahn.Modellier using the wax-up of the planned restoration; milling and subsequent digitisation of the bar
- Digital adjustment of the planned maxillary zirconia structure and verification in the PS1 virtual articulator
- Fabrication of the structure with the Zirkonzahn milling unit
- Minimal colour accentuation and sintering
- Manual individualisation for age-appropriate aesthetics: colouring with ICE Stains 3D by Enrico Steger and layering of the gingival aspect with ICE Ceramics
- Check and anodisation of the bar, cementation

Result:

Prettau[®] bridge with optimal aesthetics and function thanks to the strong cooperation within the treatment team. After placement in the mouth, the patient remarks with excitement: "This bridge is my best Christmas gift."









100 % MONOLITHIC DESIGN, VENEERED ONLY IN THE GINGIVAL AREA

DT Federico Presicci, Dental Laboratory Antonio Lazetera, Savona, Italy Dr Massimo Callerio, Cassolnovo, Italy



Zirkonzahn

NEW! PRETTAU[®] 2 DISPERSIVE[®]

- Pre-coloured zirconia with a natural colour gradient, optimal flexural strength and particularly high translucency
- No limitations! Ideal for monolithic full arch bridges, but also for single crowns, inlays, onlays, veneers, bars and multi-unit bridges (reduced or monolithic)
- No ceramic chipping (thanks to monolithic design), no abrasion of the antagonist
- Can be characterised individually for each patient with Colour Liquid Prettau[®] 2 Aquarell intensive colours, ICE Ceramics and ICE 3D Stains by Enrico Steger



HUMAN ZIRCONIUM TECHNOLOGY

Zirkonzahn Worldwide – Tel +39 0474 066 680 – info@zirkonzahn.com – www.zirkonzahn.com













NEW! M2 MILLING UNIT COMFORT LINE

Now with the Teleskoper Orbit for material blanks with Ø 95, 98, 106 and 125 mm as well as for glass-ceramic or titanium abutments; high-precision repositioning in µm range