



### PRETTAU® 2 DISPERSIVE®

The Zirkonzahn Culture



# FULL ARCH IN PRETTAU® 2 DISPERSIVE® OVER TITANIUM BAR WITH EXTERNAL STOP

Initial situation: Male patient, middle-aged, maxillary edentulism with four implants, full denture temporary

restoration, minimal space due to occlusal and anterior implant position, vertical dimension

should be maintained

Patient data: Photos, master implant models, situ scans

Planned restoration: Bridge in Prettau® 2 Dispersive® over a titanium bar with external stop; zirconia structure

can thus be designed more delicately

Realisation: - Transfer patient data into virtual articulator

- Titanium bar design in Zirkonzahn. Modellier based on situ scans of full denture; titanium bar (11–21) individualised using the Free-Form function

- Export bar design as a wax-up into Zirkonzahn. Modifier. Complete virtual tooth setup with natural tooth libraries

- Adapt virtual Prettau® Bridge to the bar in Zirkonzahn.Modellier; verification in the PSI virtual articulator


- Mill titanium bar and zirconia bridge in pre-coloured Prettau® 2 Dispersive® in one milling process with the M2 Dual Wet Heavy Metal milling unit

- Minimal colour accentuation and sintering of the zirconia structure


- Finalise monolithic bridge with ICE Stains 3D

- Layer gingival aspect with ICE Ceramics

- Check result on milled implant model with ScanAnalogs; bond bar into zirconia structure









## 100% MONOLITHIC DESIGN, VENEERED ONLY IN THE GINGIVAL AREA

Zirkonzahn Education Center Brunico, South Tyrol, Italy











#### NEW! PRETTAU® 2 DISPERSIVE®

- Pre-coloured zirconia with a natural colour gradient, optimal flexural strength and particularly high translucency
- No limitations! Ideal for monolithic full arch bridges, but also for single crowns, inlays, onlays, veneers, bars and multiunit bridges (reduced or monolithic)
- No ceramic chipping (thanks to monolithic design), no abrasion of the antagonist
- Can be characterised individually for each patient with Colour Liquid Prettau® 2 Aquarell intensive colours, ICE Ceramics and ICE Stains 3D



#### **HUMAN ZIRCONIUM TECHNOLOGY**

Zirkonzahn Worldwide - Tel +39 0474 066 680 - info@zirkonzahn.com - www.zirkonzahn.com

#### NEW! M2 DUAL WET HEAVY METAL MILLING UNIT

ONE MACHINE, TWO MILLING CHAMBERS

FOR SEQUENTIAL WET AND DRY PROCESSING WITHOUT IN-BETWEEN CLEANING



WET DRY

