



PRETTAU[®] 2

The Zirkonzahn Culture

Zirkonzahn

COMPLETE PRETTAU® 2 ZIRCONIA RESTORATION ON ANODISED TITANIUM BARS


Available data:

Photos, 3D facial scans, digital situation of the mouth, X-ray data, digital master models, cross-mounted scans of immediate provisional

Planned restoration: Realisation:

- Full-arch maxillary and mandibular bridges made of Prettau® 2 on anodised titanium bars
- Digital mounting of the oral situation with Face Hunter 3D facial scanner and PlaneSystem[®]
- Design of the first virtual set-up with Zirkonzahn. Modifier
- Creation of prototypes in Multistratum[®] Flexible resin on titanium bases
- After an extended wear period: scans of the prototypes for the final restoration
- Design of titanium bars in Zirkonzahn. Modellier on the basis of the wax-ups: the lingual and palatal sides of the bars were partially adopted from the prototypes design in order to incorporate the original shape into the final restoration
- Adjusting upper and lower wax-ups on titanium bars, positioning the structures in Zirkonzahn.Nesting and milling in the M2 Dual Teleskoper milling unit
 - Colour infiltration with Colour Liquids Prettau[®] 2 Aquarell and sintering at 1600 °C
- Temporary cementation of Prettau[®] 2 structures on titanium bars for further verification of aesthetic and functional parameters
- Minimal veneering of the gingival area with Fresco Gingiva and characterisation with ICE Stains 3D by Enrico Steger
- Anodisation and bonding of the titanium bars in the zirconia structures










100% FUNCTION AND AESTHETICS A SYMBIOSIS OF TITANIUM AND ZIRCONIA

MDT Antonio Corradini – Zirkonzahn Education Center Brunico, South Tyrol, Italy Dr. med. dent. Manrique Fonseca MDT Alexander Bassermann









Zirkonzahn

MORE INFORMATION



NEW! PRETTAU® 2

- Particularly highly translucent zirconia with an excellent flexural strength
- No limitations! Especially suitable for full arch restorations (full anatomical or reduced for ceramic veneering), but also for single crowns, inlays, onlays, veneers, bars and multi-unit bridges
- No ceramic chipping (thanks to monolithic design), no abrasion of the antagonist
- Can be characterised individually for each patient with Colour Liquids Prettau[®] 2 Aquarell, ICE Ceramics, Fresco Ceramics and ICE Stains 3D by Enrico Steger
- Available also pre-coloured with a natural colour gradient (Prettau[®] 2 Dispersive[®])



HUMAN ZIRCONIUM TECHNOLOGY

Zirkonzahn Worldwide – Tel +39 0474 066 680 – info@zirkonzahn.com – www.zirkonzahn.com

NEW! 1 ORBIT – 4 DIAMETERS 125, 106, 98 AND 95 mm

M2 MILLING UNIT COMFORT LINE, WITH EXTRA LARGE TELESKOPER ORBIT. FULLY AUTOMATIC, FLEXIBLE, VIBRATION-FREE

LARGER MILLING AREA WITH Ø 95 mm

2 x bite splints







3 x bite splints

