

ALPHA-BIO TEC® SPI/DFI/ATID

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE AB-SP-ABCH 3.75 (Raw-Abutment® customizable hex and screw)	BSAF1881	20	-	10-12
SET ZZ BASE AB-SP-CCTB 3.75 (Conical cemented titan base non hex and screw)	BSAD1831	20	-	10-12
SET ZZ BASE AB-SP-CCTB TC (Conical cemented titan base non hex and screw)	BSAD1832	25	-	10-12
SET ZZ BASE AB-SP-PCTBH 3.75 (Parallel cemented titan base hex and screw)	BSAB1881	20	-	10-12
SET ZZ BASE AB-SP-SM 3.75 (Scanmarker and screw)	ZBAD1831	-	-	10-12
SET ZZ BASE AB-SP-SM TC (Scanmarker and screw)	ZBAD1832	-	-	10-12
SET ZZ BASE AB-SP-WSM 3.75 (White Scanmarker and screw)	ZBCD1831	-	-	10-12
SET ZZ BASE AB-SP-WSM TC (White Scanmarker and screw)	ZBCD1832	-	-	10-12
ZZ BASE AB-SP-MI 3.75 (Model implant)	MIAG1801	-	-	10-12
ZZ BASE AB-SP-MI TC (Model implant)	MIAG1802	-	-	10-12
ZZ BASE AB-SP-SC 3.75 (Screw)	BSAE1801	20	-	10-12
ZZ BASE AB-SP-SC TC (Screw)	BSAE1802	25	-	10-12
ZZ BASE AB-SP-SC-G 3.75 (Screw Gold)	BSCE1801	20	-	10-12
ZZ BASE AB-SP-ZSC 3.75 (Zirconia screw)	BSAE1831	20	-	10-12
ZZ BASE AB-SP-ZSC TC (Zirconia screw)	BSAE1832	25	-	10-12

ANTHOGYR AXIOM®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE AG-AX-ABCH REG/PX (Raw-Abutment® customizable hex and screw)	BSAF8381	25	-	10-12
SET ZZ BASE AG-AX-CCTB ACA (Conical cemented titan base non hex and screw)	BSAD8333	25	-	10-12
SET ZZ BASE AG-AX-CCTB MUA 4.0 (Conical cemented titan base non hex and screw)	BSAD8334	15	-	10-12
SET ZZ BASE AG-AX-CCTB MUA 4.8 (Conical cemented titan base non hex and screw)	BSAD8335	15	-	10-12
SET ZZ BASE AG-AX-CCTB REG/PX (Conical cemented titan base non hex and screw)	BSAD8331	25	-	10-12
SET ZZ BASE AG-AX-CCTB SCA (Conical cemented titan base non hex and screw)	BSAD8332	25	-	10-12
SET ZZ BASE AG-AX-NTB REG/PX (Narrow titan base non hex and screw)	BSAP8331	25	-	10-12
SET ZZ BASE AG-AX-NTBH REG/PX (Narrow titan base hex and screw)	BSAP8381	25	-	10-12
SET ZZ BASE AG-AX-PCTBH REG/PX (Parallel cemented titan base hex and screw)	BSAB8381	25	-	10-12
SET ZZ BASE AG-AX-SM ACA (Scanmarker and screw)	ZBAD8333	-	-	10-12
SET ZZ BASE AG-AX-SM MUA 4.0 (Scanmarker and screw)	ZBAD8334	-	-	10-12
SET ZZ BASE AG-AX-SM MUA 4.8 (Scanmarker and screw)	ZBAD8335	-	-	10-12
SET ZZ BASE AG-AX-SM REG/PX (Scanmarker and screw)	ZBAD8331	-	-	10-12
SET ZZ BASE AG-AX-SM SCA (Scanmarker and screw)	ZBAD8332	-	-	10-12
SET ZZ BASE AG-AX-WSM ACA (White Scanmarker and screw)	ZBCD8333	-	-	10-12
SET ZZ BASE AG-AX-WSM MUA 4.0 (White Scanmarker and screw)	ZBCD8334	-	-	10-12
SET ZZ BASE AG-AX-WSM MUA 4.8 (White Scanmarker and screw)	ZBCD8335	-	-	10-12
SET ZZ BASE AG-AX-WSM REG/PX (White Scanmarker and screw)	ZBCD8331	-	-	10-12
SET ZZ BASE AG-AX-WSM SCA (White Scanmarker and screw)	ZBCD8332	-	-	10-12
ZZ BASE AG-AX-SC ACA (Screw)	BSAE8303	25	-	10-12
ZZ BASE AG-AX-SC MUA 4.0 (Screw)	BSAE8304	15	-	10-12
ZZ BASE AG-AX-SC MUA 4.8 (Screw)	BSAE8305	15	-	10-12
ZZ BASE AG-AX-SC REG/PX (Screw)	BSAE8301	25	-	10-12
ZZ BASE AG-AX-SC SCA (Screw)	BSAE8302	25	-	10-12
ZZ BASE AG-AX-ZSC ACA (Zirconia Screw)	BSAE8333	25	-	10-12
ZZ BASE AG-AX-ZSC MUA 4.0 (Zirconia Screw)	BSAE8334	15	-	10-12
ZZ BASE AG-AX-ZSC MUA 4.8 (Zirconia Screw)	BSAE8335	15	-	10-12
ZZ BASE AG-AX-ZSC REG/PX (Zirconia Screw)	BSAE8331	25	-	10-12
ZZ BASE AG-AX-ZSC SCA (Zirconia Screw)	BSAE8332	25	-	10-12

ASTRA TECH IMPLANT SYSTEM™ EV

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE AS-EV-ABCH 3.0 (Raw-Abutment® customizable hex and screw)	BSAF7681	25	-	10-12
SET ZZ BASE AS-EV-ABCH 3.6 (Raw-Abutment® customizable hex and screw)	BSAF7682	25	-	10-12
SET ZZ BASE AS-EV-ABCH 4.2 (Raw-Abutment® customizable hex and screw)	BSAF7683	25	-	10-12
SET ZZ BASE AS-EV-ABCH 4.8 (Raw-Abutment® customizable hex and screw)	BSAF7684	25	-	10-12
SET ZZ BASE AS-EV-ABCH 5.4 (Raw-Abutment® customizable hex and screw)	BSAF7685	25	-	10-12
SET ZZ BASE AS-EV-ABCH-D14 3.0 (Raw-Abutment® customizable hex D14 and screw)	BSAF7681D14	25	-	10-12
SET ZZ BASE AS-EV-ABCH-D14 3.6 (Raw-Abutment® customizable hex D14 and screw)	BSAF7682D14	25	-	10-12
SET ZZ BASE AS-EV-ABCH-D14 4.2 (Raw-Abutment® customizable hex D14 and screw)	BSAF7683D14	25	-	10-12
SET ZZ BASE AS-EV-ABCH-D14 4.8 (Raw-Abutment® customizable hex D14 and screw)	BSAF7684D14	25	-	10-12
SET ZZ BASE AS-EV-ABCH-D14 5.4 (Raw-Abutment® customizable hex D14 and screw)	BSAF7685D14	25	-	10-12
SET ZZ BASE AS-EV-CCTB 3.0 (Conical cemented titan base non hex and screw)	BSAD7631	25	-	10-12
SET ZZ BASE AS-EV-CCTB 3.6 (Conical cemented titan base non hex and screw)	BSAD7632	25	-	10-12
SET ZZ BASE AS-EV-CCTB 4.2 (Conical cemented titan base non hex and screw)	BSAD7633	25	-	10-12
SET ZZ BASE AS-EV-CCTB 4.8 (Conical cemented titan base non hex and screw)	BSAD7634	25	-	10-12
SET ZZ BASE AS-EV-CCTB 5.4 (Conical cemented titan base non hex and screw)	BSAD7635	25	-	10-12
SET ZZ BASE AS-EV-CCTB P4.2 (Conical cemented titan base non hex and screw)	BSAD7637	25	-	10-12
SET ZZ BASE AS-EV-CCTB P4.8 (Conical cemented titan base non hex and screw)	BSAD7638	25	-	10-12
SET ZZ BASE AS-EV-CCTB UA (Conical cemented titan base non hex and screw)	BSAD7636	15	-	10-12
SET ZZ BASE AS-EV-CCTB-L15 3.6 (Conical cemented titan base non hex L15 and screw)	BSAD7632L15	25	-	10-12
SET ZZ BASE AS-EV-CCTB-L15 4.2 (Conical cemented titan base non hex L15 and screw)	BSAD7633L15	25	-	10-12
SET ZZ BASE AS-EV-CCTB-L15 4.8 (Conical cemented titan base non hex L15 and screw)	BSAD7634L15	25	-	10-12
SET ZZ BASE AS-EV-CCTB-L15 5.4 (Conical cemented titan base non hex L15 and screw)	BSAD7635L15	25	-	10-12
SET ZZ BASE AS-EV-CCTB-L20 3.0 (Conical cemented titan base non hex L20 and screw)	BSAD7631L20	25	-	10-12
SET ZZ BASE AS-EV-CCTB-L20 3.6 (Conical cemented titan base non hex L20 and screw)	BSAD7632L20	25	-	10-12
SET ZZ BASE AS-EV-CCTB-L20 4.8 (Conical cemented titan base non hex L20 and screw)	BSAD7634L20	25	-	10-12
SET ZZ BASE AS-EV-CCTB-L20 5.4 (Conical cemented titan base non hex L20 and screw)	BSAD7635L20	25	-	10-12
SET ZZ BASE AS-EV-CCTB-L30 3.0 (Conical cemented titan base non hex L30 and screw)	BSAD7631L30	25	-	10-12
SET ZZ BASE AS-EV-CCTB-L30 3.6 (Conical cemented titan base non hex L30 and screw)	BSAD7632L30	25	-	10-12
SET ZZ BASE AS-EV-CCTB-L30 4.2 (Conical cemented titan base non hex L30 and screw)	BSAD7633L30	25	-	10-12
SET ZZ BASE AS-EV-CCTB-L30 4.8 (Conical cemented titan base non hex L30 and screw)	BSAD7634L30	25	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

ASTRA TECH IMPLANT SYSTEM™ EV

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE AS-EV-CCTB-L30 5.4 (Conical cemented titan base non hex L30 and screw)	BSAD7635L30	25	-	10-12
SET ZZ BASE AS-EV-CCTB-L40 3.0 (Conical cemented titan base non hex L40 and screw)	BSAD7631L40	25	-	10-12
SET ZZ BASE AS-EV-CCTB-L40 3.6 (Conical cemented titan base non hex L40 and screw)	BSAD7632L40	25	-	10-12
SET ZZ BASE AS-EV-CCTB-L40 4.2 (Conical cemented titan base non hex L40 and screw)	BSAD7633L40	25	-	10-12
SET ZZ BASE AS-EV-CCTB-L40 4.8 (Conical cemented titan base non hex L40 and screw)	BSAD7634L40	25	-	10-12
SET ZZ BASE AS-EV-CCTB-L40 5.4 (Conical cemented titan base non hex L40 and screw)	BSAD7635L40	25	-	10-12
SET ZZ BASE AS-EV-NTB 3.0 (Narrow titan base non hex and screw)	BSAP7631	25	-	10-12
SET ZZ BASE AS-EV-NTB 3.6 (Narrow titan base non hex and screw)	BSAP7632	25	-	10-12
SET ZZ BASE AS-EV-NTB 4.2 (Narrow titan base non hex and screw)	BSAP7633	25	-	10-12
SET ZZ BASE AS-EV-NTB 4.8 (Narrow titan base non hex and screw)	BSAP7634	25	-	10-12
SET ZZ BASE AS-EV-NTBH 3.0 (Narrow titan base hex and screw)	BSAP7681	25	-	10-12
SET ZZ BASE AS-EV-NTBH 3.6 (Narrow titan base hex and screw)	BSAP7682	25	-	10-12
SET ZZ BASE AS-EV-NTBH 4.2 (Narrow titan base hex and screw)	BSAP7683	25	-	10-12
SET ZZ BASE AS-EV-NTBH 4.8 (Narrow titan base hex and screw)	BSAP7684	25	-	10-12
SET ZZ BASE AS-EV-PCTBH 3.0 (Parallel cemented titan base hex and screw)	BSAB7681	25	-	10-12
SET ZZ BASE AS-EV-PCTBH 3.6 (Parallel cemented titan base hex and screw)	BSAB7682	25	-	10-12
SET ZZ BASE AS-EV-PCTBH 4.2 (Parallel cemented titan base hex and screw)	BSAB7683	25	-	10-12
SET ZZ BASE AS-EV-PCTBH 4.8 (Parallel cemented titan base hex and screw)	BSAB7684	25	-	10-12
SET ZZ BASE AS-EV-PCTBH 5.4 (Parallel cemented titan base hex and screw)	BSAB7685	25	-	10-12
SET ZZ BASE AS-EV-PCTBH-L15 3.0 (Parallel cemented titan base hex L15 and screw)	BSAB7691L15	25	-	10-12
SET ZZ BASE AS-EV-PCTBH-L15 3.6 (Parallel cemented titan base hex L15 and screw)	BSAB7692L15	25	-	10-12
SET ZZ BASE AS-EV-PCTBH-L15 4.2 (Parallel cemented titan base hex L15 and screw)	BSAB7693L15	25	-	10-12
SET ZZ BASE AS-EV-PCTBH-L15 4.8 (Parallel cemented titan base hex L15 and screw)	BSAB7694L15	25	-	10-12
SET ZZ BASE AS-EV-PCTBH-L15 5.4 (Parallel cemented titan base hex L15 and screw)	BSAB7695L15	25	-	10-12
SET ZZ BASE AS-EV-PCTBH-L20 3.0 (Parallel cemented titan base hex L20 and screw)	BSAB7691L20	25	-	10-12
SET ZZ BASE AS-EV-PCTBH-L20 3.6 (Parallel cemented titan base hex L20 and screw)	BSAB7692L20	25	-	10-12
SET ZZ BASE AS-EV-PCTBH-L20 4.2 (Parallel cemented titan base hex L20 and screw)	BSAB7693L20	25	-	10-12
SET ZZ BASE AS-EV-PCTBH-L20 4.8 (Parallel cemented titan base hex L20 and screw)	BSAB7694L20	25	-	10-12
SET ZZ BASE AS-EV-PCTBH-L20 5.4 (Parallel cemented titan base hex L20 and screw)	BSAB7695L20	25	-	10-12
SET ZZ BASE AS-EV-PCTBH-L30 3.0 (Parallel cemented titan base hex L30 and screw)	BSAB7691L30	25	-	10-12

ASTRA TECH IMPLANT SYSTEM™ EV

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE AS-EV-PCTBH-L30 3.6 (Parallel cemented titan base hex L30 and screw)	BSAB7692L30	25	-	10-12
SET ZZ BASE AS-EV-PCTBH-L30 4.2 (Parallel cemented titan base hex L30 and screw)	BSAB7693L30	25	-	10-12
SET ZZ BASE AS-EV-PCTBH-L30 4.8 (Parallel cemented titan base hex L30 and screw)	BSAB7694L30	25	-	10-12
SET ZZ BASE AS-EV-PCTBH-L30 5.4 (Parallel cemented titan base hex L30 and screw)	BSAB7695L30	25	-	10-12
SET ZZ BASE AS-EV-PCTBH-L40 3.0 (Parallel cemented titan base hex L40 and screw)	BSAB7691L40	25	-	10-12
SET ZZ BASE AS-EV-PCTBH-L40 3.6 (Parallel cemented titan base hex L40 and screw)	BSAB7692L40	25	-	10-12
SET ZZ BASE AS-EV-PCTBH-L40 4.2 (Parallel cemented titan base hex L40 and screw)	BSAB7693L40	25	-	10-12
SET ZZ BASE AS-EV-PCTBH-L40 4.8 (Parallel cemented titan base hex L40 and screw)	BSAB7694L40	25	-	10-12
SET ZZ BASE AS-EV-PCTBH-L40 5.4 (Parallel cemented titan base hex L40 and screw)	BSAB7695L40	25	-	10-12
SET ZZ BASE AS-EV-SA 3.0 (Scananalog and screw)	ZBDE7631	15	-	10-12
SET ZZ BASE AS-EV-SA 3.6 (Scananalog and screw)	ZBDE7632	15	-	10-12
SET ZZ BASE AS-EV-SA 4.2 (Scananalog and screw)	ZBDE7633	15	-	10-12
SET ZZ BASE AS-EV-SA 4.8 (Scananalog and screw)	ZBDE7634	15	-	10-12
SET ZZ BASE AS-EV-SA 5.4 (Scananalog and screw)	ZBDE7635	15	-	10-12
SET ZZ BASE AS-EV-SM 3.0 (Scanmarker and screw)	ZBAD7631	-	-	10-12
SET ZZ BASE AS-EV-SM 3.6 (Scanmarker and screw)	ZBAD7632	-	-	10-12
SET ZZ BASE AS-EV-SM 4.2 (Scanmarker and screw)	ZBAD7633	-	-	10-12
SET ZZ BASE AS-EV-SM 4.8 (Scanmarker and screw)	ZBAD7634	-	-	10-12
SET ZZ BASE AS-EV-SM 5.4 (Scanmarker and screw)	ZBAD7635	-	-	10-12
SET ZZ BASE AS-EV-SM P4.2 (Scanmarker and screw)	ZBAD7637	-	-	10-12
SET ZZ BASE AS-EV-SM P4.8 (Scanmarker and screw)	ZBAD7638	-	-	10-12
SET ZZ BASE AS-EV-SM UA (Scanmarker and screw)	ZBAD7636	-	-	10-12
SET ZZ BASE AS-EV-WSM 3.0 (White Scanmarker and screw)	ZBCD7631	-	-	10-12
SET ZZ BASE AS-EV-WSM 3.6 (White Scanmarker and screw)	ZBCD7632	-	-	10-12
SET ZZ BASE AS-EV-WSM 4.2 (White Scanmarker and screw)	ZBCD7633	-	-	10-12
SET ZZ BASE AS-EV-WSM 4.8 (White Scanmarker and screw)	ZBCD7634	-	-	10-12
SET ZZ BASE AS-EV-WSM 5.4 (White Scanmarker and screw)	ZBCD7635	-	-	10-12
SET ZZ BASE AS-EV-WSM P4.2 (White Scanmarker and screw)	ZBCD7637	-	-	10-12
SET ZZ BASE AS-EV-WSM P4.8 (White Scanmarker and screw)	ZBCD7638	-	-	10-12
SET ZZ BASE AS-EV-WSM UA (White Scanmarker and screw)	ZBCD7636	-	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

ASTRA TECH IMPLANT SYSTEM™ EV

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
ZZ BASE AS-EV-MI 3.0 (Model implant)	MIAG7601	-	-	10-12
ZZ BASE AS-EV-MI 3.6 (Model implant)	MIAG7602	-	-	10-12
ZZ BASE AS-EV-MI 4.2 (Model implant)	MIAG7603	-	-	10-12
ZZ BASE AS-EV-MI 4.8 (Model implant)	MIAG7604	-	-	10-12
ZZ BASE AS-EV-MI 5.4 (Model implant)	MIAG7605	-	-	10-12
ZZ BASE AS-EV-MI P4.2 (Model implant)	MIAG7607	-	-	10-12
ZZ BASE AS-EV-MI P4.8 (Model implant)	MIAG7608	-	-	10-12
ZZ BASE AS-EV-MI UA (Model implant)	MIAG7606	-	-	10-12
ZZ BASE AS-EV-SC 3.0 (Screw)	BSAE7601	25	-	10-12
ZZ BASE AS-EV-SC 3.6 (Screw)	BSAE7602	25	-	10-12
ZZ BASE AS-EV-SC 4.2 (Screw)	BSAE7603	25	-	10-12
ZZ BASE AS-EV-SC 4.8 (Screw)	BSAE7604	25	-	10-12
ZZ BASE AS-EV-SC 5.4 (Screw)	BSAE7605	25	-	10-12
ZZ BASE AS-EV-SC P4.2 (Screw)	BSAE7607	25	-	10-12
ZZ BASE AS-EV-SC P4.8 (Screw)	BSAE7608	25	-	10-12
ZZ BASE AS-EV-SC UA (Screw)	BSAE7606	15	-	10-12
ZZ BASE AS-EV-ZSC 3.0 (Zirconia Screw)	BSAE7631	25	-	10-12
ZZ BASE AS-EV-ZSC 3.6 (Zirconia Screw)	BSAE7632	25	-	10-12
ZZ BASE AS-EV-ZSC 4.2 (Zirconia Screw)	BSAE7633	25	-	10-12
ZZ BASE AS-EV-ZSC 4.8 (Zirconia Screw)	BSAE7634	25	-	10-12
ZZ BASE AS-EV-ZSC 5.4 (Zirconia Screw)	BSAE7635	25	-	10-12
ZZ BASE AS-EV-ZSC P4.2 (Zirconia Screw)	BSAE7637	25	-	10-12
ZZ BASE AS-EV-ZSC P4.8 (Zirconia Screw)	BSAE7638	25	-	10-12
ZZ BASE AS-EV-ZSC UA (Zirconia Screw)	BSAE7636	15	-	10-12

ASTRA TECH OSSEOSPEED™

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE AS-OS-ABCH 3.0 (Raw-Abutment® customizable hex and screw)	BSAF0181	15	-	10-12
SET ZZ BASE AS-OS-ABCH 3.5/4.0 (Raw-Abutment® customizable hex and screw)	BSAF0182	20	-	10-12
SET ZZ BASE AS-OS-ABCH 4.5/5.0 (Raw-Abutment® customizable hex and screw)	BSAF0183	25	-	10-12
SET ZZ BASE AS-OS-ABCH-D14 3.0 (Raw-Abutment® customizable hex D14 and screw)	BSAF0181D14	15	-	10-12
SET ZZ BASE AS-OS-ABCH-D14 3.5/4.0 (Raw-Abutment® customizable hex D14 and screw)	BSAF0182D14	20	-	10-12
SET ZZ BASE AS-OS-ABCH-D14 4.5/5.0 (Raw-Abutment® customizable hex D14 and screw)	BSAF0183D14	25	-	10-12
SET ZZ BASE AS-OS-CCTB 3.0 (Conical cemented titan base non hex and screw)	BSAD0131	15	-	10-12
SET ZZ BASE AS-OS-CCTB 3.5/4.0 (Conical cemented titan base non hex and screw)	BSAD0132	20	-	10-12
SET ZZ BASE AS-OS-CCTB 4.5/5.0 (Conical cemented titan base non hex and screw)	BSAD0133	25	-	10-12
SET ZZ BASE AS-OS-CCTB-L10 3.0 (Conical cemented titan base non hex L10 and screw)	BSAD0131L10	15	-	10-12
SET ZZ BASE AS-OS-CCTB-L10 3.5/4.0 (Conical cemented titan base non hex L10 and screw)	BSAD0132L10	20	-	10-12
SET ZZ BASE AS-OS-CCTB-L10 4.5/5.0 (Conical cemented titan base non hex L10 and screw)	BSAD0133L10	25	-	10-12
SET ZZ BASE AS-OS-CCTB-L20 3.0 (Conical cemented titan base non hex L20 and screw)	BSAD0131L20	15	-	10-12
SET ZZ BASE AS-OS-CCTB-L20 3.5/4.0 (Conical cemented titan base non hex L20 and screw)	BSAD0132L20	20	-	10-12
SET ZZ BASE AS-OS-CCTB-L20 4.5/5.0 (Conical cemented titan base non hex L20 and screw)	BSAD0133L20	25	-	10-12
SET ZZ BASE AS-OS-CCTB-L30 3.0 (Conical cemented titan base non hex L30 and screw)	BSAD0131L30	15	-	10-12
SET ZZ BASE AS-OS-CCTB-L30 3.5/4.0 (Conical cemented titan base non hex L30 and screw)	BSAD0132L30	20	-	10-12
SET ZZ BASE AS-OS-CCTB-L30 4.5/5.0 (Conical cemented titan base non hex L30 and screw)	BSAD0133L30	25	-	10-12
SET ZZ BASE AS-OS-CCTB-L40 3.5/4.0 (Conical cemented titan base non hex L40 and screw)	BSAD0132L40	20	-	10-12
SET ZZ BASE AS-OS-CCTB-L40 4.5/5.0 (Conical cemented titan base non hex L40 and screw)	BSAD0133L40	25	-	10-12
SET ZZ BASE AS-OS-NTB 3.0 (Narrow titan base non hex and screw)	BSAP0131	15	-	10-12
SET ZZ BASE AS-OS-NTB 3.5/4.0 (Narrow titan base non hex and screw)	BSAP0132	20	-	10-12
SET ZZ BASE AS-OS-NTBH 3.0 (Narrow titan base hex and screw)	BSAP0181	15	-	10-12
SET ZZ BASE AS-OS-NTBH 3.5/4.0 (Narrow titan base hex and screw)	BSAP0182	20	-	10-12
SET ZZ BASE AS-OS-PCTBH 3.0 (Parallel cemented titan base hex and screw)	BSAB0181	15	-	10-12
SET ZZ BASE AS-OS-PCTBH 3.5/4.0 (Parallel cemented titan base hex and screw)	BSAB0192	20	-	10-12
SET ZZ BASE AS-OS-PCTBH 4.5/5.0 (Parallel cemented titan base hex and screw)	BSAB0193	25	-	10-12
SET ZZ BASE AS-OS-PCTBH-L20 3.0 (Parallel cemented titan base hex L20 and screw)	BSAB0181L20	15	-	10-12
SET ZZ BASE AS-OS-PCTBH-L20 3.5/4.0 (Parallel cemented titan base hex L20 and screw)	BSAB0182L20	20	-	10-12
SET ZZ BASE AS-OS-PCTBH-L20 4.5/5.0 (Parallel cemented titan base hex L20 and screw)	BSAB0183L20	25	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

ASTRA TECH OSSEOSPEED™

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE AS-OS-PCTBH-L30 3.0 (Parallel cemented titan base hex L30 and screw)	BSAB0181L30	15	-	10-12
SET ZZ BASE AS-OS-PCTBH-L30 3.5/4.0 (Parallel cemented titan base hex L30 and screw)	BSAB0182L30	20	-	10-12
SET ZZ BASE AS-OS-PCTBH-L30 4.5/5.0 (Parallel cemented titan base hex L30 and screw)	BSAB0183L30	25	-	10-12
SET ZZ BASE AS-OS-PCTBH-L40 3.5/4.0 (Parallel cemented titan base hex L40 and screw)	BSAB0182L40	20	-	10-12
SET ZZ BASE AS-OS-PCTBH-L40 4.5/5.0 (Parallel cemented titan base hex L40 and screw)	BSAB0183L40	25	-	10-12
SET ZZ BASE AS-OS-SA 3.0 (Scanalog and screw)	ZBDE0131	15	-	10-12
SET ZZ BASE AS-OS-SA 3.5/4.0 (Scanalog and screw)	ZBDE0132	15	-	10-12
SET ZZ BASE AS-OS-SA 4.5/5.0 (Scanalog and screw)	ZBDE0133	15	-	10-12
SET ZZ BASE AS-OS-SM 3.0 (Scanmarker and screw)	ZBAD0131	-	-	10-12
SET ZZ BASE AS-OS-SM 3.5/4.0 (Scanmarker and screw)	ZBAD0132	-	-	10-12
SET ZZ BASE AS-OS-SM 4.5/5.0 (Scanmarker and screw)	ZBAD0133	-	-	10-12
SET ZZ BASE AS-OS-WSM 3.0 (White Scanmarker and screw)	ZBCD0131	-	-	10-12
SET ZZ BASE AS-OS-WSM 3.5/4.0 (White Scanmarker and screw)	ZBCD0132	-	-	10-12
SET ZZ BASE AS-OS-WSM 4.5/5.0 (White Scanmarker and screw)	ZBCD0133	-	-	10-12
SET ZZ BASE AS-OS-ZZMUA-L15 3.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0131L15	25	15	10-12
SET ZZ BASE AS-OS-ZZMUA-L15 3.5/4.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0132L15	25	20	10-12
SET ZZ BASE AS-OS-ZZMUA-L15 4.5/5.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0133L15	25	25	10-12
SET ZZ BASE AS-OS-ZZMUA-L25 3.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0131L25	25	15	10-12
SET ZZ BASE AS-OS-ZZMUA-L25 3.5/4.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0132L25	25	20	10-12
SET ZZ BASE AS-OS-ZZMUA-L25 4.5/5.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0133L25	25	25	10-12
SET ZZ BASE AS-OS-ZZMUA-L35 3.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0131L35	25	15	10-12
SET ZZ BASE AS-OS-ZZMUA-L35 3.5/4.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0132L35	25	20	10-12
SET ZZ BASE AS-OS-ZZMUA-L35 4.5/5.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0133L35	25	25	10-12
SET ZZ BASE AS-OS-ZZMUA-L45 3.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0131L45	25	15	10-12
SET ZZ BASE AS-OS-ZZMUA-L45 3.5/4.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0132L45	25	20	10-12
SET ZZ BASE AS-OS-ZZMUA-L45 4.5/5.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0133L45	25	25	10-12
ZZ BASE AS-OS-MI 3.0 (Model implant)	MIAG0101	-	-	10-12
ZZ BASE AS-OS-MI 3.5/4.0 (Model implant)	MIAG0102	-	-	10-12
ZZ BASE AS-OS-MI 4.5/5.0 (Model implant)	MIAG0103	-	-	10-12
ZZ BASE AS-OS-SC 3.0 (Screw)	BSAE0101	15	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

ASTRA TECH OSSEOSPEED™

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
ZZ BASE AS-OS-SC 3.5/4.0 (Screw)	BSAE0102	20	-	10-12
ZZ BASE AS-OS-SC 4.5/5.0 (Screw)	BSAE0103	25	-	10-12
ZZ BASE AS-OS-SC-G 3.0 (Screw Gold)	BSCE0101	15	-	10-12
ZZ BASE AS-OS-SC-G 3.5/4.0 (Screw Gold)	BSCE0102	20	-	10-12
ZZ BASE AS-OS-SC-G 4.5/5.0 (Screw Gold)	BSCE0103	25	-	10-12
ZZ BASE AS-OS-ZSC 3.5/4.0 (Zirconia screw)	BSAE0132	20	-	10-12
ZZ BASE AS-OS-ZSC 4.5/5.0 (Zirconia screw)	BSAE0133	25	-	10-12

ASTRA TECH UNIABUTMENT®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE AS-UA-CCTB 20° (Conical cemented titan base non hex and screw)	BSAD0731	15	-	10-12
SET ZZ BASE AS-UA-CCTB 45° (Conical cemented titan base non hex and screw)	BSAD0732	15	-	10-12
SET ZZ BASE AS-UA-SM 20° (Scanmarker and screw)	ZBAD0731	-	-	10-12
SET ZZ BASE AS-UA-SM 45° (Scanmarker and screw)	ZBAD0732	-	-	10-12
SET ZZ BASE AS-UA-WSM 20° (White Scanmarker and screw)	ZBCD0731	-	-	10-12
SET ZZ BASE AS-UA-WSM 45° (White Scanmarker and screw)	ZBCD0732	-	-	10-12
ZZ BASE AS-UA-MI 20° (Model implant)	MIAG0701	-	-	10-12
ZZ BASE AS-UA-MI 45° (Model implant)	MIAG0702	-	-	10-12
ZZ BASE AS-UA-SC 20° (Screw)	BSAE0701	15	-	10-12
ZZ BASE AS-UA-SC 45° (Screw)	BSAE0702	15	-	10-12
ZZ BASE AS-UA-SC-G 20° (Screw Gold)	BSCE0701	15	-	10-12
ZZ BASE AS-UA-SC-G 45° (Screw Gold)	BSCE0702	15	-	10-12
ZZ BASE AS-UA-ZSC 20° (Zirconia screw)	BSAE0731	15	-	10-12
ZZ BASE AS-UA-ZSC 45° (Zirconia screw)	BSAE0732	15	-	10-12

A-Z® VL

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE AZ-VL-ABCH-ID 3.2 (Raw-Abutment® customizable hex ID and screw)	BSAF5694	-	-	10-12
SET ZZ BASE AZ-VL-ABCH-ID 3.8 (Raw-Abutment® customizable hex ID and screw)	BSAF5691	-	-	10-12
SET ZZ BASE AZ-VL-ABCH-ID 4.4 (Raw-Abutment® customizable hex ID and screw)	BSAF5692	-	-	10-12
SET ZZ BASE AZ-VL-ABCH-ID 5.5 (Raw-Abutment® customizable hex ID and screw)	BSAF5693	-	-	10-12
SET ZZ BASE AZ-VL-CCTB 3.2 (Conical cemented titan base non hex and screw)	BSAD5634	-	-	10-12
SET ZZ BASE AZ-VL-CCTB 3.8 (Conical cemented titan base non hex and screw)	BSAD5631	-	-	10-12
SET ZZ BASE AZ-VL-CCTB 4.4 (Conical cemented titan base non hex and screw)	BSAD5632	-	-	10-12
SET ZZ BASE AZ-VL-CCTB 5.5 (Conical cemented titan base non hex and screw)	BSAD5633	-	-	10-12
SET ZZ BASE AZ-VL-PCTBH 3.2 (Parallel cemented titan base hex and screw)	BSAB5684	-	-	10-12
SET ZZ BASE AZ-VL-PCTBH 3.8 (Parallel cemented titan base hex and screw)	BSAB5681	-	-	10-12
SET ZZ BASE AZ-VL-PCTBH 4.4 (Parallel cemented titan base hex and screw)	BSAB5682	-	-	10-12
SET ZZ BASE AZ-VL-PCTBH 5.5 (Parallel cemented titan base hex and screw)	BSAB5683	-	-	10-12
SET ZZ BASE AZ-VL-SM 3.2 (Scanmarker and screw)	ZBAD5634	-	-	10-12
SET ZZ BASE AZ-VL-SM 3.8 (Scanmarker and screw)	ZBAD5631	-	-	10-12
SET ZZ BASE AZ-VL-SM 4.4 (Scanmarker and screw)	ZBAD5632	-	-	10-12
SET ZZ BASE AZ-VL-SM 5.5 (Scanmarker and screw)	ZBAD5633	-	-	10-12
SET ZZ BASE AZ-VL-WSM 3.2 (White Scanmarker and screw)	ZBCD5634	-	-	10-12
SET ZZ BASE AZ-VL-WSM 3.8 (White Scanmarker and screw)	ZBCD5631	-	-	10-12
SET ZZ BASE AZ-VL-WSM 4.4 (White Scanmarker and screw)	ZBCD5632	-	-	10-12
SET ZZ BASE AZ-VL-WSM 5.5 (White Scanmarker and screw)	ZBCD5633	-	-	10-12
ZZ BASE AZ-VL-SC 3.2 (Screw)	BSAE5604	-	-	10-12
ZZ BASE AZ-VL-SC 3.8 (Screw)	BSAE5601	-	-	10-12
ZZ BASE AZ-VL-SC 4.4 (Screw)	BSAE5602	-	-	10-12
ZZ BASE AZ-VL-SC 5.5 (Screw)	BSAE5603	-	-	10-12
ZZ BASE AZ-VL-ZSC 3.2 (Zirconia screw)	BSAE5634	-	-	10-12
ZZ BASE AZ-VL-ZSC 3.8 (Zirconia screw)	BSAE5631	-	-	10-12
ZZ BASE AZ-VL-ZSC 4.4 (Zirconia screw)	BSAE5632	-	-	10-12
ZZ BASE AZ-VL-ZSC 5.5 (Zirconia screw)	BSAE5633	-	-	10-12

BEGO MINI

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE BE-MI-ABCH-ID 2.7-3.1 (Raw-Abutment® customizable hex ID and screw)	BSAF7591	25	-	10-12
SET ZZ BASE BE-MI-CCTB 2.7-3.1 (Conical cemented titan base non hex and screw)	BSAD7531	25	-	10-12
SET ZZ BASE BE-MI-PCTBH 2.7-3.1 (Parallel cemented titan base hex and screw)	BSAB7581	25	-	10-12
SET ZZ BASE BE-MI-SM 2.7-3.1 (Scanmarker and screw)	ZBAD7531	-	-	10-12
ZZ BASE BE-MI-SC 2.7-3.1 (Screw)	BSAE7501	25	-	10-12
ZZ BASE BE-MI-ZSC 2.7-3.1 (Zirconia Screw)	BSAE7531	25	-	10-12

BEGO SEMADOS® S/RI/SC/SCX/RS/RXS-LINE®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE BE-SE-S/RI-ABCH 3.25 (Raw-Abutment® customizable hex and screw)	BSAF5081	30	-	10-12
SET ZZ BASE BE-SE-S/RI-ABCH 3.75 (Raw-Abutment® customizable hex and screw)	BSAF5082	30	-	10-12
SET ZZ BASE BE-SE-S/RI-ABCH 4.1 (Raw-Abutment® customizable hex and screw)	BSAF5083	30	-	10-12
SET ZZ BASE BE-SE-S/RI-ABCH 4.5 (Raw-Abutment® customizable hex and screw)	BSAF5084	30	-	10-12
SET ZZ BASE BE-SE-S/RI-ABCH 5.5 (Raw-Abutment® customizable hex and screw)	BSAF5085	30	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB 3.0 (Conical cemented titan base non hex and screw)	BSAD5037	20	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB 3.25 (Conical cemented titan base non hex and screw)	BSAD5031	30	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB 3.75 (Conical cemented titan base non hex and screw)	BSAD5032	30	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB 4.1 (Conical cemented titan base non hex and screw)	BSAD5033	30	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB 4.5 (Conical cemented titan base non hex and screw)	BSAD5034	30	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB 5.5 (Conical cemented titan base non hex and screw)	BSAD5035	30	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB MP (Conical cemented titan base non hex and screw)	BSAD5036	30	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB MP (Conical cemented titan base non hex and screw)	BSAD5046	30	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L10 4.5 (Conical cemented titan base non hex L10 and screw)	BSAD5034L10	30	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L20 4.1 (Conical cemented titan base non hex L20 and screw)	BSAD5033L20	30	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L20 4.5 (Conical cemented titan base non hex L20 and screw)	BSAD5034L20	30	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L20 5.5 (Conical cemented titan base non hex L20 and screw)	BSAD5035L20	30	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L30 3.25 (Conical cemented titan base non hex L30 and screw)	BSAD5031L30	30	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L30 3.75 (Conical cemented titan base non hex L30 and screw)	BSAD5032L30	30	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L30 4.1 (Conical cemented titan base non hex L30 and screw)	BSAD5033L30	30	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L30 4.5 (Conical cemented titan base non hex L30 and screw)	BSAD5034L30	30	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L30 5.5 (Conical cemented titan base non hex L30 and screw)	BSAD5035L30	30	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L40 3.25 (Conical cemented titan base non hex L40 and screw)	BSAD5031L40	30	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L40 3.75 (Conical cemented titan base non hex L40 and screw)	BSAD5032L40	30	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L40 4.1 (Conical cemented titan base non hex L40 and screw)	BSAD5033L40	30	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L40 4.5 (Conical cemented titan base non hex L40 and screw)	BSAD5034L40	30	-	10-12
SET ZZ BASE BE-SE-S/RI-CCTB-L40 5.5 (Conical cemented titan base non hex L40 and screw)	BSAD5035L40	30	-	10-12
SET ZZ BASE BE-SE-S/RI-NTB 3.0 (Narrow titan base non hex and screw)	BSAP5037	20	-	10-12
SET ZZ BASE BE-SE-S/RI-NTB 3.25 (Narrow titan base non hex and screw)	BSAP5031	30	-	10-12
SET ZZ BASE BE-SE-S/RI-NTB 3.75 (Narrow titan base non hex and screw)	BSAP5032	30	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

BEGO SEMADOS® S/RI/SC/SCX/RS/RXS-LINE®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE BE-SE-S/RI-NTB 4.1 (Narrow titan base non hex and screw)	BSAP5033	30	-	10-12
SET ZZ BASE BE-SE-S/RI-NTB 4.5 (Narrow titan base non hex and screw)	BSAP5034	30	-	10-12
SET ZZ BASE BE-SE-S/RI-NTBH 3.0 (Narrow titan base hex and screw)	BSAP5087	20	-	10-12
SET ZZ BASE BE-SE-S/RI-NTBH 3.25 (Narrow titan base hex and screw)	BSAP5081	30	-	10-12
SET ZZ BASE BE-SE-S/RI-NTBH 3.75 (Narrow titan base hex and screw)	BSAP5082	30	-	10-12
SET ZZ BASE BE-SE-S/RI-NTBH 4.1 (Narrow titan base hex and screw)	BSAP5083	30	-	10-12
SET ZZ BASE BE-SE-S/RI-NTBH 4.5 (Narrow titan base hex and screw)	BSAP5084	30	-	10-12
SET ZZ BASE BE-SE-S/RI-PCTBH 3.0 (Parallel cemented titan base hex and screw)	BSAB5087	20	-	10-12
SET ZZ BASE BE-SE-S/RI-PCTBH 3.25 (Parallel cemented titan base hex and screw)	BSAB5081	30	-	10-12
SET ZZ BASE BE-SE-S/RI-PCTBH 3.75 (Parallel cemented titan base hex and screw)	BSAB5082	30	-	10-12
SET ZZ BASE BE-SE-S/RI-PCTBH 4.1 (Parallel cemented titan base hex and screw)	BSAB5083	30	-	10-12
SET ZZ BASE BE-SE-S/RI-PCTBH 4.5 (Parallel cemented titan base hex and screw)	BSAB5084	30	-	10-12
SET ZZ BASE BE-SE-S/RI-PCTBH 5.5 (Parallel cemented titan base hex and screw)	BSAB5085	30	-	10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L20 3.25 (Parallel cemented titan base hex L20 and screw)	BSAB5091L20	30	-	10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L20 3.75 (Parallel cemented titan base hex L20 and screw)	BSAB5092L20	30	-	10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L20 4.1 (Parallel cemented titan base hex L20 and screw)	BSAB5093L20	30	-	10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L20 4.5 (Parallel cemented titan base hex L20 and screw)	BSAB5094L20	30	-	10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L20 5.5 (Parallel cemented titan base hex L20 and screw)	BSAB5095L20	30	-	10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L30 3.25 (Parallel cemented titan base hex L30 and screw)	BSAB5091L30	30	-	10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L30 3.75 (Parallel cemented titan base hex L30 and screw)	BSAB5092L30	30	-	10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L30 4.1 (Parallel cemented titan base hex L30 and screw)	BSAB5093L30	30	-	10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L30 4.5 (Parallel cemented titan base hex L30 and screw)	BSAB5094L30	30	-	10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L30 5.5 (Parallel cemented titan base hex L30 and screw)	BSAB5095L30	30	-	10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L40 3.25 (Parallel cemented titan base hex L40 and screw)	BSAB5091L40	30	-	10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L40 3.75 (Parallel cemented titan base hex L40 and screw)	BSAB5092L40	30	-	10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L40 4.1 (Parallel cemented titan base hex L40 and screw)	BSAB5093L40	30	-	10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L40 4.5 (Parallel cemented titan base hex L40 and screw)	BSAB5094L40	30	-	10-12
SET ZZ BASE BE-SE-S/RI-PCTBH-L40 5.5 (Parallel cemented titan base hex L40 and screw)	BSAB5095L40	30	-	10-12
SET ZZ BASE BE-SE-S/RI-SM 3.0 (Scanmarker and screw)	ZBAD5037	-	-	10-12
SET ZZ BASE BE-SE-S/RI-SM 3.25 (Scanmarker and screw)	ZBAD5031	-	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

BEGO SEMADOS® S/RI/SC/SCX/RS/R SX-LINE®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE BE-SE-S/RI-SM 3.75 (Scanmarker and screw)	ZBAD5032	-	-	10-12
SET ZZ BASE BE-SE-S/RI-SM 4.1 (Scanmarker and screw)	ZBAD5033	-	-	10-12
SET ZZ BASE BE-SE-S/RI-SM 4.5 (Scanmarker and screw)	ZBAD5034	-	-	10-12
SET ZZ BASE BE-SE-S/RI-SM 5.5 (Scanmarker and screw)	ZBAD5035	-	-	10-12
SET ZZ BASE BE-SE-S/RI-SM MP (Scanmarker and screw)	ZBAD5036	-	-	10-12
SET ZZ BASE BE-SE-S/RI-WSM 3.25 (White Scanmarker and screw)	ZBCD5031	-	-	10-12
SET ZZ BASE BE-SE-S/RI-WSM 3.75 (White Scanmarker and screw)	ZBCD5032	-	-	10-12
SET ZZ BASE BE-SE-S/RI-WSM 4.1 (White Scanmarker and screw)	ZBCD5033	-	-	10-12
SET ZZ BASE BE-SE-S/RI-WSM 4.5 (White Scanmarker and screw)	ZBCD5034	-	-	10-12
SET ZZ BASE BE-SE-S/RI-WSM 5.5 (White Scanmarker and screw)	ZBCD5035	-	-	10-12
SET ZZ BASE BE-SE-S/RI-WSM MP (White Scanmarker and screw)	ZBCD5036	-	-	10-12
ZZ BASE BE-SE-S/RI-SC 3.0 (Screw)	BSAE5007	20	-	10-12
ZZ BASE BE-SE-S/RI-SC 3.25 (Screw)	BSAE5001	30	-	10-12
ZZ BASE BE-SE-S/RI-SC 3.75 (Screw)	BSAE5002	30	-	10-12
ZZ BASE BE-SE-S/RI-SC 4.1 (Screw)	BSAE5003	30	-	10-12
ZZ BASE BE-SE-S/RI-SC 4.5 (Screw)	BSAE5004	30	-	10-12
ZZ BASE BE-SE-S/RI-SC 5.5 (Screw)	BSAE5005	30	-	10-12
ZZ BASE BE-SE-S/RI-SC MP (Screw)	BSAE5006	30	-	10-12
ZZ BASE BE-SE-S/RI-ZSC 3.0 (Zirconia screw)	BSAE5037	20	-	10-12
ZZ BASE BE-SE-S/RI-ZSC 3.25 (Zirconia screw)	BSAE5031	30	-	10-12
ZZ BASE BE-SE-S/RI-ZSC 3.75 (Zirconia screw)	BSAE5032	30	-	10-12
ZZ BASE BE-SE-S/RI-ZSC 4.1 (Zirconia screw)	BSAE5033	30	-	10-12
ZZ BASE BE-SE-S/RI-ZSC 4.5 (Zirconia screw)	BSAE5034	30	-	10-12
ZZ BASE BE-SE-S/RI-ZSC 5.5 (Zirconia screw)	BSAE5035	30	-	10-12
ZZ BASE BE-SE-S/RI-ZSC MP (Zirconia screw)	BSAE5036	30	-	10-12

BIOHORIZONS® EXTERNAL

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE BH-EX-ABCH-ID 3.5 (Raw-Abutment® customizable hex ID and screw)	BSAF3891	30	-	10-12
SET ZZ BASE BH-EX-ABCH-ID 4.0 (Raw-Abutment® customizable hex ID and screw)	BSAF3892	30	-	10-12
SET ZZ BASE BH-EX-ABCH-ID 5.0 (Raw-Abutment® customizable hex ID and screw)	BSAF3893	30	-	10-12
SET ZZ BASE BH-EX-ABCH-ID 6.0 (Raw-Abutment® customizable hex ID and screw)	BSAF3894	30	-	10-12
SET ZZ BASE BH-EX-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD3831	30	-	10-12
SET ZZ BASE BH-EX-CCTB 4.0 (Conical cemented titan base non hex and screw)	BSAD3832	30	-	10-12
SET ZZ BASE BH-EX-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD3833	30	-	10-12
SET ZZ BASE BH-EX-CCTB 6.0 (Conical cemented titan base non hex and screw)	BSAD3834	30	-	10-12
SET ZZ BASE BH-EX-NTB 3.5 (Narrow titan base non hex and screw)	BSAP3831	30	-	10-12
SET ZZ BASE BH-EX-NTBH 3.5 (Narrow titan base hex and screw)	BSAP3881	30	-	10-12
SET ZZ BASE BH-EX-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB3881	30	-	10-12
SET ZZ BASE BH-EX-PCTBH 4.0 (Parallel cemented titan base hex and screw)	BSAB3882	30	-	10-12
SET ZZ BASE BH-EX-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB3883	30	-	10-12
SET ZZ BASE BH-EX-PCTBH 6.0 (Parallel cemented titan base hex and screw)	BSAB3884	30	-	10-12
SET ZZ BASE BH-EX-SM 3.5 (Scanmarker and screw)	ZBAD3831	-	-	10-12
SET ZZ BASE BH-EX-SM 4.0 (Scanmarker and screw)	ZBAD3832	-	-	10-12
SET ZZ BASE BH-EX-SM 5.0 (Scanmarker and screw)	ZBAD3833	-	-	10-12
SET ZZ BASE BH-EX-SM 6.0 (Scanmarker and screw)	ZBAD3834	-	-	10-12
SET ZZ BASE BH-EX-WSM 3.5 (White Scanmarker and screw)	ZBCD3831	-	-	10-12
SET ZZ BASE BH-EX-WSM 4.0 (White Scanmarker and screw)	ZBCD3832	-	-	10-12
SET ZZ BASE BH-EX-WSM 5.0 (White Scanmarker and screw)	ZBCD3833	-	-	10-12
SET ZZ BASE BH-EX-WSM 6.0 (White Scanmarker and screw)	ZBCD3834	-	-	10-12
ZZ BASE BH-EX-SC 3.5 (Screw)	BSAE3801	30	-	10-12
ZZ BASE BH-EX-SC 4.0 (Screw)	BSAE3802	30	-	10-12
ZZ BASE BH-EX-SC 5.0 (Screw)	BSAE3803	30	-	10-12
ZZ BASE BH-EX-SC 6.0 (Screw)	BSAE3804	30	-	10-12
ZZ BASE BH-EX-ZSC 3.5 (Zirconia screw)	BSAE3831	30	-	10-12
ZZ BASE BH-EX-ZSC 4.0 (Zirconia screw)	BSAE3832	30	-	10-12
ZZ BASE BH-EX-ZSC 5.0 (Zirconia screw)	BSAE3833	30	-	10-12
ZZ BASE BH-EX-ZSC 6.0 (Zirconia screw)	BSAE3834	30	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

BIOHORIZONS® INTERNAL

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE BH-INT-ABCH 3.0 (Raw-Abutment® customizable hex and screw)	BSAF3781	30	-	10-12
SET ZZ BASE BH-INT-ABCH 3.5 (Raw-Abutment® customizable hex and screw)	BSAF3782	30	-	10-12
SET ZZ BASE BH-INT-ABCH 4.5 (Raw-Abutment® customizable hex and screw)	BSAF3783	30	-	10-12
SET ZZ BASE BH-INT-ABCH 5.7 (Raw-Abutment® customizable hex and screw)	BSAF3784	30	-	10-12
SET ZZ BASE BH-INT-CCTB 3.0 (Conical cemented titan base non hex and screw)	BSAD3731	30	-	10-12
SET ZZ BASE BH-INT-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD3732	30	-	10-12
SET ZZ BASE BH-INT-CCTB 4.5 (Conical cemented titan base non hex and screw)	BSAD3733	30	-	10-12
SET ZZ BASE BH-INT-CCTB 5.7 (Conical cemented titan base non hex and screw)	BSAD3734	30	-	10-12
SET ZZ BASE BH-INT-NTB 3.0 (Narrow titan base non hex and screw)	BSAP3731	30	-	10-12
SET ZZ BASE BH-INT-NTB 3.5 (Narrow titan base non hex and screw)	BSAP3732	30	-	10-12
SET ZZ BASE BH-INT-NTBH 3.0 (Narrow titan base hex and screw)	BSAP3781	30	-	10-12
SET ZZ BASE BH-INT-NTBH 3.5 (Narrow titan base hex and screw)	BSAP3782	30	-	10-12
SET ZZ BASE BH-INT-PCTBH 3.0 (Parallel cemented titan base hex and screw)	BSAB3781	30	-	10-12
SET ZZ BASE BH-INT-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB3782	30	-	10-12
SET ZZ BASE BH-INT-PCTBH 4.5 (Parallel cemented titan base hex and screw)	BSAB3783	30	-	10-12
SET ZZ BASE BH-INT-PCTBH 5.7 (Parallel cemented titan base hex and screw)	BSAB3784	30	-	10-12
SET ZZ BASE BH-INT-SM 3.0 (Scanmarker and screw)	ZBAD3731	-	-	10-12
SET ZZ BASE BH-INT-SM 3.5 (Scanmarker and screw)	ZBAD3732	-	-	10-12
SET ZZ BASE BH-INT-SM 4.5 (Scanmarker and screw)	ZBAD3733	-	-	10-12
SET ZZ BASE BH-INT-SM 5.7 (Scanmarker and screw)	ZBAD3734	-	-	10-12
SET ZZ BASE BH-INT-WSM 3.0 (White Scanmarker and screw)	ZBCD3731	-	-	10-12
SET ZZ BASE BH-INT-WSM 3.5 (White Scanmarker and screw)	ZBCD3732	-	-	10-12
SET ZZ BASE BH-INT-WSM 4.5 (White Scanmarker and screw)	ZBCD3733	-	-	10-12
SET ZZ BASE BH-INT-WSM 5.7 (White Scanmarker and screw)	ZBCD3734	-	-	10-12
ZZ BASE BH-INT-SC 3.0 (Screw)	BSAE3701	30	-	10-12
ZZ BASE BH-INT-SC 3.5 (Screw)	BSAE3702	30	-	10-12
ZZ BASE BH-INT-SC 4.5 (Screw)	BSAE3703	30	-	10-12
ZZ BASE BH-INT-SC 5.7 (Screw)	BSAE3704	30	-	10-12
ZZ BASE BH-INT-ZSC 3.0 (Zirconia screw)	BSAE3731	30	-	10-12
ZZ BASE BH-INT-ZSC 3.5 (Zirconia screw)	BSAE3732	30	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

BIOHORIZONS® INTERNAL

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
ZZ BASE BH-INT-ZSC 4.5 (Zirconia screw)	BSAE3733	30	-	10-12
ZZ BASE BH-INT-ZSC 5.7 (Zirconia screw)	BSAE3734	30	-	10-12

BIOHORIZONS® MULTI UNIT ABUTMENT

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE BH-MUA-CCTB (Conical cemented titan base non hex and screw)	BSAD7731	15	-	10-12
SET ZZ BASE BH-MUA-SM (Scanmarker and screw)	ZBAD7731	-	-	10-12
ZZ BASE BH-MUA-SC (Screw)	BSAE7701	15	-	10-12
ZZ BASE BH-MUA-ZSC (Zirconia screw)	BSAE7731	15	-	10-12

BIOMET 3I™ LOW PROFILE ABUTMENT

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE B-LPA-CCTB (Conical cemented titan base non hex and screw)	BSAD6431	20	-	10-12
SET ZZ BASE B-LPA-CCTB (Conical cemented titan base non hex and screw)	BSAD6441	20	-	10-12
SET ZZ BASE B-LPA-SM (Scanmarker and screw)	ZBAD6431	-	-	10-12
SET ZZ BASE B-LPA-WSM (White Scanmarker and screw)	ZBCD6431	-	-	10-12
ZZ BASE B-LPA-MI (Model implant)	MIAG6401	-	-	10-12
ZZ BASE B-LPA-SC (Screw)	BSAE6401	20	-	10-12
ZZ BASE B-LPA-ZSC (Zirconia Screw)	BSAE6431	20	-	10-12

BIOMET 3I™ OSSEOTITE® CERTAIN®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE B-C-ABCH-ID 3.4 (Raw-Abutment® customizable hex ID and screw)	BSAF0891	20	-	10-12
SET ZZ BASE B-C-ABCH-ID 4.1 (Raw-Abutment® customizable hex ID and screw)	BSAF0892	20	-	10-12
SET ZZ BASE B-C-ABCH-ID 5.0 (Raw-Abutment® customizable hex ID and screw)	BSAF0893	20	-	10-12
SET ZZ BASE B-C-ABCH-ID 6.0 (Raw-Abutment® customizable hex ID and screw)	BSAF0894	20	-	10-12
SET ZZ BASE B-C-CCTB 3.4 (Conical cemented titan base non hex and screw)	BSAD0831	20	-	10-12
SET ZZ BASE B-C-CCTB 3.4 (Conical cemented titan base non hex and thick screw)	BSAD0841	20	-	10-12
SET ZZ BASE B-C-CCTB 4.1 (Conical cemented titan base non hex and screw)	BSAD0832	20	-	10-12
SET ZZ BASE B-C-CCTB 4.1 (Conical cemented titan base non hex and thick screw)	BSAD0842	20	-	10-12
SET ZZ BASE B-C-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD0833	20	-	10-12
SET ZZ BASE B-C-CCTB 5.0 (Conical cemented titan base non hex and thick screw)	BSAD0843	20	-	10-12
SET ZZ BASE B-C-CCTB 6.0 (Conical cemented titan base non hex and screw)	BSAD0834	20	-	10-12
SET ZZ BASE B-C-CCTB 6.0 (Conical cemented titan base non hex and thick screw)	BSAD0844	20	-	10-12
SET ZZ BASE B-C-CCTB-L10 3.4 (Conical cemented titan base non hex L10 and thick screw)	BSAD0831L10	20	-	10-12
SET ZZ BASE B-C-CCTB-L10 4.1 (Conical cemented titan base non hex L10 and thick screw)	BSAD0832L10	20	-	10-12
SET ZZ BASE B-C-CCTB-L10 5.0 (Conical cemented titan base non hex L10 and thick screw)	BSAD0833L10	20	-	10-12
SET ZZ BASE B-C-CCTB-L10 6.0 (Conical cemented titan base non hex L10 and thick screw)	BSAD0834L10	20	-	10-12
SET ZZ BASE B-C-CCTB-L20 3.4 (Conical cemented titan base non hex L20 and thick screw)	BSAD0831L20	20	-	10-12
SET ZZ BASE B-C-CCTB-L20 4.1 (Conical cemented titan base non hex L20 and thick screw)	BSAD0832L20	20	-	10-12
SET ZZ BASE B-C-CCTB-L20 5.0 (Conical cemented titan base non hex L20 and thick screw)	BSAD0833L20	20	-	10-12
SET ZZ BASE B-C-CCTB-L20 6.0 (Conical cemented titan base non hex L20 and thick screw)	BSAD0834L20	20	-	10-12
SET ZZ BASE B-C-CCTB-L30 3.4 (Conical cemented titan base non hex L30 and thick screw)	BSAD0831L30	20	-	10-12
SET ZZ BASE B-C-CCTB-L30 4.1 (Conical cemented titan base non hex L30 and thick screw)	BSAD0832L30	20	-	10-12
SET ZZ BASE B-C-CCTB-L30 5.0 (Conical cemented titan base non hex L30 and thick screw)	BSAD0833L30	20	-	10-12
SET ZZ BASE B-C-CCTB-L30 6.0 (Conical cemented titan base non hex L30 and thick screw)	BSAD0834L30	20	-	10-12
SET ZZ BASE B-C-CCTB-L40 3.4 (Conical cemented titan base non hex L40 and thick screw)	BSAD0831L40	20	-	10-12
SET ZZ BASE B-C-CCTB-L40 4.1 (Conical cemented titan base non hex L40 and thick screw)	BSAD0832L40	20	-	10-12
SET ZZ BASE B-C-CCTB-L40 5.0 (Conical cemented titan base non hex L40 and thick screw)	BSAD0833L40	20	-	10-12
SET ZZ BASE B-C-CCTB-L40 6.0 (Conical cemented titan base non hex L40 and thick screw)	BSAD0834L40	20	-	10-12
SET ZZ BASE B-C-NTB 3.4 (Narrow titan base non hex and thick screw)	BSAP0831	20	-	10-12
SET ZZ BASE B-C-NTB 4.1 (Narrow titan base non hex and thick screw)	BSAP0832	20	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

BIOMET 3I™ OSSEOTITE® CERTAIN®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE B-C-NTBH 3.4 (Narrow titan base hex and screw)	BSAP0881	20	-	10-12
SET ZZ BASE B-C-NTBH 4.1 (Narrow titan base hex and screw)	BSAP0882	20	-	10-12
SET ZZ BASE B-C-PCTBH 3.4 (Parallel cemented titan base hex and screw)	BSAB0881	20	-	10-12
SET ZZ BASE B-C-PCTBH 4.1 (Parallel cemented titan base hex and screw)	BSAB0882	20	-	10-12
SET ZZ BASE B-C-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB0883	20	-	10-12
SET ZZ BASE B-C-PCTBH 6.0 (Parallel cemented titan base hex and screw)	BSAB0884	20	-	10-12
SET ZZ BASE B-C-PCTBH-L10 3.4 (Parallel cemented titan base hex L10 and screw)	BSAB0891L10	20	-	10-12
SET ZZ BASE B-C-PCTBH-L10 4.1 (Parallel cemented titan base hex L10 and screw)	BSAB0892L10	20	-	10-12
SET ZZ BASE B-C-PCTBH-L10 5.0 (Parallel cemented titan base hex L10 and screw)	BSAB0893L10	20	-	10-12
SET ZZ BASE B-C-PCTBH-L10 6.0 (Parallel cemented titan base hex L10 and screw)	BSAB0894L10	20	-	10-12
SET ZZ BASE B-C-PCTBH-L20 3.4 (Parallel cemented titan base hex L20 and screw)	BSAB0891L20	20	-	10-12
SET ZZ BASE B-C-PCTBH-L20 4.1 (Parallel cemented titan base hex L20 and screw)	BSAB0892L20	20	-	10-12
SET ZZ BASE B-C-PCTBH-L20 5.0 (Parallel cemented titan base hex L20 and screw)	BSAB0893L20	20	-	10-12
SET ZZ BASE B-C-PCTBH-L20 6.0 (Parallel cemented titan base hex L20 and screw)	BSAB0894L20	20	-	10-12
SET ZZ BASE B-C-PCTBH-L30 3.4 (Parallel cemented titan base hex L30 and screw)	BSAB0891L30	20	-	10-12
SET ZZ BASE B-C-PCTBH-L30 4.1 (Parallel cemented titan base hex L30 and screw)	BSAB0892L30	20	-	10-12
SET ZZ BASE B-C-PCTBH-L30 5.0 (Parallel cemented titan base hex L30 and screw)	BSAB0893L30	20	-	10-12
SET ZZ BASE B-C-PCTBH-L30 6.0 (Parallel cemented titan base hex L30 and screw)	BSAB0894L30	20	-	10-12
SET ZZ BASE B-C-PCTBH-L40 3.4 (Parallel cemented titan base hex L40 and screw)	BSAB0891L40	20	-	10-12
SET ZZ BASE B-C-PCTBH-L40 4.1 (Parallel cemented titan base hex L40 and screw)	BSAB0892L40	20	-	10-12
SET ZZ BASE B-C-PCTBH-L40 5.0 (Parallel cemented titan base hex L40 and screw)	BSAB0893L40	20	-	10-12
SET ZZ BASE B-C-PCTBH-L40 6.0 (Parallel cemented titan base hex L40 and screw)	BSAB0894L40	20	-	10-12
SET ZZ BASE B-C-SA 3.4 (Scananalog and screw)	ZBDE0831	15	-	10-12
SET ZZ BASE B-C-SA 4.1 (Scananalog and screw)	ZBDE0832	15	-	10-12
SET ZZ BASE B-C-SA 5.0 (Scananalog and screw)	ZBDE0833	15	-	10-12
SET ZZ BASE B-C-SM 3.4 (Scanmarker and screw)	ZBAD0831	-	-	10-12
SET ZZ BASE B-C-SM 4.1 (Scanmarker and screw)	ZBAD0832	-	-	10-12
SET ZZ BASE B-C-SM 5.0 (Scanmarker and screw)	ZBAD0833	-	-	10-12
SET ZZ BASE B-C-SM 6.0 (Scanmarker and screw)	ZBAD0834	-	-	10-12
SET ZZ BASE B-C-WSM 3.4 (White Scanmarker and screw)	ZBCD0831	-	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

BIOMET 3I™ OSSEOTITE® CERTAIN®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE B-C-WSM 4.1 (White Scanmarker and screw)	ZBCD0832	-	-	10-12
SET ZZ BASE B-C-WSM 5.0 (White Scanmarker and screw)	ZBCD0833	-	-	10-12
SET ZZ BASE B-C-WSM 6.0 (White Scanmarker and screw)	ZBCD0834	-	-	10-12
SET ZZ BASE B-C-ZZMUA-L15 3.4 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0831L15	25	20	10-12
SET ZZ BASE B-C-ZZMUA-L15 4.1 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0832L15	25	20	10-12
SET ZZ BASE B-C-ZZMUA-L15 5.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0833L15	25	20	10-12
SET ZZ BASE B-C-ZZMUA-L25 3.4 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0831L25	25	20	10-12
SET ZZ BASE B-C-ZZMUA-L25 4.1 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0832L25	25	20	10-12
SET ZZ BASE B-C-ZZMUA-L25 5.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0833L25	25	20	10-12
SET ZZ BASE B-C-ZZMUA-L35 3.4 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0831L35	25	20	10-12
SET ZZ BASE B-C-ZZMUA-L35 4.1 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0832L35	25	20	10-12
SET ZZ BASE B-C-ZZMUA-L35 5.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0833L35	25	20	10-12
SET ZZ BASE B-C-ZZMUA-L45 3.4 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0831L45	25	20	10-12
SET ZZ BASE B-C-ZZMUA-L45 4.1 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0832L45	25	20	10-12
ZZ BASE B-C-MI 3.4 (Model implant)	MIAG0801	-	-	10-12
ZZ BASE B-C-MI 4.1 (Model implant)	MIAG0802	-	-	10-12
ZZ BASE B-C-MI 5.0 (Model implant)	MIAG0803	-	-	10-12
ZZ BASE B-C-MI 6.0 (Model implant)	MIAG0804	-	-	10-12
ZZ BASE B-C-SC 3.4 (Screw)	BSAE0801	20	-	10-12
ZZ BASE B-C-SC 4.1 (Screw)	BSAE0802	20	-	10-12
ZZ BASE B-C-SC 5.0 (Screw)	BSAE0803	20	-	10-12
ZZ BASE B-C-SC 6.0 (Screw)	BSAE0804	20	-	10-12
ZZ BASE B-C-SC-G 3.4 (Screw Gold)	BSCE0801	20	-	10-12
ZZ BASE B-C-SC-G 4.1 (Screw Gold)	BSCE0802	20	-	10-12
ZZ BASE B-C-SC-G 5.0 (Screw Gold)	BSCE0803	20	-	10-12
ZZ BASE B-C-SC-G 6.0 (Screw Gold)	BSCE0804	20	-	10-12
ZZ BASE B-C-TSC 3.4 (Thick screw)	BSAE0811	20	-	10-12
ZZ BASE B-C-TSC 4.1 (Thick screw)	BSAE0812	20	-	10-12
ZZ BASE B-C-TSC 5.0 (Thick screw)	BSAE0813	20	-	10-12
ZZ BASE B-C-TSC 6.0 (Thick screw)	BSAE0814	20	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

BIOMET 3I™ OSSEOTITE® CERTAIN®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
ZZ BASE B-C-ZSC 3.4 (Zirconia screw)	BSAE0831	20	-	10-12
ZZ BASE B-C-ZSC 4.1 (Zirconia screw)	BSAE0832	20	-	10-12
ZZ BASE B-C-ZSC 5.0 (Zirconia screw)	BSAE0833	20	-	10-12
ZZ BASE B-C-ZSC 6.0 (Zirconia screw)	BSAE0834	20	-	10-12

BIOMET 3I™ OSSEOTITE® EXTERNAL HEX

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE B-O-ABCH-ID 3.4 (Raw-Abutment® customizable hex ID and screw)	BSAF0991	35	-	10-12
SET ZZ BASE B-O-ABCH-ID 4.1 (Raw-Abutment® customizable hex ID and screw)	BSAF0992	35	-	10-12
SET ZZ BASE B-O-ABCH-ID 5.0 (Raw-Abutment® customizable hex ID and screw)	BSAF0993	35	-	10-12
SET ZZ BASE B-O-ABCH-ID 6.0 (Raw-Abutment® customizable hex ID and screw)	BSAF0994	35	-	10-12
SET ZZ BASE B-O-CCTB 3.4 (Conical cemented titan base non hex and screw)	BSAD0931	35	-	10-12
SET ZZ BASE B-O-CCTB 4.1 (Conical cemented titan base non hex and screw)	BSAD0932	35	-	10-12
SET ZZ BASE B-O-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD0933	35	-	10-12
SET ZZ BASE B-O-CCTB 6.0 (Conical cemented titan base non hex and screw)	BSAD0934	35	-	10-12
SET ZZ BASE B-O-CCTB-L10 3.4 (Conical cemented titan base non hex L10 and screw)	BSAD0931L10	35	-	10-12
SET ZZ BASE B-O-CCTB-L10 4.1 (Conical cemented titan base non hex L10 and screw)	BSAD0932L10	35	-	10-12
SET ZZ BASE B-O-CCTB-L10 5.0 (Conical cemented titan base non hex L10 and screw)	BSAD0933L10	35	-	10-12
SET ZZ BASE B-O-CCTB-L10 6.0 (Conical cemented titan base non hex L10 and screw)	BSAD0934L10	35	-	10-12
SET ZZ BASE B-O-CCTB-L20 3.4 (Conical cemented titan base non hex L20 and screw)	BSAD0931L20	35	-	10-12
SET ZZ BASE B-O-CCTB-L20 4.1 (Conical cemented titan base non hex L20 and screw)	BSAD0932L20	35	-	10-12
SET ZZ BASE B-O-CCTB-L20 5.0 (Conical cemented titan base non hex L20 and screw)	BSAD0933L20	35	-	10-12
SET ZZ BASE B-O-CCTB-L20 6.0 (Conical cemented titan base non hex L20 and screw)	BSAD0934L20	35	-	10-12
SET ZZ BASE B-O-CCTB-L30 3.4 (Conical cemented titan base non hex L30 and screw)	BSAD0931L30	35	-	10-12
SET ZZ BASE B-O-CCTB-L30 4.1 (Conical cemented titan base non hex L30 and screw)	BSAD0932L30	35	-	10-12
SET ZZ BASE B-O-CCTB-L30 5.0 (Conical cemented titan base non hex L30 and screw)	BSAD0933L30	35	-	10-12
SET ZZ BASE B-O-CCTB-L30 6.0 (Conical cemented titan base non hex L30 and screw)	BSAD0934L30	35	-	10-12
SET ZZ BASE B-O-CCTB-L40 3.4 (Conical cemented titan base non hex L40 and screw)	BSAD0931L40	35	-	10-12
SET ZZ BASE B-O-CCTB-L40 4.1 (Conical cemented titan base non hex L40 and screw)	BSAD0932L40	35	-	10-12
SET ZZ BASE B-O-CCTB-L40 5.0 (Conical cemented titan base non hex L40 and screw)	BSAD0933L40	35	-	10-12
SET ZZ BASE B-O-CCTB-L40 6.0 (Conical cemented titan base non hex L40 and screw)	BSAD0934L40	35	-	10-12
SET ZZ BASE B-O-NTB 3.4 (Narrow titan base non hex and screw)	BSAP0931	35	-	10-12
SET ZZ BASE B-O-NTB 4.1 (Narrow titan base non hex and screw)	BSAP0932	35	-	10-12
SET ZZ BASE B-O-NTBH 3.4 (Narrow titan base hex and screw)	BSAP0981	35	-	10-12
SET ZZ BASE B-O-NTBH 4.1 (Narrow titan base hex and screw)	BSAP0982	35	-	10-12
SET ZZ BASE B-O-PCTBH 3.4 (Parallel cemented titan base hex and screw)	BSAB0981	35	-	10-12
SET ZZ BASE B-O-PCTBH 4.1 (Parallel cemented titan base hex and screw)	BSAB0982	35	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

BIOMET 3I™ OSSEOTITE® EXTERNAL HEX

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE B-O-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB0983	35	-	10-12
SET ZZ BASE B-O-PCTBH 6.0 (Parallel cemented titan base hex and screw)	BSAB0984	35	-	10-12
SET ZZ BASE B-O-PCTBH-L10 3.4 (Parallel cemented titan base hex L10 and screw)	BSAB0991L10	35	-	10-12
SET ZZ BASE B-O-PCTBH-L10 4.1 (Parallel cemented titan base hex L10 and screw)	BSAB0992L10	35	-	10-12
SET ZZ BASE B-O-PCTBH-L10 5.0 (Parallel cemented titan base hex L10 and screw)	BSAB0993L10	35	-	10-12
SET ZZ BASE B-O-PCTBH-L10 6.0 (Parallel cemented titan base hex L10 and screw)	BSAB0994L10	35	-	10-12
SET ZZ BASE B-O-PCTBH-L20 3.4 (Parallel cemented titan base hex L20 and screw)	BSAB0991L20	35	-	10-12
SET ZZ BASE B-O-PCTBH-L20 4.1 (Parallel cemented titan base hex L20 and screw)	BSAB0992L20	35	-	10-12
SET ZZ BASE B-O-PCTBH-L20 5.0 (Parallel cemented titan base hex L20 and screw)	BSAB0993L20	35	-	10-12
SET ZZ BASE B-O-PCTBH-L20 6.0 (Parallel cemented titan base hex L20 and screw)	BSAB0994L20	35	-	10-12
SET ZZ BASE B-O-PCTBH-L30 3.4 (Parallel cemented titan base hex L30 and screw)	BSAB0991L30	35	-	10-12
SET ZZ BASE B-O-PCTBH-L30 4.1 (Parallel cemented titan base hex L30 and screw)	BSAB0992L30	35	-	10-12
SET ZZ BASE B-O-PCTBH-L30 5.0 (Parallel cemented titan base hex L30 and screw)	BSAB0993L30	35	-	10-12
SET ZZ BASE B-O-PCTBH-L30 6.0 (Parallel cemented titan base hex L30 and screw)	BSAB0994L30	35	-	10-12
SET ZZ BASE B-O-PCTBH-L40 3.4 (Parallel cemented titan base hex L40 and screw)	BSAB0991L40	35	-	10-12
SET ZZ BASE B-O-PCTBH-L40 4.1 (Parallel cemented titan base hex L40 and screw)	BSAB0992L40	35	-	10-12
SET ZZ BASE B-O-PCTBH-L40 5.0 (Parallel cemented titan base hex L40 and screw)	BSAB0993L40	35	-	10-12
SET ZZ BASE B-O-PCTBH-L40 6.0 (Parallel cemented titan base hex L40 and screw)	BSAB0994L40	35	-	10-12
SET ZZ BASE B-O-SM 3.4 (Scanmarker and screw)	ZBAD0931	-	-	10-12
SET ZZ BASE B-O-SM 4.1 (Scanmarker and screw)	ZBAD0932	-	-	10-12
SET ZZ BASE B-O-SM 5.0 (Scanmarker and screw)	ZBAD0933	-	-	10-12
SET ZZ BASE B-O-SM 6.0 (Scanmarker and screw)	ZBAD0934	-	-	10-12
SET ZZ BASE B-O-WSM 3.4 (White Scanmarker and screw)	ZBCD0931	-	-	10-12
SET ZZ BASE B-O-WSM 4.1 (White Scanmarker and screw)	ZBCD0932	-	-	10-12
SET ZZ BASE B-O-WSM 5.0 (White Scanmarker and screw)	ZBCD0933	-	-	10-12
SET ZZ BASE B-O-WSM 6.0 (White Scanmarker and screw)	ZBCD0934	-	-	10-12
ZZ BASE B-O-MI 3.4 (Model implant)	MIAG0901	-	-	10-12
ZZ BASE B-O-MI 4.1 (Model implant)	MIAG0902	-	-	10-12
ZZ BASE B-O-MI 5.0 (Model implant)	MIAG0903	-	-	10-12
ZZ BASE B-O-MI 6.0 (Model implant)	MIAG0904	-	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

BIOMET 3I™ OSSEOTITE® EXTERNAL HEX

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
ZZ BASE B-O-SC 3.4 (Screw)	BSAE0901	35	-	10-12
ZZ BASE B-O-SC 4.1 (Screw)	BSAE0902	35	-	10-12
ZZ BASE B-O-SC 5.0 (Screw)	BSAE0903	35	-	10-12
ZZ BASE B-O-SC 6.0 (Screw)	BSAE0904	35	-	10-12
ZZ BASE B-O-ZSC 3.4 (Zirconia screw)	BSAE0931	35	-	10-12
ZZ BASE B-O-ZSC 4.1 (Zirconia screw)	BSAE0932	35	-	10-12
ZZ BASE B-O-ZSC 5.0 (Zirconia screw)	BSAE0933	35	-	10-12
ZZ BASE B-O-ZSC 6.0 (Zirconia screw)	BSAE0934	35	-	10-12

BIOTECH KONTACT®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE BI-KO-ABCH 3.0 (Raw-Abutment® customizable hex and screw)	BSAF6281	20	-	10-12
SET ZZ BASE BI-KO-ABCH 3.6/4.2/4.8/5.4 (Raw-Abutment® customizable hex and screw)	BSAF6282	20	-	10-12
SET ZZ BASE BI-KO-CCTB 3.0 (Conical cemented titan base non hex and screw)	BSAD6231	20	-	10-12
SET ZZ BASE BI-KO-CCTB 3.6/4.2/4.8/5.4 (Conical cemented titan base non hex and screw)	BSAD6232	20	-	10-12
SET ZZ BASE BI-KO-CCTB 4.0 (Conical cemented titan base non hex and screw)	BSAD6233	20	-	10-12
SET ZZ BASE BI-KO-CCTB 4.9 (Conical cemented titan base non hex and screw)	BSAD6234	20	-	10-12
SET ZZ BASE BI-KO-PCTBH 3.0 (Parallel cemented titan base hex and screw)	BSAB6281	20	-	10-12
SET ZZ BASE BI-KO-PCTBH 3.6/4.2/4.8/5.4 (Parallel cemented titan base hex and screw)	BSAB6282	20	-	10-12
SET ZZ BASE BI-KO-SM 3.0 (Scanmarker and screw)	ZBAD6231	-	-	10-12
SET ZZ BASE BI-KO-SM 3.6/4.2/4.8/5.4 (Scanmarker and screw)	ZBAD6232	-	-	10-12
SET ZZ BASE BI-KO-SM 4.0 (Scanmarker and screw)	ZBAD6233	-	-	10-12
SET ZZ BASE BI-KO-SM 4.9 (Scanmarker and screw)	ZBAD6234	-	-	10-12
SET ZZ BASE BI-KO-WSM 3.0 (White Scanmarker and screw)	ZBCD6231	-	-	10-12
SET ZZ BASE BI-KO-WSM 3.6/4.2/4.8/5.4 (White Scanmarker and screw)	ZBCD6232	-	-	10-12
SET ZZ BASE BI-KO-WSM 4.0 (White Scanmarker and screw)	ZBCD6233	-	-	10-12
SET ZZ BASE BI-KO-WSM 4.9 (White Scanmarker and screw)	ZBCD6234	-	-	10-12
ZZ BASE BI-KO-SC 3.0 (Screw)	BSAE6201	20	-	10-12
ZZ BASE BI-KO-SC 3.6/4.2/4.8/5.4 (Screw)	BSAE6202	20	-	10-12
ZZ BASE BI-KO-SC 4.0 (Screw)	BSAE6203	20	-	10-12
ZZ BASE BI-KO-SC 4.9 (Screw)	BSAE6204	20	-	10-12
ZZ BASE BI-KO-ZSC 3.0 (Zirconia screw)	BSAE6231	20	-	10-12
ZZ BASE BI-KO-ZSC 3.6/4.2/4.8/5.4 (Zirconia screw)	BSAE6232	20	-	10-12
ZZ BASE BI-KO-ZSC 4.0 (Zirconia screw)	BSAE6233	20	-	10-12
ZZ BASE BI-KO-ZSC 4.9 (Zirconia screw)	BSAE6234	20	-	10-12

BRAINBASE CORPORATION® MYTIS ARROW IMPLAN

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE BB-MA-ABCH B (Raw-Abutment® customizable hex and screw)	BSAF9481	30	-	10-12
SET ZZ BASE BB-MA-ABCH C (Raw-Abutment® customizable hex and screw)	BSAF9482	30	-	10-12
SET ZZ BASE BB-MA-ABCH E (Raw-Abutment® customizable hex and screw)	BSAF9483	30	-	10-12
SET ZZ BASE BB-MA-ABCH EW/BW (Raw-Abutment® customizable hex and screw)	BSAF9484	30	-	10-12
SET ZZ BASE BB-MA-CCTB B (Conical cemented titan base non hex and screw)	BSAD9431	30	-	10-12
SET ZZ BASE BB-MA-CCTB C (Conical cemented titan base non hex and screw)	BSAD9432	30	-	10-12
SET ZZ BASE BB-MA-CCTB E (Conical cemented titan base non hex and screw)	BSAD9433	30	-	10-12
SET ZZ BASE BB-MA-CCTB EW/BW (Conical cemented titan base non hex and screw)	BSAD9434	30	-	10-12
SET ZZ BASE BB-MA-PCTBH B (Parallel cemented titan base hex and screw)	BSAB9481	30	-	10-12
SET ZZ BASE BB-MA-PCTBH C (Parallel cemented titan base hex and screw)	BSAB9482	30	-	10-12
SET ZZ BASE BB-MA-PCTBH E (Parallel cemented titan base hex and screw)	BSAB9483	30	-	10-12
SET ZZ BASE BB-MA-PCTBH EW/BW (Parallel cemented titan base hex and screw)	BSAB9484	30	-	10-12
SET ZZ BASE BB-MA-SM B (Scanmarker and screw)	ZBAD9431	30	-	10-12
SET ZZ BASE BB-MA-SM C (Scanmarker and screw)	ZBAD9432	30	-	10-12
SET ZZ BASE BB-MA-SM E (Scanmarker and screw)	ZBAD9433	30	-	10-12
SET ZZ BASE BB-MA-SM EW/BW (Scanmarker and screw)	ZBAD9434	30	-	10-12
SET ZZ BASE BB-MA-WSM B (White Scanmarker and screw)	ZBCD9431	30	-	10-12
SET ZZ BASE BB-MA-WSM C (White Scanmarker and screw)	ZBCD9432	30	-	10-12
SET ZZ BASE BB-MA-WSM E (White Scanmarker and screw)	ZBCD9433	30	-	10-12
SET ZZ BASE BB-MA-WSM EW/BW (White Scanmarker and screw)	ZBCD9434	30	-	10-12
ZZ BASE BB-MA-SC B (Screw)	BSAE9401	30	-	10-12
ZZ BASE BB-MA-SC C (Screw)	BSAE9402	30	-	10-12
ZZ BASE BB-MA-SC E (Screw)	BSAE9403	30	-	10-12
ZZ BASE BB-MA-SC EW/BW (Screw)	BSAE9404	30	-	10-12
ZZ BASE BB-MA-ZSC B (Zirconia Screw)	BSAE9431	30	-	10-12
ZZ BASE BB-MA-ZSC C (Zirconia Screw)	BSAE9432	30	-	10-12
ZZ BASE BB-MA-ZSC E (Zirconia Screw)	BSAE9433	30	-	10-12
ZZ BASE BB-MA-ZSC EW/BW (Zirconia Screw)	BSAE9434	30	-	10-12

BREDENT SKY CLASSIC/BLUESKY

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE SY-CB-ABCH-ID 3.5 (Raw-Abutment® customizable hex ID and screw)	BSAF2292	25	-	10-12
SET ZZ BASE SY-CB-ABCH-ID 4.0/4.5/5.5 (Raw-Abutment® customizable hex ID and screw)	BSAF2291	25	-	10-12
SET ZZ BASE SY-CB-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD2232	25	-	10-12
SET ZZ BASE SY-CB-CCTB 4.0/4.5/5.5 (Conical cemented titan base non hex and screw)	BSAD2231	25	-	10-12
SET ZZ BASE SY-CB-NTB 3.5 (Narrow titan base non hex and screw)	BSAP2232	25	-	10-12
SET ZZ BASE SY-CB-NTBH 3.5 (Narrow titan base hex and screw)	BSAP2282	25	-	10-12
SET ZZ BASE SY-CB-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB2282	25	-	10-12
SET ZZ BASE SY-CB-PCTBH 4.0/4.5/5.5 (Parallel cemented titan base hex and screw)	BSAB2281	25	-	10-12
SET ZZ BASE SY-CB-SM 3.5 (Scanmarker and screw)	ZBAD2232	-	-	10-12
SET ZZ BASE SY-CB-SM 4.0/4.5/5.5 (Scanmarker and screw)	ZBAD2231	-	-	10-12
SET ZZ BASE SY-CB-WSM 3.5 (White Scanmarker and screw)	ZBCD2232	-	-	10-12
SET ZZ BASE SY-CB-WSM 4.0/4.5/5.5 (White Scanmarker and screw)	ZBCD2231	-	-	10-12
SET ZZ BASE SY-CB-ZZMUA-L15 3.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH2232L15	25	25	10-12
SET ZZ BASE SY-CB-ZZMUA-L15 4.0/4.5/5.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH2231L15	25	25	10-12
SET ZZ BASE SY-CB-ZZMUA-L25 3.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH2232L25	25	25	10-12
SET ZZ BASE SY-CB-ZZMUA-L25 4.0/4.5/5.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH2231L25	25	25	10-12
SET ZZ BASE SY-CB-ZZMUA-L35 3.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH2232L35	25	25	10-12
SET ZZ BASE SY-CB-ZZMUA-L35 4.0/4.5/5.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH2231L35	25	25	10-12
SET ZZ BASE SY-CB-ZZMUA-L45 3.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH2232L45	25	25	10-12
SET ZZ BASE SY-CB-ZZMUA-L45 4.0/4.5/5.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH2231L45	25	25	10-12
ZZ BASE SY-CB-SC 3.5 (Screw)	BSAE2202	25	-	10-12
ZZ BASE SY-CB-SC 4.0/4.5/5.5 (Screw)	BSAE2201	25	-	10-12
ZZ BASE SY-CB-ZSC 3.5 (Zirconia screw)	BSAE2232	25	-	10-12
ZZ BASE SY-CB-ZSC 4.0/4.5/5.5 (Zirconia screw)	BSAE2231	25	-	10-12
ZZ BASE SY-CB-ZSC-G 4.0/4.5/5.5 (Zirconia Screw Gold)	BSCE2231	25	-	10-12

BREDENT SKY FAST & FIXED

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE SY-FF-CCTB (Conical cemented titan base non hex and screw)	BSAD2331	20	-	10-12
SET ZZ BASE SY-FF-SM (Scanmarker and screw)	ZBAD2331	-	-	10-12
ZZ BASE SY-FF-SC (Screw)	BSAE2301	20	-	10-12
ZZ BASE SY-FF-ZSC (Zirconia screw)	BSAE2331	20	-	10-12

BREDENT SKY UNI.CONE

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE SY-UC-CCTB (Conical cemented titan base non hex and screw)	BSAD8931	18	-	10-12
SET ZZ BASE SY-UC-SM (Scanmarker and screw)	ZBAD8931	-	-	10-12
SET ZZ BASE SY-UC-WSM (White Scanmarker and screw)	ZBCD8931	-	-	10-12
ZZ BASE SY-UC-MI (Model implant)	MIAG8901	-	-	10-12
ZZ BASE SY-UC-SC (Screw)	BSAE8901	18	-	10-12
ZZ BASE SY-UC-ZSC (Zirconia screw)	BSAE8931	18	-	10-12

BTI® CONICAL SPACER

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE BT-MCS-CCTB 4.1/5.5 (Conical cemented titan base non hex and screw)	BSAD3333	10	-	10-12
SET ZZ BASE BT-MCS-SM 4.1/5.5 (Scanmarker and screw)	ZBAD3333	-	-	10-12
SET ZZ BASE BT-MCS-WSM 4.1/5.5 (White Scanmarker and screw)	ZBCD3333	-	-	10-12
ZZ BASE BT-MCS-SC 4.1/5.5 (Screw)	BSAE3303	10	-	10-12
ZZ BASE BT-MCS-ZSC 4.1/5.5 (Zirconia screw)	BSAE3333	10	-	10-12

BTI® EXTERNA®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE BT-E-ABCH-ID 4.1 (Raw-Abutment® customizable hex ID and screw)	BSAF3291	20	-	10-12
SET ZZ BASE BT-E-ABCH-ID 5.5 (Raw-Abutment® customizable hex ID and screw)	BSAF3292	20	-	10-12
SET ZZ BASE BT-E-CCTB 4.1 (Conical cemented titan base non hex and screw)	BSAD3231	20	-	10-12
SET ZZ BASE BT-E-CCTB 5.5 (Conical cemented titan base non hex and screw)	BSAD3232	20	-	10-12
SET ZZ BASE BT-E-PCTBH 4.1 (Parallel cemented titan base hex and screw)	BSAB3281	20	-	10-12
SET ZZ BASE BT-E-PCTBH 5.5 (Parallel cemented titan base hex and screw)	BSAB3282	20	-	10-12
SET ZZ BASE BT-E-SM 4.1 (Scanmarker and screw)	ZBAD3231	-	-	10-12
SET ZZ BASE BT-E-SM 5.5 (Scanmarker and screw)	ZBAD3232	-	-	10-12
SET ZZ BASE BT-E-WSM 4.1 (White Scanmarker and screw)	ZBCD3231	-	-	10-12
SET ZZ BASE BT-E-WSM 5.5 (White Scanmarker and screw)	ZBCD3232	-	-	10-12
ZZ BASE BT-E-SC 4.1 (Screw)	BSAE3201	20	-	10-12
ZZ BASE BT-E-SC 5.5 (Screw)	BSAE3202	20	-	10-12
ZZ BASE BT-E-SC-G 4.1 (Screw Gold)	BSCE3201	20	-	10-12
ZZ BASE BT-E-SC-G 5.5 (Screw Gold)	BSCE3202	20	-	10-12
ZZ BASE BT-E-ZSC 4.1 (Zirconia screw)	BSAE3231	20	-	10-12
ZZ BASE BT-E-ZSC 5.5 (Zirconia screw)	BSAE3232	20	-	10-12
ZZ BASE BT-E-ZSC-G 4.1 (Zirconia Screw Gold)	BSCE3231	20	-	10-12
ZZ BASE BT-E-ZSC-G 5.5 (Zirconia Screw Gold)	BSCE3232	20	-	10-12

BTI® INTERNA®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE BT-I-ABCH-ID 4.1 (Raw-Abutment® customizable hex ID and screw)	BSAF3691	20	-	10-12
SET ZZ BASE BT-I-ABCH-ID 5.5 (Raw-Abutment® customizable hex ID and screw)	BSAF3692	20	-	10-12
SET ZZ BASE BT-I-CCTB 4.1 (Conical cemented titan base non hex and screw)	BSAD3631	20	-	10-12
SET ZZ BASE BT-I-CCTB 5.5 (Conical cemented titan base non hex and screw)	BSAD3632	20	-	10-12
SET ZZ BASE BT-I-NTB 4.1 (Narrow titan base non hex and screw)	BSAP3631	20	-	10-12
SET ZZ BASE BT-I-NTBH 4.1 (Narrow titan base hex and screw)	BSAP3681	20	-	10-12
SET ZZ BASE BT-I-PCTBH 4.1 (Parallel cemented titan base hex and screw)	BSAB3681	20	-	10-12
SET ZZ BASE BT-I-PCTBH 5.5 (Parallel cemented titan base hex and screw)	BSAB3682	20	-	10-12
SET ZZ BASE BT-I-SM 4.1 (Scanmarker and screw)	ZBAD3631	-	-	10-12
SET ZZ BASE BT-I-SM 5.5 (Scanmarker and screw)	ZBAD3632	-	-	10-12
SET ZZ BASE BT-I-WSM 4.1 (White Scanmarker and screw)	ZBCD3631	-	-	10-12
SET ZZ BASE BT-I-WSM 5.5 (White Scanmarker and screw)	ZBCD3632	-	-	10-12
ZZ BASE BT-I-SC 4.1 (Screw)	BSAE3601	20	-	10-12
ZZ BASE BT-I-SC 5.5 (Screw)	BSAE3602	20	-	10-12
ZZ BASE BT-I-ZSC 4.1 (Zirconia screw)	BSAE3631	20	-	10-12
ZZ BASE BT-I-ZSC 5.5 (Zirconia screw)	BSAE3632	20	-	10-12

BTI® MULTI-IM

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE BT-MIM-CCTB 3.5/4.1 (Conical cemented titan base non hex and screw)	BSAD3331	10	-	10-12
SET ZZ BASE BT-MIM-CCTB 5.5 (Conical cemented titan base non hex and screw)	BSAD3332	10	-	10-12
SET ZZ BASE BT-MIM-SM 3.5/4.1 (Scanmarker and screw)	ZBAD3331	-	-	10-12
SET ZZ BASE BT-MIM-SM 5.5 (Scanmarker and screw)	ZBAD3332	-	-	10-12
SET ZZ BASE BT-MIM-WSM 3.5/4.1 (White Scanmarker and screw)	ZBCD3331	-	-	10-12
SET ZZ BASE BT-MIM-WSM 5.5 (White Scanmarker and screw)	ZBCD3332	-	-	10-12
ZZ BASE BT-MIM-SC 3.5/4.1 (Screw)	BSAE3301	10	-	10-12
ZZ BASE BT-MIM-SC 5.5 (Screw)	BSAE3302	10	-	10-12
ZZ BASE BT-MIM-ZSC 3.5/4.1 (Zirconia screw)	BSAE3331	10	-	10-12
ZZ BASE BT-MIM-ZSC 5.5 (Zirconia screw)	BSAE3332	10	-	10-12

BTI® MULTI-IM ANGLED

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE BT-MIMA-CCTB 3.5/4.1 (Conical cemented titan base non hex and screw)	BSAD3334	10	-	10-12
SET ZZ BASE BT-MIMA-SM 3.5/4.1 (Scanmarker and screw)	ZBAD3334	-	-	10-12
SET ZZ BASE BT-MIMA-WSM 3.5/4.1 (White Scanmarker and screw)	ZBCD3334	-	-	10-12
ZZ BASE BT-MIMA-SC 3.5/4.1 (Screw)	BSAE3304	10	-	10-12
ZZ BASE BT-MIMA-ZSC 3.5/4.1 (Zirconia screw)	BSAE3334	10	-	10-12

BTI® TINY®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE BT-T-ABCH-ID 3.5 (Raw-Abutment® customizable hex ID and screw)	BSAF3191	10	-	10-12
SET ZZ BASE BT-T-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD3131	10	-	10-12
SET ZZ BASE BT-T-NTB 3.5 (Narrow titan base non hex and screw)	BSAP3131	10	-	10-12
SET ZZ BASE BT-T-NTBH 3.5 (Narrow titan base hex and screw)	BSAP3181	10	-	10-12
SET ZZ BASE BT-T-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB3181	10	-	10-12
SET ZZ BASE BT-T-SM 3.5 (Scanmarker and screw)	ZBAD3131	-	-	10-12
SET ZZ BASE BT-T-WSM 3.5 (White Scanmarker and screw)	ZBCD3131	-	-	10-12
ZZ BASE BT-T-SC 3.5 (Screw)	BSAE3101	10	-	10-12
ZZ BASE BT-T-SC/ZSC-G 3.5 (Screw Gold)	BSCE3101	10	-	10-12
ZZ BASE BT-T-ZSC 3.5 (Zirconia screw)	BSAE3131	10	-	10-12

CAMLOG® BAR ABUTMENTS (COMFOUR®)

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE C-BA-CCTB 3.3/3.8/4.3 (Conical cemented titan base non hex and screw)	BSAD2131	20	-	10-12
SET ZZ BASE C-BA-CCTB 5.0/6.0 (Conical cemented titan base non hex and screw)	BSAD2132	35	-	10-12
SET ZZ BASE C-BA-SM 3.3/3.8/4.3 (Scanmarker and screw)	ZBAD2131	-	-	10-12
SET ZZ BASE C-BA-SM 5.0/6.0 (Scanmarker and screw)	ZBAD2132	-	-	10-12
ZZ BASE C-BA-SC 3.3/3.8/4.3 (Screw)	BSAE2101	20	-	10-12
ZZ BASE C-BA-SC 5.0/6.0 (Screw)	BSAE2102	35	-	10-12
ZZ BASE C-BA-SC-G 3.3/3.8/4.3 (Screw Gold)	BSCE2101	20	-	10-12
ZZ BASE C-BA-SC-G 5.0/6.0 (Screw Gold)	BSCE2102	35	-	10-12
ZZ BASE C-BA-ZSC 3.3/3.8/4.3 (Zirconia screw)	BSAE2131	20	-	10-12
ZZ BASE C-BA-ZSC 5.0/6.0 (Zirconia screw)	BSAE2132	35	-	10-12

CAMLOG® CONELOG®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE C-CL-ABCH 3.3 (Raw-Abutment® customizable hex and screw)	BSAF4781	20	-	10-12
SET ZZ BASE C-CL-ABCH 3.8 (Raw-Abutment® customizable hex and screw)	BSAF4782	20	-	10-12
SET ZZ BASE C-CL-ABCH 4.3 (Raw-Abutment® customizable hex and screw)	BSAF4783	20	-	10-12
SET ZZ BASE C-CL-ABCH 5.0 (Raw-Abutment® customizable hex and screw)	BSAF4784	20	-	10-12
SET ZZ BASE C-CL-CCTB 3.3 (Conical cemented titan base non hex and screw)	BSAD4731	20	-	10-12
SET ZZ BASE C-CL-CCTB 3.8 (Conical cemented titan base non hex and screw)	BSAD4732	20	-	10-12
SET ZZ BASE C-CL-CCTB 4.3 (Conical cemented titan base non hex and screw)	BSAD4733	20	-	10-12
SET ZZ BASE C-CL-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD4734	20	-	10-12
SET ZZ BASE C-CL-NTB 3.3 (Narrow titan base non hex and screw)	BSAP4731	20	-	10-12
SET ZZ BASE C-CL-NTB 3.8 (Narrow titan base non hex and screw)	BSAP4732	20	-	10-12
SET ZZ BASE C-CL-NTB 4.3 (Narrow titan base non hex and screw)	BSAP4733	20	-	10-12
SET ZZ BASE C-CL-NTBH 3.3 (Narrow titan base hex and screw)	BSAP4781	20	-	10-12
SET ZZ BASE C-CL-NTBH 3.8 (Narrow titan base hex and screw)	BSAP4782	20	-	10-12
SET ZZ BASE C-CL-NTBH 4.3 (Narrow titan base hex and screw)	BSAP4783	20	-	10-12
SET ZZ BASE C-CL-PCTBH 3.3 (Parallel cemented titan base hex and screw)	BSAB4781	20	-	10-12
SET ZZ BASE C-CL-PCTBH 3.8 (Parallel cemented titan base hex and screw)	BSAB4782	20	-	10-12
SET ZZ BASE C-CL-PCTBH 4.3 (Parallel cemented titan base hex and screw)	BSAB4783	20	-	10-12
SET ZZ BASE C-CL-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB4784	20	-	10-12
SET ZZ BASE C-CL-SM 3.3 (Scanmarker and screw)	ZBAD4731	-	-	10-12
SET ZZ BASE C-CL-SM 3.8 (Scanmarker and screw)	ZBAD4732	-	-	10-12
SET ZZ BASE C-CL-SM 4.3 (Scanmarker and screw)	ZBAD4733	-	-	10-12
SET ZZ BASE C-CL-SM 5.0 (Scanmarker and screw)	ZBAD4734	-	-	10-12
SET ZZ BASE C-CL-WSM 3.3 (White Scanmarker and screw)	ZBCD4731	-	-	10-12
SET ZZ BASE C-CL-WSM 3.8 (White Scanmarker and screw)	ZBCD4732	-	-	10-12
SET ZZ BASE C-CL-WSM 4.3 (White Scanmarker and screw)	ZBCD4733	-	-	10-12
SET ZZ BASE C-CL-WSM 5.0 (White Scanmarker and screw)	ZBCD4734	-	-	10-12
ZZ BASE C-CL-MI 3.3 (Model implant)	MIAG4701	-	-	10-12
ZZ BASE C-CL-MI 3.8 (Model implant)	MIAG4702	-	-	10-12
ZZ BASE C-CL-MI 4.3 (Model implant)	MIAG4703	-	-	10-12
ZZ BASE C-CL-MI 5.0 (Model implant)	MIAG4704	-	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

CAMLOG® CONELOG®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
ZZ BASE C-CL-SC 3.3 (Screw)	BSAE4701	20	-	10-12
ZZ BASE C-CL-SC 3.8 (Screw)	BSAE4702	20	-	10-12
ZZ BASE C-CL-SC 4.3 (Screw)	BSAE4703	20	-	10-12
ZZ BASE C-CL-SC 5.0 (Screw)	BSAE4704	20	-	10-12
ZZ BASE C-CL-ZSC 3.3 (Zirconia screw)	BSAE4731	20	-	10-12
ZZ BASE C-CL-ZSC 3.8 (Zirconia screw)	BSAE4732	20	-	10-12
ZZ BASE C-CL-ZSC 4.3 (Zirconia screw)	BSAE4733	20	-	10-12
ZZ BASE C-CL-ZSC 5.0 (Zirconia screw)	BSAE4734	20	-	10-12

CAMLOG® J-TYPE/K-TYPE

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE C-J-SM 3.3 (Scanmarker and screw)	ZBAD2031	-	-	10-12
SET ZZ BASE C-J-SM 3.8 (Scanmarker and screw)	ZBAD2032	-	-	10-12
SET ZZ BASE C-J-SM 4.3 (Scanmarker and screw)	ZBAD2033	-	-	10-12
SET ZZ BASE C-J-SM 5.0 (Scanmarker and screw)	ZBAD2034	-	-	10-12
SET ZZ BASE C-J-ABCH-D14-ID 3.3 (Raw-Abutment® customizable hex D14 ID and screw)	BSAF2091D14	20	-	10-12
SET ZZ BASE C-J-ABCH-D14-ID 3.8 (Raw-Abutment® customizable hex D14 ID and screw)	BSAF2092D14	20	-	10-12
SET ZZ BASE C-J-ABCH-D14-ID 4.3 (Raw-Abutment® customizable hex D14 ID and screw)	BSAF2093D14	20	-	10-12
SET ZZ BASE C-J-ABCH-D14-ID 5.0 (Raw-Abutment® customizable hex D14 ID and screw)	BSAF2094D14	30	-	10-12
SET ZZ BASE C-J-ABCH-D14-ID 6.0 (Raw-Abutment® customizable hex D14 ID and screw)	BSAF2095D14	30	-	10-12
SET ZZ BASE C-J-ABCH-ID 3.3 (Raw-Abutment® customizable hex ID and screw)	BSAF2091	20	-	10-12
SET ZZ BASE C-J-ABCH-ID 3.8 (Raw-Abutment® customizable hex ID and screw)	BSAF2092	20	-	10-12
SET ZZ BASE C-J-ABCH-ID 4.3 (Raw-Abutment® customizable hex ID and screw)	BSAF2093	20	-	10-12
SET ZZ BASE C-J-ABCH-ID 5.0 (Raw-Abutment® customizable hex ID and screw)	BSAF2094	30	-	10-12
SET ZZ BASE C-J-ABCH-ID 6.0 (Raw-Abutment® customizable hex ID and screw)	BSAF2095	30	-	10-12
SET ZZ BASE C-J-CCTB 3.3 (Conical cemented titan base non hex and screw)	BSAD2031	20	-	10-12
SET ZZ BASE C-J-CCTB 3.8 (Conical cemented titan base non hex and screw)	BSAD2032	20	-	10-12
SET ZZ BASE C-J-CCTB 4.3 (Conical cemented titan base non hex and screw)	BSAD2033	20	-	10-12
SET ZZ BASE C-J-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD2034	30	-	10-12
SET ZZ BASE C-J-CCTB 6.0 (Conical cemented titan base non hex and screw)	BSAD2035	30	-	10-12
SET ZZ BASE C-J-CCTB-L10 3.3 (Conical cemented titan base non hex L10 and screw)	BSAD2031L10	20	-	10-12
SET ZZ BASE C-J-CCTB-L10 3.8 (Conical cemented titan base non hex L10 and screw)	BSAD2032L10	20	-	10-12
SET ZZ BASE C-J-CCTB-L10 4.3 (Conical cemented titan base non hex L10 and screw)	BSAD2033L10	20	-	10-12
SET ZZ BASE C-J-CCTB-L10 5.0 (Conical cemented titan base non hex L10 and screw)	BSAD2034L10	30	-	10-12
SET ZZ BASE C-J-CCTB-L10 6.0 (Conical cemented titan base non hex L10 and screw)	BSAD2035L10	30	-	10-12
SET ZZ BASE C-J-CCTB-L20 3.3 (Conical cemented titan base non hex L20 and screw)	BSAD2031L20	20	-	10-12
SET ZZ BASE C-J-CCTB-L20 3.8 (Conical cemented titan base non hex L20 and screw)	BSAD2032L20	20	-	10-12
SET ZZ BASE C-J-CCTB-L20 4.3 (Conical cemented titan base non hex L20 and screw)	BSAD2033L20	20	-	10-12
SET ZZ BASE C-J-CCTB-L20 5.0 (Conical cemented titan base non hex L20 and screw)	BSAD2034L20	30	-	10-12
SET ZZ BASE C-J-CCTB-L20 6.0 (Conical cemented titan base non hex L20 and screw)	BSAD2035L20	30	-	10-12
SET ZZ BASE C-J-CCTB-L30 3.3 (Conical cemented titan base non hex L30 and screw)	BSAD2031L30	20	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

CAMLOG® J-TYPE/K-TYPE

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE C-J-CCTB-L30 3.8 (Conical cemented titan base non hex L30 and screw)	BSAD2032L30	20	-	10-12
SET ZZ BASE C-J-CCTB-L30 4.3 (Conical cemented titan base non hex L30 and screw)	BSAD2033L30	20	-	10-12
SET ZZ BASE C-J-CCTB-L30 5.0 (Conical cemented titan base non hex L30 and screw)	BSAD2034L30	30	-	10-12
SET ZZ BASE C-J-CCTB-L30 6.0 (Conical cemented titan base non hex L30 and screw)	BSAD2035L30	30	-	10-12
SET ZZ BASE C-J-NTB 3.3 (Narrow titan base non hex and screw)	BSAP2031	20	-	10-12
SET ZZ BASE C-J-NTB 3.8 (Narrow titan base non hex and screw)	BSAP2032	20	-	10-12
SET ZZ BASE C-J-NTBH 3.3 (Narrow titan base hex and screw)	BSAP2081	20	-	10-12
SET ZZ BASE C-J-NTBH 3.8 (Narrow titan base hex and screw)	BSAP2082	20	-	10-12
SET ZZ BASE C-J-PCTBH 3.3 (Parallel cemented titan base hex and screw)	BSAB2081	20	-	10-12
SET ZZ BASE C-J-PCTBH 3.8 (Parallel cemented titan base hex and screw)	BSAB2082	20	-	10-12
SET ZZ BASE C-J-PCTBH 4.3 (Parallel cemented titan base hex and screw)	BSAB2083	20	-	10-12
SET ZZ BASE C-J-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB2084	30	-	10-12
SET ZZ BASE C-J-PCTBH 6.0 (Parallel cemented titan base hex and screw)	BSAB2085	30	-	10-12
SET ZZ BASE C-J-PCTBH-L10 3.8 (Parallel cemented titan base hex L10 and screw)	BSAB2092L10	20	-	10-12
SET ZZ BASE C-J-PCTBH-L10 4.3 (Parallel cemented titan base hex L10 and screw)	BSAB2093L10	20	-	10-12
SET ZZ BASE C-J-PCTBH-L10 5.0 (Parallel cemented titan base hex L10 and screw)	BSAB2094L10	30	-	10-12
SET ZZ BASE C-J-PCTBH-L10 6.0 (Parallel cemented titan base hex L10 and screw)	BSAB2095L10	30	-	10-12
SET ZZ BASE C-J-PCTBH-L20 3.3 (Parallel cemented titan base hex L20 and screw)	BSAB2091L20	20	-	10-12
SET ZZ BASE C-J-PCTBH-L20 3.8 (Parallel cemented titan base hex L20 and screw)	BSAB2092L20	20	-	10-12
SET ZZ BASE C-J-PCTBH-L20 4.3 (Parallel cemented titan base hex L20 and screw)	BSAB2093L20	20	-	10-12
SET ZZ BASE C-J-PCTBH-L20 5.0 (Parallel cemented titan base hex L20 and screw)	BSAB2094L20	30	-	10-12
SET ZZ BASE C-J-PCTBH-L20 6.0 (Parallel cemented titan base hex L20 and screw)	BSAB2095L20	30	-	10-12
SET ZZ BASE C-J-PCTBH-L30 3.3 (Parallel cemented titan base hex L30 and screw)	BSAB2091L30	20	-	10-12
SET ZZ BASE C-J-PCTBH-L30 3.8 (Parallel cemented titan base hex L30 and screw)	BSAB2092L30	20	-	10-12
SET ZZ BASE C-J-PCTBH-L30 4.3 (Parallel cemented titan base hex L30 and screw)	BSAB2093L30	20	-	10-12
SET ZZ BASE C-J-PCTBH-L30 5.0 (Parallel cemented titan base hex L30 and screw)	BSAB2094L30	30	-	10-12
SET ZZ BASE C-J-PCTBH-L30 6.0 (Parallel cemented titan base hex L30 and screw)	BSAB2095L30	30	-	10-12
SET ZZ BASE C-J-SA 3.3 (Scanalog and screw)	ZBDE2031	15	-	10-12
SET ZZ BASE C-J-SA 3.8 (Scanalog and screw)	ZBDE2032	15	-	10-12
SET ZZ BASE C-J-SA 4.3 (Scanalog and screw)	ZBDE2033	15	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

CAMLOG® J-TYPE/K-TYPE

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE C-J-SA 5.0 (Scanalog and screw)	ZBDE2034	15	-	10-12
SET ZZ BASE C-J-SA 6.0 (Scanalog and screw)	ZBDE2035	15	-	10-12
SET ZZ BASE C-J-SM 6.0 (Scanmarker and screw)	ZBAD2035	-	-	10-12
SET ZZ BASE C-J-WSM 3.3 (White Scanmarker and screw)	ZBCD2031	-	-	10-12
SET ZZ BASE C-J-WSM 3.8 (White Scanmarker and screw)	ZBCD2032	-	-	10-12
SET ZZ BASE C-J-WSM 4.3 (White Scanmarker and screw)	ZBCD2033	-	-	10-12
SET ZZ BASE C-J-WSM 5.0 (White Scanmarker and screw)	ZBCD2034	-	-	10-12
SET ZZ BASE C-J-WSM 6.0 (White Scanmarker and screw)	ZBCD2035	-	-	10-12
SET ZZ BASE C-J-ZZMUA-L05 3.8 (Zirkonzahn Multi Unit Abutment non hex L05 and Screw)	BSAH2032L05	25	20	10-12
SET ZZ BASE C-J-ZZMUA-L05 4.3 (Zirkonzahn Multi Unit Abutment non hex L05 and Screw)	BSAH2033L05	25	20	10-12
SET ZZ BASE C-J-ZZMUA-L05 5.0 (Zirkonzahn Multi Unit Abutment non hex L05 and Screw)	BSAH2034L05	25	30	10-12
SET ZZ BASE C-J-ZZMUA-L05 6.0 (Zirkonzahn Multi Unit Abutment non hex L05 and Screw)	BSAH2035L05	25	30	10-12
SET ZZ BASE C-J-ZZMUA-L15 3.3 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH2031L15	25	20	10-12
SET ZZ BASE C-J-ZZMUA-L15 3.8 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH2032L15	25	20	10-12
SET ZZ BASE C-J-ZZMUA-L15 4.3 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH2033L15	25	20	10-12
SET ZZ BASE C-J-ZZMUA-L15 5.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH2034L15	25	30	10-12
SET ZZ BASE C-J-ZZMUA-L15 6.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH2035L15	25	30	10-12
SET ZZ BASE C-J-ZZMUA-L25 3.3 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH2031L25	25	20	10-12
SET ZZ BASE C-J-ZZMUA-L25 3.8 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH2032L25	25	20	10-12
SET ZZ BASE C-J-ZZMUA-L25 4.3 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH2033L25	25	20	10-12
SET ZZ BASE C-J-ZZMUA-L25 5.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH2034L25	25	30	10-12
SET ZZ BASE C-J-ZZMUA-L25 6.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH2035L25	25	30	10-12
SET ZZ BASE C-J-ZZMUA-L35 3.3 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH2031L35	25	20	10-12
SET ZZ BASE C-J-ZZMUA-L35 3.8 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH2032L35	25	20	10-12
SET ZZ BASE C-J-ZZMUA-L35 4.3 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH2033L35	25	20	10-12
SET ZZ BASE C-J-ZZMUA-L35 5.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH2034L35	25	30	10-12
SET ZZ BASE C-J-ZZMUA-L35 6.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH2035L35	25	30	10-12
SET ZZ BASE C-J-ZZMUA-L45 3.3 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH2031L45	25	20	10-12
SET ZZ BASE C-J-ZZMUA-L45 3.8 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH2032L45	25	20	10-12
SET ZZ BASE C-J-ZZMUA-L45 4.3 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH2033L45	25	20	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

CAMLOG® J-TYPE/K-TYPE

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE C-J-ZZMUA-L45 5.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH2034L45	25	30	10-12
SET ZZ BASE C-J-ZZMUA-L45 6.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH2035L45	25	30	10-12
ZZ BASE C-J-MI 3.3 (Model implant)	MIAG2001	-	-	10-12
ZZ BASE C-J-MI 3.8 (Model implant)	MIAG2002	-	-	10-12
ZZ BASE C-J-MI 4.3 (Model implant)	MIAG2003	-	-	10-12
ZZ BASE C-J-MI 5.0 (Model implant)	MIAG2004	-	-	10-12
ZZ BASE C-J-MI 6.0 (Model implant)	MIAG2005	-	-	10-12
ZZ BASE C-J-SC 3.3 (Screw)	BSAE2001	20	-	10-12
ZZ BASE C-J-SC 3.8 (Screw)	BSAE2002	20	-	10-12
ZZ BASE C-J-SC 4.3 (Screw)	BSAE2003	20	-	10-12
ZZ BASE C-J-SC 5.0 (Screw)	BSAE2004	30	-	10-12
ZZ BASE C-J-SC 6.0 (Screw)	BSAE2005	30	-	10-12
ZZ BASE C-J-SC-G 3.3 (Screw Gold)	BSCE2001	20	-	10-12
ZZ BASE C-J-SC-G 3.8 (Screw Gold)	BSCE2002	20	-	10-12
ZZ BASE C-J-SC-G 4.3 (Screw Gold)	BSCE2003	20	-	10-12
ZZ BASE C-J-SC-G 5.0 (Screw Gold)	BSCE2004	30	-	10-12
ZZ BASE C-J-SC-G 6.0 (Screw Gold)	BSCE2005	30	-	10-12
ZZ BASE C-J-ZSC 3.3 (Zirconia screw)	BSAE2031	20	-	10-12
ZZ BASE C-J-ZSC 3.8 (Zirconia screw)	BSAE2032	20	-	10-12
ZZ BASE C-J-ZSC 4.3 (Zirconia screw)	BSAE2033	20	-	10-12
ZZ BASE C-J-ZSC 5.0 (Zirconia screw)	BSAE2034	30	-	10-12
ZZ BASE C-J-ZSC 6.0 (Zirconia screw)	BSAE2035	30	-	10-12
ZZ BASE C-J-ZSC-G 3.3 (Zirconia Screw Gold)	BSCE2031	20	-	10-12
ZZ BASE C-J-ZSC-G 3.8 (Zirconia Screw Gold)	BSCE2032	20	-	10-12
ZZ BASE C-J-ZSC-G 4.3 (Zirconia Screw Gold)	BSCE2033	20	-	10-12
ZZ BASE C-J-ZSC-G 5.0 (Zirconia Screw Gold)	BSCE2034	30	-	10-12
ZZ BASE C-J-ZSC-G 6.0 (Zirconia Screw Gold)	BSCE2035	30	-	10-12

CAMLOG® VARIO SR

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE C-VSR-CCTB 3.8/4.3 (Conical cemented titan base non hex and screw)	BSAD5431	20	-	10-12
SET ZZ BASE C-VSR-CCTB 5.0/6.0 (Conical cemented titan base non hex and screw)	BSAD5432	20	-	10-12
SET ZZ BASE C-VSR-SM 3.8/4.3 (Scanmarker and screw)	ZBAD5431	-	-	10-12
SET ZZ BASE C-VSR-SM 5.0/6.0 (Scanmarker and screw)	ZBAD5432	-	-	10-12
ZZ BASE C-VSR-SC 3.8/4.3 (Screw)	BSAE5401	20	-	10-12
ZZ BASE C-VSR-SC 5.0/6.0 (Screw)	BSAE5402	20	-	10-12
ZZ BASE C-VSR-ZSC 3.8/4.3 (Zirconia screw)	BSAE5431	20	-	10-12
ZZ BASE C-VSR-ZSC 5.0/6.0 (Zirconia screw)	BSAE5432	20	-	10-12

CONMET HEX

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE CO-H-CCTB HEX 2.2 (Conical cemented titan base non hex and screw)	BSAD7831	20	-	10-12
SET ZZ BASE CO-H-CCTB HEX 2.7 (Conical cemented titan base non hex and screw)	BSAD7832	20	-	10-12
SET ZZ BASE CO-H-PCTBH HEX 2.2 (Parallel cemented titan base hex and screw)	BSAB7881	20	-	10-12
SET ZZ BASE CO-H-PCTBH HEX 2.7 (Parallel cemented titan base hex and screw)	BSAB7882	20	-	10-12
SET ZZ BASE CO-H-SM HEX 2.2 (Scanmarker and screw)	ZBAD7831	-	-	10-12
SET ZZ BASE CO-H-SM HEX 2.7 (Scanmarker and screw)	ZBAD7832	-	-	10-12
ZZ BASE CO-H-SC HEX 2.2 (Screw)	BSAE7801	20	-	10-12
ZZ BASE CO-H-SC HEX 2.7 (Screw)	BSAE7802	20	-	10-12
ZZ BASE CO-H-ZSC HEX 2.2 (Zirconia Screw)	BSAE7831	20	-	10-12
ZZ BASE CO-H-ZSC HEX 2.7 (Zirconia Screw)	BSAE7832	20	-	10-12

CONNECT®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE CN-MC-CCTB (Conical cemented titan base non hex and screw)	BSAD8831	20	-	10-12
SET ZZ BASE CN-MC-PCTBH (Parallel cemented titan base hex and screw)	BSAB8881	20	-	10-12
SET ZZ BASE CN-MC-SM (Scanmarker and screw)	ZBAD8831	-	-	10-12
SET ZZ BASE CN-ST-CCTB (Conical cemented titan base non hex and screw)	BSAD8832	35	-	10-12
SET ZZ BASE CN-ST-PCTBH (Parallel cemented titan base hex and screw)	BSAB8882	35	-	10-12
SET ZZ BASE CN-ST-SM (Scanmarker and screw)	ZBAD8832	-	-	10-12
ZZ BASE CN-MC-SC (Screw)	BSAE8801	20	-	10-12
ZZ BASE CN-MC-ZSC (Zirconia Screw)	BSAE8831	20	-	10-12
ZZ BASE CN-ST-SC (Screw)	BSAE8802	35	-	10-12
ZZ BASE CN-ST-ZSC (Zirconia Screw)	BSAE8832	35	-	10-12

COWELLMEDI INNO

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE CM-IN-ABCH 3.5-6.0 (Raw-Abutment® customizable hex and screw)	BSAF6981	30	-	10-12
SET ZZ BASE CM-IN-SM 3.5-6.0 (Scanmarker and screw)	ZBAD6931	-	-	10-12

DENTIUM IMPLANTIUM®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE D-I-ABCH 3.4 (Raw-Abutment® customizable hex and screw)	BSAF4881	30	-	10-12
SET ZZ BASE D-I-ABCH 3.8 (Raw-Abutment® customizable hex and screw)	BSAF4882	30	-	10-12
SET ZZ BASE D-I-ABCH 4.3 (Raw-Abutment® customizable hex and screw)	BSAF4883	30	-	10-12
SET ZZ BASE D-I-ABCH 4.8 (Raw-Abutment® customizable hex and screw)	BSAF4884	30	-	10-12
SET ZZ BASE D-I-ABCH-D14 3.4 (Raw-Abutment® customizable hex D14 and screw)	BSAF4881D14	30	-	10-12
SET ZZ BASE D-I-ABCH-D14 3.8 (Raw-Abutment® customizable hex D14 and screw)	BSAF4882D14	30	-	10-12
SET ZZ BASE D-I-ABCH-D14 4.3 (Raw-Abutment® customizable hex D14 and screw)	BSAF4883D14	30	-	10-12
SET ZZ BASE D-I-ABCH-D14 4.8 (Raw-Abutment® customizable hex D14 and screw)	BSAF4884D14	30	-	10-12
SET ZZ BASE D-I-CCTB 3.4 (Conical cemented titan base non hex and screw)	BSAD4831	30	-	10-12
SET ZZ BASE D-I-CCTB 3.8 (Conical cemented titan base non hex and screw)	BSAD4832	30	-	10-12
SET ZZ BASE D-I-CCTB 4.3 (Conical cemented titan base non hex and screw)	BSAD4833	30	-	10-12
SET ZZ BASE D-I-CCTB 4.8 (Conical cemented titan base non hex and screw)	BSAD4834	30	-	10-12
SET ZZ BASE D-I-NTB 3.4 (Narrow titan base non hex and screw)	BSAP4831	30	-	10-12
SET ZZ BASE D-I-NTB 3.8 (Narrow titan base non hex and screw)	BSAP4832	30	-	10-12
SET ZZ BASE D-I-NTB 4.3 (Narrow titan base non hex and screw)	BSAP4833	30	-	10-12
SET ZZ BASE D-I-NTB 4.8 (Narrow titan base non hex and screw)	BSAP4834	30	-	10-12
SET ZZ BASE D-I-NTBH 3.4 (Narrow titan base hex and screw)	BSAP4881	30	-	10-12
SET ZZ BASE D-I-NTBH 3.8 (Narrow titan base hex and screw)	BSAP4882	30	-	10-12
SET ZZ BASE D-I-NTBH 4.3 (Narrow titan base hex and screw)	BSAP4883	30	-	10-12
SET ZZ BASE D-I-NTBH 4.8 (Narrow titan base hex and screw)	BSAP4884	30	-	10-12
SET ZZ BASE D-I-PCTBH 3.4 (Parallel cemented titan base hex and screw)	BSAB4881	30	-	10-12
SET ZZ BASE D-I-PCTBH 3.8 (Parallel cemented titan base hex and screw)	BSAB4882	30	-	10-12
SET ZZ BASE D-I-PCTBH 4.3 (Parallel cemented titan base hex and screw)	BSAB4883	30	-	10-12
SET ZZ BASE D-I-PCTBH 4.8 (Parallel cemented titan base hex and screw)	BSAB4884	30	-	10-12
SET ZZ BASE D-I-SM 3.4 (Scanmarker and screw)	ZBAD4831	-	-	10-12
SET ZZ BASE D-I-SM 3.8 (Scanmarker and screw)	ZBAD4832	-	-	10-12
SET ZZ BASE D-I-SM 4.3 (Scanmarker and screw)	ZBAD4833	-	-	10-12
SET ZZ BASE D-I-SM 4.8 (Scanmarker and screw)	ZBAD4834	-	-	10-12
SET ZZ BASE D-I-WSM 3.4 (White Scanmarker and screw)	ZBCD4831	-	-	10-12
SET ZZ BASE D-I-WSM 3.8 (White Scanmarker and screw)	ZBCD4832	-	-	10-12

DENTIUM IMPLANTIUM®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE D-I-WSM 4.3 (White Scanmarker and screw)	ZBCD4833	30	-	10-12
SET ZZ BASE D-I-WSM 4.8 (White Scanmarker and screw)	ZBCD4834	30	-	10-12
ZZ BASE D-I-MI 3.4 (Model implant)	MIAG4801	-	-	10-12
ZZ BASE D-I-MI 3.8 (Model implant)	MIAG4802	-	-	10-12
ZZ BASE D-I-MI 4.3 (Model implant)	MIAG4803	-	-	10-12
ZZ BASE D-I-MI 4.8 (Model implant)	MIAG4804	-	-	10-12
ZZ BASE D-I-SC 3.4 (Screw)	BSAE4801	30	-	10-12
ZZ BASE D-I-SC 3.8 (Screw)	BSAE4802	30	-	10-12
ZZ BASE D-I-SC 4.3 (Screw)	BSAE4803	30	-	10-12
ZZ BASE D-I-SC 4.8 (Screw)	BSAE4804	30	-	10-12
ZZ BASE D-I-ZSC 3.4 (Zirconia screw)	BSAE4831	30	-	10-12
ZZ BASE D-I-ZSC 3.8 (Zirconia screw)	BSAE4832	30	-	10-12
ZZ BASE D-I-ZSC 4.3 (Zirconia screw)	BSAE4833	30	-	10-12
ZZ BASE D-I-ZSC 4.8 (Zirconia screw)	BSAE4834	30	-	10-12

DENTIUM SCREW ABUTMENT

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE D-SA-CCTB 4.5 (Conical cemented titan base non hex and screw)	BSAD4931	30	-	10-12
SET ZZ BASE D-SA-CCTB 5.5 (Conical cemented titan base non hex and screw)	BSAD4932	30	-	10-12
SET ZZ BASE D-SA-SM 4.5 (Scanmarker and screw)	ZBAD4931	-	-	10-12
SET ZZ BASE D-SA-SM 5.5 (Scanmarker and screw)	ZBAD4932	-	-	10-12
ZZ BASE D-SA-SC 4.5 (Screw)	BSAE4901	30	-	10-12
ZZ BASE D-SA-SC 5.5 (Screw)	BSAE4902	30	-	10-12
ZZ BASE D-SA-ZSC 4.5 (Zirconia screw)	BSAE4931	30	-	10-12
ZZ BASE D-SA-ZSC 5.5 (Zirconia screw)	BSAE4932	30	-	10-12

DENTSPLY IMPLANTS ANKYLOS®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE F-DA-ABC 3.5 (Raw-Abutment® customizable non hex and screw)	BSAF3031	-	-	10-12
SET ZZ BASE F-DA-ABC 4.5 (Raw-Abutment® customizable non hex and screw)	BSAF3032	-	-	10-12
SET ZZ BASE F-DA-ABC 5.5 (Raw-Abutment® customizable non hex and screw)	BSAF3033	-	-	10-12
SET ZZ BASE F-DA-ABC 7.0 (Raw-Abutment® customizable non hex and screw)	BSAF3034	-	-	10-12
SET ZZ BASE F-DA-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD3031	-	-	10-12
SET ZZ BASE F-DA-CCTB 4.5 (Conical cemented titan base non hex and screw)	BSAD3032	-	-	10-12
SET ZZ BASE F-DA-CCTB 5.5 (Conical cemented titan base non hex and screw)	BSAD3033	-	-	10-12
SET ZZ BASE F-DA-CCTB 7.0 (Conical cemented titan base non hex and screw)	BSAD3034	-	-	10-12
SET ZZ BASE F-DA-NTB 3.5 (Narrow titan base non hex and screw)	BSAP3031	-	-	10-12
SET ZZ BASE F-DA-NTB 4.5 (Narrow titan base non hex and screw)	BSAP3032	-	-	10-12
SET ZZ BASE F-DA-SM 3.5 (Scanmarker and screw)	ZBAD3031	-	-	10-12
SET ZZ BASE F-DA-SM 4.5 (Scanmarker and screw)	ZBAD3032	-	-	10-12
SET ZZ BASE F-DA-SM 5.5 (Scanmarker and screw)	ZBAD3033	-	-	10-12
SET ZZ BASE F-DA-SM 7.0 (Scanmarker and screw)	ZBAD3034	-	-	10-12
SET ZZ BASE F-DA-WSM 3.5 (White Scanmarker and screw)	ZBCD3031	-	-	10-12
SET ZZ BASE F-DA-WSM 4.5 (White Scanmarker and screw)	ZBCD3032	-	-	10-12
SET ZZ BASE F-DA-WSM 5.5 (White Scanmarker and screw)	ZBCD3033	-	-	10-12
SET ZZ BASE F-DA-WSM 7.0 (White Scanmarker and screw)	ZBCD3034	-	-	10-12
SET ZZ BASE F-DA-ZZMUA-L15 3.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH3031L15	25	35	10-12
SET ZZ BASE F-DA-ZZMUA-L15 4.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH3032L15	25	35	10-12
SET ZZ BASE F-DA-ZZMUA-L15 5.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH3033L15	25	35	10-12
SET ZZ BASE F-DA-ZZMUA-L15 7.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH3034L15	25	35	10-12
SET ZZ BASE F-DA-ZZMUA-L25 3.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH3031L25	25	35	10-12
SET ZZ BASE F-DA-ZZMUA-L25 4.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH3032L25	25	35	10-12
SET ZZ BASE F-DA-ZZMUA-L25 5.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH3033L25	25	35	10-12
SET ZZ BASE F-DA-ZZMUA-L25 7.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH3034L25	25	35	10-12
SET ZZ BASE F-DA-ZZMUA-L35 3.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH3031L35	25	35	10-12
SET ZZ BASE F-DA-ZZMUA-L35 4.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH3032L35	25	35	10-12
SET ZZ BASE F-DA-ZZMUA-L35 5.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH3033L35	25	35	10-12
SET ZZ BASE F-DA-ZZMUA-L35 7.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH3034L35	25	35	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

DENTSPLY IMPLANTS ANKYLOS®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE F-DA-ZZMUA-L45 3.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH3031L45	25	35	10-12
SET ZZ BASE F-DA-ZZMUA-L45 4.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH3032L45	25	35	10-12
SET ZZ BASE F-DA-ZZMUA-L45 5.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH3033L45	25	35	10-12
SET ZZ BASE F-DA-ZZMUA-L45 7.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH3034L45	25	35	10-12
ZZ BASE F-DA-SC 3.5 (Screw)	BSAE3001	-	-	10-12
ZZ BASE F-DA-SC 4.5 (Screw)	BSAE3002	-	-	10-12
ZZ BASE F-DA-SC 5.5 (Screw)	BSAE3003	-	-	10-12
ZZ BASE F-DA-SC 7.0 (Screw)	BSAE3004	-	-	10-12

DENTSPLY IMPLANTS ANKYLOS® BALANCE BASE A

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE F-DA-CCTB BBAN (Conical cemented titan base non hex and screw)	BSAD8231	25	-	10-12
SET ZZ BASE F-DA-SM BBAN (Scanmarker and screw)	ZBAD8231	-	-	10-12
SET ZZ BASE F-DA-WSM BBAN (White Scanmarker and screw)	ZBCD8231	-	-	10-12
ZZ BASE F-DA-SC BBAN (Screw)	BSAE8201	25	-	10-12
ZZ BASE F-DA-ZSC BBAN (Zirconia Screw)	BSAE8231	25	-	10-12

DENTSPLY IMPLANTS XIVE®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE F-DX-ABCH-ID 3.4 (Raw-Abutment® customizable hex ID and screw)	BSAF1091	25	-	10-12
SET ZZ BASE F-DX-ABCH-ID 3.8 (Raw-Abutment® customizable hex ID and screw)	BSAF1092	25	-	10-12
SET ZZ BASE F-DX-ABCH-ID 4.5 (Raw-Abutment® customizable hex ID and screw)	BSAF1093	25	-	10-12
SET ZZ BASE F-DX-ABCH-ID 5.5 (Raw-Abutment® customizable hex ID and screw)	BSAF1094	25	-	10-12
SET ZZ BASE F-DX-CCTB 3.0 (Conical cemented titan base non hex and screw)	BSAD1035	24	-	10-12
SET ZZ BASE F-DX-CCTB 3.4 (Conical cemented titan base non hex and screw)	BSAD1031	25	-	10-12
SET ZZ BASE F-DX-CCTB 3.8 (Conical cemented titan base non hex and screw)	BSAD1032	25	-	10-12
SET ZZ BASE F-DX-CCTB 4.5 (Conical cemented titan base non hex and screw)	BSAD1033	25	-	10-12
SET ZZ BASE F-DX-CCTB 5.5 (Conical cemented titan base non hex and screw)	BSAD1034	25	-	10-12
SET ZZ BASE F-DX-NTB 3.0 (Narrow titan base non hex and screw)	BSAP1035	24	-	10-12
SET ZZ BASE F-DX-NTB 3.4 (Narrow titan base non hex and screw)	BSAP1031	25	-	10-12
SET ZZ BASE F-DX-NTB 3.8 (Narrow titan base non hex and screw)	BSAP1032	25	-	10-12
SET ZZ BASE F-DX-NTBH 3.0 (Narrow titan base hex and screw)	BSAP1085	24	-	10-12
SET ZZ BASE F-DX-NTBH 3.4 (Narrow titan base hex and screw)	BSAP1081	25	-	10-12
SET ZZ BASE F-DX-NTBH 3.8 (Narrow titan base hex and screw)	BSAP1082	25	-	10-12
SET ZZ BASE F-DX-PCTBH 3.0 (Parallel cemented titan base hex and screw)	BSAB1085	24	-	10-12
SET ZZ BASE F-DX-PCTBH 3.4 (Parallel cemented titan base hex and screw)	BSAB1081	25	-	10-12
SET ZZ BASE F-DX-PCTBH 3.8 (Parallel cemented titan base hex and screw)	BSAB1082	25	-	10-12
SET ZZ BASE F-DX-PCTBH 4.5 (Parallel cemented titan base hex and screw)	BSAB1083	25	-	10-12
SET ZZ BASE F-DX-PCTBH 5.5 (Parallel cemented titan base hex and screw)	BSAB1084	25	-	10-12
SET ZZ BASE F-DX-SM 3.0 (Scanmarker and screw)	ZBAD1035	-	-	10-12
SET ZZ BASE F-DX-SM 3.4 (Scanmarker and screw)	ZBAD1031	-	-	10-12
SET ZZ BASE F-DX-SM 3.8 (Scanmarker and screw)	ZBAD1032	-	-	10-12
SET ZZ BASE F-DX-SM 4.5 (Scanmarker and screw)	ZBAD1033	-	-	10-12
SET ZZ BASE F-DX-SM 5.5 (Scanmarker and screw)	ZBAD1034	-	-	10-12
SET ZZ BASE F-DX-WSM 3.0 (White Scanmarker and screw)	ZBCD1035	-	-	10-12
SET ZZ BASE F-DX-WSM 3.4 (White Scanmarker and screw)	ZBCD1031	-	-	10-12
SET ZZ BASE F-DX-WSM 3.8 (White Scanmarker and screw)	ZBCD1032	-	-	10-12
SET ZZ BASE F-DX-WSM 4.5 (White Scanmarker and screw)	ZBCD1033	-	-	10-12
SET ZZ BASE F-DX-WSM 5.5 (White Scanmarker and screw)	ZBCD1034	-	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

DENTSPLY IMPLANTS XIVE®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
ZZ BASE F-DX-MI 3.0 (Model implant)	MIAG1005	-	-	10-12
ZZ BASE F-DX-MI 3.4 (Model implant)	MIAG1001	-	-	10-12
ZZ BASE F-DX-MI 3.8 (Model implant)	MIAG1002	-	-	10-12
ZZ BASE F-DX-MI 4.5 (Model implant)	MIAG1003	-	-	10-12
ZZ BASE F-DX-MI 5.5 (Model implant)	MIAG1004	-	-	10-12
ZZ BASE F-DX-SC 3.0 (Screw)	BSAE1005	24	-	10-12
ZZ BASE F-DX-SC 3.4 (Screw)	BSAE1001	25	-	10-12
ZZ BASE F-DX-SC 3.8 (Screw)	BSAE1002	25	-	10-12
ZZ BASE F-DX-SC 4.5 (Screw)	BSAE1003	25	-	10-12
ZZ BASE F-DX-SC 5.5 (Screw)	BSAE1004	25	-	10-12
ZZ BASE F-DX-ZSC 3.0 (Zirconia screw)	BSAE1035	24	-	10-12
ZZ BASE F-DX-ZSC 3.4 (Zirconia screw)	BSAE1031	25	-	10-12
ZZ BASE F-DX-ZSC 3.8 (Zirconia screw)	BSAE1032	25	-	10-12
ZZ BASE F-DX-ZSC 4.5 (Zirconia screw)	BSAE1033	25	-	10-12
ZZ BASE F-DX-ZSC 5.5 (Zirconia screw)	BSAE1034	25	-	10-12

DENTSPLY IMPLANTS XIVE® MP/TG

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE F-MP/TG-CCTB 3.4 (Conical cemented titan base non hex and screw)	BSAD4031	24	-	10-12
SET ZZ BASE F-MP/TG-CCTB 3.8 (Conical cemented titan base non hex and screw)	BSAD4032	24	-	10-12
SET ZZ BASE F-MP/TG-CCTB 4.5 (Conical cemented titan base non hex and screw)	BSAD4033	24	-	10-12
SET ZZ BASE F-MP/TG-CCTB 5.5 (Conical cemented titan base non hex and screw)	BSAD4034	24	-	10-12
SET ZZ BASE F-MP/TG-SM 3.4 (Scanmarker and screw)	ZBAD4031	-	-	10-12
SET ZZ BASE F-MP/TG-SM 3.8 (Scanmarker and screw)	ZBAD4032	-	-	10-12
SET ZZ BASE F-MP/TG-SM 4.5 (Scanmarker and screw)	ZBAD4033	-	-	10-12
SET ZZ BASE F-MP/TG-SM 5.5 (Scanmarker and screw)	ZBAD4034	-	-	10-12
ZZ BASE F-MP/TG-SC 3.4 (Screw)	BSAE4001	24	-	10-12
ZZ BASE F-MP/TG-SC 3.8 (Screw)	BSAE4002	24	-	10-12
ZZ BASE F-MP/TG-SC 4.5 (Screw)	BSAE4003	24	-	10-12
ZZ BASE F-MP/TG-SC 5.5 (Screw)	BSAE4004	24	-	10-12
ZZ BASE F-MP/TG-ZSC 3.4 (Zirconia screw)	BSAE4031	24	-	10-12
ZZ BASE F-MP/TG-ZSC 3.8 (Zirconia screw)	BSAE4032	24	-	10-12
ZZ BASE F-MP/TG-ZSC 4.5 (Zirconia screw)	BSAE4033	24	-	10-12
ZZ BASE F-MP/TG-ZSC 5.5 (Zirconia screw)	BSAE4034	24	-	10-12

DYNA® OCTALOCK®/HELIX®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE DYN-OL-CCTB 3.6-4.0 (Conical cemented titan base non hex and screw)	BSAD4131	35	-	10-12
SET ZZ BASE DYN-OL-PCTBH 3.6-4.0 (Parallel cemented titan base hex and screw)	BSAB4181	35	-	10-12
SET ZZ BASE DYN-OL-SM 3.6-4.0 (Scanmarker and screw)	ZBAD4131	-	-	10-12
ZZ BASE DYN-OL-SC 3.6-4.0 (Screw)	BSAE4101	35	-	10-12
ZZ BASE DYN-OL-ZSC 3.6-4.0 (Zirconia screw)	BSAE4131	35	-	10-12

FAIRIMPLANT™ FAIRTWO™

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE FI-FT-ABCH L (Raw-Abutment® customizable hex and screw)	BSAF9582	20	-	10-12
SET ZZ BASE FI-FT-ABCH S (Raw-Abutment® customizable hex and screw)	BSAF9581	20	-	10-12
SET ZZ BASE FI-FT-CCTB L (Conical cemented titan base non hex and screw)	BSAD9532	20	-	10-12
SET ZZ BASE FI-FT-CCTB S (Conical cemented titan base non hex and screw)	BSAD9531	20	-	10-12
SET ZZ BASE FI-FT-PCTBH L (Parallel cemented titan base hex and screw)	BSAB9582	20	-	10-12
SET ZZ BASE FI-FT-PCTBH S (Parallel cemented titan base hex and screw)	BSAB9581	20	-	10-12
SET ZZ BASE FI-FT-SM L (Scanmarker and screw)	ZBAD9532	-	-	10-12
SET ZZ BASE FI-FT-SM S (Scanmarker and screw)	ZBAD9531	-	-	10-12
SET ZZ BASE FI-FT-WSM L (White Scanmarker and screw)	ZBCD9532	-	-	10-12
SET ZZ BASE FI-FT-WSM S (White Scanmarker and screw)	ZBCD9531	-	-	10-12
ZZ BASE FI-FT-SC L (Screw)	BSAE9502	20	-	10-12
ZZ BASE FI-FT-SC S (Screw)	BSAE9501	20	-	10-12
ZZ BASE FI-FT-ZSC L (Zirconia Screw)	BSAE9532	20	-	10-12
ZZ BASE FI-FT-ZSC S (Zirconia Screw)	BSAE9531	20	-	10-12

ICX-MULTI® KONZEPT

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE I-MK-CCTB (Conical cemented titan base non hex and screw)	BSAD10031	27	-	10-12
SET ZZ BASE I-MK-SM (Scanmarker and screw)	ZBAD10031	-	-	10-12
SET ZZ BASE I-MK-WSM (White Scanmarker and screw)	ZBCD10031	-	-	10-12
ZZ BASE I-MK-SC (Screw)	BSAE10001	27	-	10-12
ZZ BASE I-MK-ZSC (Zirconia screw)	BSAE10031	27	-	10-12

ICX-TEMPLANT®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE I-T-SM 3.45-4.80 (Scanmarker and screw)	ZBAD1931	-	-	10-12
SET ZZ BASE I-T-ABCH 3.45-4.8 (Raw-Abutment® customizable hex and screw)	BSAF1981	25	-	10-12
SET ZZ BASE I-T-ABCH-D14 3.45-4.8 (Raw-Abutment® customizable hex D14 and screw)	BSAF1981D14	25	-	10-12
SET ZZ BASE I-T-CCTB 3.45-4.8 (Conical cemented titan base non hex and screw)	BSAD1931	25	-	10-12
SET ZZ BASE I-T-NTB 3.45-4.8 (Narrow titan base non hex and screw)	BSAP1931	25	-	10-12
SET ZZ BASE I-T-NTBH 3.45-4.8 (Narrow titan base hex and screw)	BSAP1981	25	-	10-12
SET ZZ BASE I-T-PCTBH 3.45-4.80 (Parallel cemented titan base hex and screw)	BSAB1981	25	-	10-12
SET ZZ BASE I-T-WSM 3.45-4.80 (White Scanmarker and screw)	ZBCD1931	-	-	10-12
SET ZZ BASE I-T-ZZMUA-L15 3.45-4.8 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH1931L15	25	25	10-12
SET ZZ BASE I-T-ZZMUA-L25 3.45-4.8 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH1931L25	25	25	10-12
SET ZZ BASE I-T-ZZMUA-L35 3.45-4.8 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH1931L35	25	25	10-12
SET ZZ BASE I-T-ZZMUA-L45 3.45-4.8 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH1931L45	25	25	10-12
ZZ BASE I-T-MI 3.45-4.80 (Model implant)	MIAG1901	-	-	10-12
ZZ BASE I-T-SC 3.45-4.80 (Screw)	BSAE1901	25	-	10-12
ZZ BASE I-T-ZSC 3.45-4.80 (Zirconia screw)	BSAE1931	25	-	10-12
ZZ BASE I-T-ZSC-G 3.45-4.80 (Zirconia Screw Gold)	BSCE1931	25	-	10-12

IMPLANT DIRECT LEGACY™

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE ID-LE-ABCH 3.5 (Raw-Abutment® customizable hex and screw)	BSAF5782	30	-	10-12
SET ZZ BASE ID-LE-ABCH 4.5 (Raw-Abutment® customizable hex and screw)	BSAF5783	30	-	10-12
SET ZZ BASE ID-LE-ABCH 5.7 (Raw-Abutment® customizable hex and screw)	BSAF5784	30	-	10-12
SET ZZ BASE ID-LE-ABCH-D14 3.5 (Raw-Abutment® customizable hex D14 and screw)	BSAF5782D14	30	-	10-12
SET ZZ BASE ID-LE-ABCH-D14 4.5 (Raw-Abutment® customizable hex D14 and screw)	BSAF5783D14	30	-	10-12
SET ZZ BASE ID-LE-ABCH-D14 5.7 (Raw-Abutment® customizable hex D14 and screw)	BSAF5784D14	30	-	10-12
SET ZZ BASE ID-LE-ABCH-D14-ID 3.0 (Raw-Abutment® customizable hex D14 ID and screw)	BSAF5791D14	30	-	10-12
SET ZZ BASE ID-LE-ABCH-ID 3.0 (Raw-Abutment® customizable hex ID and screw)	BSAF5791	30	-	10-12
SET ZZ BASE ID-LE-CCTB 3.0 (Conical cemented titan base non hex and screw)	BSAD5731	30	-	10-12
SET ZZ BASE ID-LE-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD5732	30	-	10-12
SET ZZ BASE ID-LE-CCTB 4.5 (Conical cemented titan base non hex and screw)	BSAD5733	30	-	10-12
SET ZZ BASE ID-LE-CCTB 5.7 (Conical cemented titan base non hex and screw)	BSAD5734	30	-	10-12
SET ZZ BASE ID-LE-NTB 3.0 (Narrow titan base non hex and screw)	BSAP5731	30	-	10-12
SET ZZ BASE ID-LE-NTB 3.5 (Narrow titan base non hex and screw)	BSAP5732	30	-	10-12
SET ZZ BASE ID-LE-NTBH 3.0 (Narrow titan base hex and screw)	BSAP5781	30	-	10-12
SET ZZ BASE ID-LE-NTBH 3.5 (Narrow titan base hex and screw)	BSAP5782	30	-	10-12
SET ZZ BASE ID-LE-PCTBH 3.0 (Parallel cemented titan base hex and screw)	BSAB5781	30	-	10-12
SET ZZ BASE ID-LE-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB5782	30	-	10-12
SET ZZ BASE ID-LE-PCTBH 4.5 (Parallel cemented titan base hex and screw)	BSAB5783	30	-	10-12
SET ZZ BASE ID-LE-PCTBH 5.7 (Parallel cemented titan base hex and screw)	BSAB5784	30	-	10-12
SET ZZ BASE ID-LE-SM 3.0 (Scanmarker and screw)	ZBAD5731	-	-	10-12
SET ZZ BASE ID-LE-SM 3.5 (Scanmarker and screw)	ZBAD5732	-	-	10-12
SET ZZ BASE ID-LE-SM 4.5 (Scanmarker and screw)	ZBAD5733	-	-	10-12
SET ZZ BASE ID-LE-SM 5.7 (Scanmarker and screw)	ZBAD5734	-	-	10-12
SET ZZ BASE ID-LE-WSM 3.0 (White Scanmarker and screw)	ZBCD5731	-	-	10-12
SET ZZ BASE ID-LE-WSM 3.5 (White Scanmarker and screw)	ZBCD5732	-	-	10-12
SET ZZ BASE ID-LE-WSM 4.5 (White Scanmarker and screw)	ZBCD5733	-	-	10-12
SET ZZ BASE ID-LE-WSM 5.7 (White Scanmarker and screw)	ZBCD5734	-	-	10-12
ZZ BASE ID-LE-SC 3.0 (Screw)	BSAE5701	30	-	10-12
ZZ BASE ID-LE-SC 3.5 (Screw)	BSAE5702	30	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

IMPLANT DIRECT LEGACY™

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
ZZ BASE ID-LE-SC 4.5 (Screw)	BSAE5703	30	-	10-12
ZZ BASE ID-LE-SC 5.7 (Screw)	BSAE5704	30	-	10-12
ZZ BASE ID-LE-ZSC 3.0 (Zirconia screw)	BSAE5731	30	-	10-12
ZZ BASE ID-LE-ZSC 3.5 (Zirconia screw)	BSAE5732	30	-	10-12
ZZ BASE ID-LE-ZSC 4.5 (Zirconia screw)	BSAE5733	30	-	10-12
ZZ BASE ID-LE-ZSC 5.7 (Zirconia screw)	BSAE5734	30	-	10-12

IMPLANT DIRECT OVERDENTURE ABUTMENT

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE ID-ODA-CCTB (Conical cemented titan base non hex and screw)	BSAD4431	30	-	10-12
SET ZZ BASE ID-ODA-SM (Scanmarker and screw)	ZBAD4431	-	-	10-12
ZZ BASE ID-ODA-SC (Screw)	BSAE4401	30	-	10-12
ZZ BASE ID-ODA-ZSC (Zirconia Screw)	BSAE4431	30	-	10-12

INTRALOCK® INTERNAL IMPLANTS

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE IL-II-CCTB Sq (Conical cemented titan base non hex and screw)	BSAD8131	30	-	10-12
SET ZZ BASE IL-II-CCTB St (Conical cemented titan base non hex and screw)	BSAD8132	30	-	10-12
SET ZZ BASE IL-II-CCTB Wd (Conical cemented titan base non hex and screw)	BSAD8133	30	-	10-12
SET ZZ BASE IL-II-PCTBH Sq (Parallel cemented titan base hex and screw)	BSAB8181	30	-	10-12
SET ZZ BASE IL-II-PCTBH St (Parallel cemented titan base hex and screw)	BSAB8182	30	-	10-12
SET ZZ BASE IL-II-PCTBH Wd (Parallel cemented titan base hex and screw)	BSAB8183	30	-	10-12
SET ZZ BASE IL-II-SM Sq (Scanmarker and screw)	ZBAD8131	-	-	10-12
SET ZZ BASE IL-II-SM St (Scanmarker and screw)	ZBAD8132	-	-	10-12
SET ZZ BASE IL-II-SM Wd (Scanmarker and screw)	ZBAD8133	-	-	10-12
SET ZZ BASE IL-II-WSM Sq (White Scanmarker and screw)	ZBCD8131	-	-	10-12
SET ZZ BASE IL-II-WSM St (White Scanmarker and screw)	ZBCD8132	-	-	10-12
SET ZZ BASE IL-II-WSM Wd (White Scanmarker and screw)	ZBCD8133	-	-	10-12
ZZ BASE IL-II-SC Sq (Screw)	BSAE8101	30	-	10-12
ZZ BASE IL-II-SC St (Screw)	BSAE8102	30	-	10-12
ZZ BASE IL-II-SC Wd (Screw)	BSAE8103	30	-	10-12
ZZ BASE IL-II-ZSC Sq (Zirconia Screw)	BSAE8131	30	-	10-12
ZZ BASE IL-II-ZSC St (Zirconia Screw)	BSAE8132	30	-	10-12
ZZ BASE IL-II-ZSC Wd (Zirconia Screw)	BSAE8133	30	-	10-12

INTRALOCK® UNIHEX™

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE IL-UH-CCTB 3.4 (Conical cemented titan base non hex and screw)	BSAD7931	30	-	10-12
SET ZZ BASE IL-UH-CCTB 4.0 (Conical cemented titan base non hex and screw)	BSAD7932	30	-	10-12
SET ZZ BASE IL-UH-CCTB 4.75 (Conical cemented titan base non hex and screw)	BSAD7933	30	-	10-12
SET ZZ BASE IL-UH-PCTBH 3.4 (Parallel cemented titan base hex and screw)	BSAB7981	30	-	10-12
SET ZZ BASE IL-UH-PCTBH 4.0 (Parallel cemented titan base hex and screw)	BSAB7982	30	-	10-12
SET ZZ BASE IL-UH-PCTBH 4.75 (Parallel cemented titan base hex and screw)	BSAB7983	30	-	10-12
SET ZZ BASE IL-UH-SM 3.4 (Scanmarker and screw)	ZBAD7931	-	-	10-12
SET ZZ BASE IL-UH-SM 4.0 (Scanmarker and screw)	ZBAD7932	-	-	10-12
SET ZZ BASE IL-UH-SM 4.75 (Scanmarker and screw)	ZBAD7933	-	-	10-12
SET ZZ BASE IL-UH-WSM 3.4 (White Scanmarker and screw)	ZBCD7931	-	-	10-12
SET ZZ BASE IL-UH-WSM 4.0 (White Scanmarker and screw)	ZBCD7932	-	-	10-12
SET ZZ BASE IL-UH-WSM 4.75 (White Scanmarker and screw)	ZBCD7933	-	-	10-12
ZZ BASE IL-UH-SC 3.4 (Screw)	BSAE7901	30	-	10-12
ZZ BASE IL-UH-SC 4.0 (Screw)	BSAE7902	30	-	10-12
ZZ BASE IL-UH-SC 4.75 (Screw)	BSAE7903	30	-	10-12
ZZ BASE IL-UH-ZSC 3.4 (Zirconia Screw)	BSAE7931	30	-	10-12
ZZ BASE IL-UH-ZSC 4.0 (Zirconia Screw)	BSAE7932	30	-	10-12
ZZ BASE IL-UH-ZSC 4.75 (Zirconia Screw)	BSAE7933	30	-	10-12

K3PRO®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE K3-ABCH Standard 3 Pro (Raw-Abutment® customizable hex and screw)	BSAF6382	-	-	10-12
SET ZZ BASE K3-CCTB Mini 2 Pro (Conical cemented titan base non hex and screw)	BSAD6331	-	-	10-12
SET ZZ BASE K3-CCTB Standard 3 Pro (Conical cemented titan base non hex and screw)	BSAD6332	-	-	10-12
SET ZZ BASE K3-CCTB-L15 Mini 2 Pro (Conical cemented titan base non hex L15 and screw)	BSAD6331L15	-	-	10-12
SET ZZ BASE K3-CCTB-L15 Standard 3 Pro (Conical cemented titan base non hex L15 and screw)	BSAD6332L15	-	-	10-12
SET ZZ BASE K3-PCTBH Mini 2 Pro (Parallel cemented titan base hex and screw)	BSAB6381	-	-	10-12
SET ZZ BASE K3-PCTBH Standard 3 Pro (Parallel cemented titan base hex and screw)	BSAB6382	-	-	10-12
SET ZZ BASE K3-PCTBH-L15 Mini 2 Pro (Parallel cemented titan base hex L15 and screw)	BSAB6391L15	-	-	10-12
SET ZZ BASE K3-PCTBH-L15 Standard 3 Pro (Parallel cemented titan base hex L15 and screw)	BSAB6392L15	-	-	10-12
SET ZZ BASE K3-SM Mini 2 Pro (Scanmarker and screw)	ZBAD6331	-	-	10-12
SET ZZ BASE K3-SM Standard 3 Pro (Scanmarker and screw)	ZBAD6332	-	-	10-12
SET ZZ BASE K3-WSM Mini 2 Pro (White Scanmarker and screw)	ZBCD6331	-	-	10-12
SET ZZ BASE K3-WSM Standard 3 Pro (White Scanmarker and screw)	ZBCD6332	-	-	10-12
ZZ BASE K3-SC Mini 2 Pro (Screw)	BSAE6301	-	-	10-12
ZZ BASE K3-SC Standard 3 Pro (Screw)	BSAE6302	-	-	10-12
ZZ BASE K3-ZSC Mini 2 Pro (Zirconia Screw)	BSAE6331	-	-	10-12
ZZ BASE K3-ZSC Standard 3 Pro (Zirconia Screw)	BSAE6332	-	-	10-12

KLOCKNER® ESSENTIAL® CONE

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE K-EC-CCTB 3.5/4.0/4.5 (Conical cemented titan base non hex and screw)	BSAD2831	30	-	10-12
SET ZZ BASE K-EC-CCTB 4.8 (Conical cemented titan base non hex and screw)	BSAD2832	30	-	10-12
SET ZZ BASE K-EC-PCTBH 3.5/4.0/4.5 (Parallel cemented titan base hex and screw)	BSAB2881	30	-	10-12
SET ZZ BASE K-EC-PCTBH 4.8 (Parallel cemented titan base hex and screw)	BSAB2882	30	-	10-12
SET ZZ BASE K-EC-SM 3.5/4.0/4.5 (Scanmarker and screw)	ZBAD2831	-	-	10-12
SET ZZ BASE K-EC-SM 4.8 (Scanmarker and screw)	ZBAD2832	-	-	10-12
ZZ BASE K-EC-SC 3.5/4.0/4.5 (Screw)	BSAE2801	30	-	10-12
ZZ BASE K-EC-SC 4.8 (Screw)	BSAE2802	30	-	10-12
ZZ BASE K-EC-ZSC 3.5/4.0/4.5 (Zirconia screw)	BSAE2831	30	-	10-12
ZZ BASE K-EC-ZSC 4.8 (Zirconia screw)	BSAE2832	30	-	10-12

KLOCKNER® NK2/SK2

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE K-NK2/SK2-CCTB 3.3/3.8/4.2 (Conical cemented titan base non hex and screw)	BSAD2931	25	-	10-12
SET ZZ BASE K-NK2/SK2-CCTB 4.7 (Conical cemented titan base non hex and screw)	BSAD2932	25	-	10-12
SET ZZ BASE K-NK2/SK2-PCTBH 3.3/3.8/4.2 (Parallel cemented titan base hex and screw)	BSAB2981	25	-	10-12
SET ZZ BASE K-NK2/SK2-PCTBH 4.7 (Parallel cemented titan base hex and screw)	BSAB2982	25	-	10-12
SET ZZ BASE K-NK2/SK2-SM 3.3/3.8/4.2 (Scanmarker and screw)	ZBAD2931	-	-	10-12
SET ZZ BASE K-NK2/SK2-SM 4.7 (Scanmarker and screw)	ZBAD2932	-	-	10-12
ZZ BASE K-NK2/SK2-SC 3.3/3.8/4.2 (Screw)	BSAE2901	25	-	10-12
ZZ BASE K-NK2/SK2-SC 4.7 (Screw)	BSAE2902	25	-	10-12
ZZ BASE K-NK2/SK2-ZSC 3.3/3.8/4.2 (Zirconia screw)	BSAE2931	25	-	10-12
ZZ BASE K-NK2/SK2-ZSC 4.7 (Zirconia screw)	BSAE2932	25	-	10-12

MEDENTIKA® MEDENTIBASE®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE MD-MB-CCTB (Conical cemented titan base non hex and screw)	BSAD8631	15	-	10-12
SET ZZ BASE MD-MB-SM (Scanmarker and screw)	ZBAD8631	-	-	10-12
SET ZZ BASE MD-MB-WSM (White Scanmarker and screw)	ZBCD8631	-	-	10-12
ZZ BASE MD-MB-MI (Model implant)	MIAG8601	-	-	10-12
ZZ BASE MD-MB-SC (Screw)	BSAE8601	15	-	10-12
ZZ BASE MD-MB-ZSC (Zirconia Screw)	BSAE8631	15	-	10-12

MEDENTIKA® M-IMPLANT

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE MD-MI-ABCH 3.0 (Raw-Abutment® customizable hex and screw)	BSAF8781	15	-	10-12
SET ZZ BASE MD-MI-ABCH 3.5-5.0 (Raw-Abutment® customizable hex and screw)	BSAF8782	25	-	10-12
SET ZZ BASE MD-MI-CCTB 3.0 (Conical cemented titan base non hex and screw)	BSAD8731	15	-	10-12
SET ZZ BASE MD-MI-CCTB 3.5-5.0 (Conical cemented titan base non hex and screw)	BSAD8732	25	-	10-12
SET ZZ BASE MD-MI-NTB 3.0 (Narrow titan base non hex and screw)	BSAP8731	15	-	10-12
SET ZZ BASE MD-MI-NTB 3.5-5.0 (Narrow titan base non hex and screw)	BSAP8732	25	-	10-12
SET ZZ BASE MD-MI-NTBH 3.0 (Narrow titan base hex and screw)	BSAP8781	15	-	10-12
SET ZZ BASE MD-MI-NTBH 3.5-5.0 (Narrow titan base hex and screw)	BSAP8782	25	-	10-12
SET ZZ BASE MD-MI-PCTBH 3.0 (Parallel cemented titan base hex and screw)	BSAB8781	25	-	10-12
SET ZZ BASE MD-MI-PCTBH 3.5-5.0 (Parallel cemented titan base hex and screw)	BSAB8782	25	-	10-12
SET ZZ BASE MD-MI-SM 3.0 (Scanmarker and screw)	ZBAD8731	-	-	10-12
SET ZZ BASE MD-MI-SM 3.5-5.0 (Scanmarker and screw)	ZBAD8732	-	-	10-12
SET ZZ BASE MD-MI-WSM 3.0 (White Scanmarker and screw)	ZBCD8731	-	-	10-12
SET ZZ BASE MD-MI-WSM 3.5-5.0 (White Scanmarker and screw)	ZBCD8732	-	-	10-12
SET ZZ BASE MD-MI-ZZMUA-L15 3.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH8731L15	25	30	10-12
SET ZZ BASE MD-MI-ZZMUA-L15 3.5-5.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH8732L15	25	30	10-12
SET ZZ BASE MD-MI-ZZMUA-L25 3.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH8731L25	25	30	10-12
SET ZZ BASE MD-MI-ZZMUA-L25 3.5-5.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH8732L25	25	30	10-12
SET ZZ BASE MD-MI-ZZMUA-L35 3.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH8731L35	25	30	10-12
SET ZZ BASE MD-MI-ZZMUA-L35 3.5-5.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH8732L35	25	30	10-12
SET ZZ BASE MD-MI-ZZMUA-L45 3.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH8731L45	25	30	10-12
SET ZZ BASE MD-MI-ZZMUA-L45 3.5-5.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH8732L45	25	30	10-12
ZZ BASE MD-MI-MI 3.0 (Model implant)	MIAG8701	-	-	10-12
ZZ BASE MD-MI-MI 3.5-5.0 (Model implant)	MIAG8702	-	-	10-12
ZZ BASE MD-MI-SC 3.0 (Screw)	BSAE8701	15	-	10-12
ZZ BASE MD-MI-SC 3.5-5.0 (Screw)	BSAE8702	25	-	10-12
ZZ BASE MD-MI-ZSC 3.0 (Zirconia Screw)	BSAE8731	15	-	10-12
ZZ BASE MD-MI-ZSC 3.5-5.0 (Zirconia Screw)	BSAE8732	25	-	10-12

MEGAGEN ANYRIDGE®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE ME-AR-ABCH 3.5-8.0 (Raw-Abutment® customizable hex and screw)	BSAF6881	30	-	10-12
SET ZZ BASE ME-AR-CCTB 3.5-8.0 (Conical cemented titan base non hex and screw)	BSAD6831	30	-	10-12
SET ZZ BASE ME-AR-NTB 3.5-8.0 (Narrow titan base non hex and screw)	BSAP6831	30	-	10-12
SET ZZ BASE ME-AR-NTBH 3.5-8.0 (Narrow titan base hex and screw)	BSAP6881	30	-	10-12
SET ZZ BASE ME-AR-PCTBH 3.5-8.0 (Parallel cemented titan base hex and screw)	BSAB6881	30	-	10-12
SET ZZ BASE ME-AR-SM 3.5-8.0 (Scanmarker and screw)	ZBAD6831	-	-	10-12
SET ZZ BASE ME-AR-WSM 3.5-8.0 (White Scanmarker and screw)	ZBCD6831	-	-	10-12
ZZ BASE ME-AR-SC 3.5-8.0 (Screw)	BSAE6801	30	-	10-12
ZZ BASE ME-AR-ZSC 3.5-8.0 (Zirconia Screw)	BSAE6831	-	-	10-12

MIS® C1

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE MS-C-ABCH Narrow (Raw-Abutment® customizable hex and screw)	BSAF5881	30	-	10-12
SET ZZ BASE MS-C-ABCH Standard (Raw-Abutment® customizable hex and screw)	BSAF5882	30	-	10-12
SET ZZ BASE MS-C-ABCH Wide (Raw-Abutment® customizable hex and screw)	BSAF5883	30	-	10-12
SET ZZ BASE MS-C-CCTB Narrow (Conical cemented titan base non hex and screw)	BSAD5831	30	-	10-12
SET ZZ BASE MS-C-CCTB Standard (Conical cemented titan base non hex and screw)	BSAD5832	30	-	10-12
SET ZZ BASE MS-C-CCTB Wide (Conical cemented titan base non hex and screw)	BSAD5833	30	-	10-12
SET ZZ BASE MS-C-NTB Narrow (Narrow titan base non hex and screw)	BSAP5831	30	-	10-12
SET ZZ BASE MS-C-NTB Standard (Narrow titan base non hex and screw)	BSAP5832	30	-	10-12
SET ZZ BASE MS-C-NTBH Narrow (Narrow titan base hex and screw)	BSAP5881	30	-	10-12
SET ZZ BASE MS-C-NTBH Standard (Narrow titan base hex and screw)	BSAP5882	30	-	10-12
SET ZZ BASE MS-C-PCTBH Narrow (Parallel cemented titan base hex and screw)	BSAB5881	30	-	10-12
SET ZZ BASE MS-C-PCTBH Standard (Parallel cemented titan base hex and screw)	BSAB5882	30	-	10-12
SET ZZ BASE MS-C-PCTBH Wide (Parallel cemented titan base hex and screw)	BSAB5883	30	-	10-12
SET ZZ BASE MS-C-SM Narrow (Scanmarker and screw)	ZBAD5831	-	-	10-12
SET ZZ BASE MS-C-SM Standard (Scanmarker and screw)	ZBAD5832	-	-	10-12
SET ZZ BASE MS-C-SM Wide (Scanmarker and screw)	ZBAD5833	-	-	10-12
SET ZZ BASE MS-C-WSM Narrow (White Scanmarker and screw)	ZBCD5831	-	-	10-12
SET ZZ BASE MS-C-WSM Standard (White Scanmarker and screw)	ZBCD5832	-	-	10-12
SET ZZ BASE MS-C-WSM Wide (White Scanmarker and screw)	ZBCD5833	-	-	10-12
ZZ BASE MS-C-MI Narrow (Model implant)	MIAG5801	-	-	10-12
ZZ BASE MS-C-MI Standard (Model implant)	MIAG5802	-	-	10-12
ZZ BASE MS-C-MI Wide (Model implant)	MIAG5803	-	-	10-12
ZZ BASE MS-C-SC Narrow (Screw)	BSAE5801	30	-	10-12
ZZ BASE MS-C-SC Standard (Screw)	BSAE5802	30	-	10-12
ZZ BASE MS-C-SC Wide (Screw)	BSAE5803	30	-	10-12
ZZ BASE MS-C-ZSC Narrow (Zirconia screw)	BSAE5831	30	-	10-12
ZZ BASE MS-C-ZSC Standard (Zirconia screw)	BSAE5832	30	-	10-12
ZZ BASE MS-C-ZSC Wide (Zirconia screw)	BSAE5833	30	-	10-12

MIS® MULTI UNIT ABUTMENT

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE MS-MUA-SM (Scanmarker and screw)	ZBAD1731	-	-	10-12
SET ZZ BASE MS-MUA-CCTB (Conical cemented titan base non hex and screw)	BSAD1731	20	-	10-12
SET ZZ BASE MS-MUA-WSM (White Scanmarker and screw)	ZBCD1731	-	-	10-12
ZZ BASE MS-MUA-SC (Screw)	BSAE1701	20	-	10-12
ZZ BASE MS-MUA-SC-G (Screw Gold)	BSCE1701	20	-	10-12
ZZ BASE MS-MUA-ZSC (Zirconia screw)	BSAE1731	20	-	10-12
ZZ BASE MS-MUA-ZSC-G (Zirconia Screw Gold)	BSCE1731	20	-	10-12

MIS® MULTI UNIT SYSTEM

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE MS-MUS-CCTB (Conical cemented titan base non hex and screw)	BSAD5931	30	-	10-12
SET ZZ BASE MS-MUS-SM (Scanmarker and screw)	ZBAD5931	-	-	10-12
SET ZZ BASE MS-MUS-WSM (White Scanmarker and screw)	ZBCD5931	-	-	10-12
ZZ BASE MS-MUS-MI (Model implant)	MIAG5901	-	-	10-12
ZZ BASE MS-MUS-SC (Screw)	BSAE5901	30	-	10-12
ZZ BASE MS-MUS-ZSC (Zirconia screw)	BSAE5931	30	-	10-12

MIS® SEVEN

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE MS-S-SM Standard (Scanmarker and screw)	ZBAD1631	-	-	10-12
SET ZZ BASE MS-S-SM Wide (Scanmarker and screw)	ZBAD1632	-	-	10-12
SET ZZ BASE MS-S-ABCH Narrow (Raw-Abutment® customizable hex and screw)	BSAF1683	30	-	10-12
SET ZZ BASE MS-S-ABCH Standard (Raw-Abutment® customizable hex and screw)	BSAF1681	30	-	10-12
SET ZZ BASE MS-S-ABCH Wide (Raw-Abutment® customizable hex and screw)	BSAF1682	30	-	10-12
SET ZZ BASE MS-S-ABCH-D14 Narrow (Raw-Abutment® customizable hex D14 and screw)	BSAF1683D14	30	-	10-12
SET ZZ BASE MS-S-ABCH-D14 Standard (Raw-Abutment® customizable hex D14 and screw)	BSAF1681D14	30	-	10-12
SET ZZ BASE MS-S-ABCH-D14 Wide (Raw-Abutment® customizable hex D14 and screw)	BSAF1682D14	30	-	10-12
SET ZZ BASE MS-S-CCTB Narrow (Conical cemented titan base non hex and screw)	BSAD1633	30	-	10-12
SET ZZ BASE MS-S-CCTB Standard (Conical cemented titan base non hex and screw)	BSAD1631	30	-	10-12
SET ZZ BASE MS-S-CCTB Wide (Conical cemented titan base non hex and screw)	BSAD1632	30	-	10-12
SET ZZ BASE MS-S-NTB NARROW (Narrow titan base non hex and screw)	BSAP1633	30	-	10-12
SET ZZ BASE MS-S-NTB STANDARD (Narrow titan base non hex and screw)	BSAP1631	30	-	10-12
SET ZZ BASE MS-S-NTBH NARROW (Narrow titan base hex and screw)	BSAP1683	30	-	10-12
SET ZZ BASE MS-S-NTBH STANDARD (Narrow titan base hex and screw)	BSAP1681	30	-	10-12
SET ZZ BASE MS-S-PCTBH Narrow (Parallel cemented titan base hex and screw)	BSAB1683	30	-	10-12
SET ZZ BASE MS-S-PCTBH Standard (Parallel cemented titan base hex and screw)	BSAB1681	30	-	10-12
SET ZZ BASE MS-S-PCTBH Wide (Parallel cemented titan base hex and screw)	BSAB1682	30	-	10-12
SET ZZ BASE MS-S-SM Narrow (Scanmarker and screw)	ZBAD1633	-	-	10-12
SET ZZ BASE MS-S-WSM Narrow (White Scanmarker and screw)	ZBCD1633	-	-	10-12
SET ZZ BASE MS-S-WSM Standard (White Scanmarker and screw)	ZBCD1631	-	-	10-12
SET ZZ BASE MS-S-WSM Wide (White Scanmarker and screw)	ZBCD1632	-	-	10-12
ZZ BASE MS-S-MI Narrow (Model implant)	MIAG1603	-	-	10-12
ZZ BASE MS-S-MI Standard (Model implant)	MIAG1601	-	-	10-12
ZZ BASE MS-S-MI Wide (Model implant)	MIAG1602	-	-	10-12
ZZ BASE MS-S-SC Narrow (Screw)	BSAE1603	30	-	10-12
ZZ BASE MS-S-SC Standard (Screw)	BSAE1601	30	-	10-12
ZZ BASE MS-S-SC Wide (Screw)	BSAE1602	30	-	10-12
ZZ BASE MS-S-ZSC Narrow (Zirconia screw)	BSAE1633	30	-	10-12
ZZ BASE MS-S-ZSC Standard (Zirconia screw)	BSAE1631	30	-	10-12

MIS® SEVEN

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
ZZ BASE MS-S-ZSC Wide (Zirconia screw)	BSAE1632	30	-	10-12

MIS® V3

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE MS-V3-ABCH Narrow (Raw-Abutment® customizable hex and screw)	BSAF9381	30	-	10-12
SET ZZ BASE MS-V3-ABCH Standard (Raw-Abutment® customizable hex and screw)	BSAF9382	30	-	10-12
SET ZZ BASE MS-V3-CCTB Narrow (Conical cemented titan base non hex and screw)	BSAD9331	30	-	10-12
SET ZZ BASE MS-V3-CCTB Standard (Conical cemented titan base non hex and screw)	BSAD9332	30	-	10-12
SET ZZ BASE MS-V3-NTB Standard (Narrow titan base non hex and screw)	BSAP9332	30	-	10-12
SET ZZ BASE MS-V3-NTBH Standard (Narrow titan base hex and screw)	BSAP9382	30	-	10-12
SET ZZ BASE MS-V3-PCTBH Narrow (Parallel cemented titan base hex and screw)	BSAB9381	30	-	10-12
SET ZZ BASE MS-V3-PCTBH Standard (Parallel cemented titan base hex and screw)	BSAB9382	30	-	10-12
SET ZZ BASE MS-V3-SM Narrow (Scanmarker and screw)	ZBAD9331	-	-	10-12
SET ZZ BASE MS-V3-SM Standard (Scanmarker and screw)	ZBAD9332	-	-	10-12
SET ZZ BASE MS-V3-WSM Narrow (White Scanmarker and screw)	ZBCD9331	-	-	10-12
SET ZZ BASE MS-V3-WSM Standard (White Scanmarker and screw)	ZBCD9332	-	-	10-12
ZZ BASE MS-V3-MI Standard (Model implant)	MIAG9302	-	-	10-12
ZZ BASE MS-V3-SC Narrow (Screw)	BSAE9301	30	-	10-12
ZZ BASE MS-V3-SC Standard (Screw)	BSAE9302	30	-	10-12
ZZ BASE MS-V3-ZSC Narrow (Zirconia Screw)	BSAE9331	30	-	10-12
ZZ BASE MS-V3-ZSC Standard (Zirconia screw)	BSAE9332	30	-	10-12

MOZO GRAU® TAPERED SCREW

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE MG-TS-CCTB 3.4 10° (Conical cemented titan base non hex and screw)	BSAD5131	20	-	10-12
SET ZZ BASE MG-TS-CCTB 4.1 10° (Conical cemented titan base non hex and screw)	BSAD5132	20	-	10-12
SET ZZ BASE MG-TS-CCTB 4.1 30° (Conical cemented titan base non hex and screw)	BSAD5133	20	-	10-12
SET ZZ BASE MG-TS-CCTB 5.0 10° (Conical cemented titan base non hex and screw)	BSAD5134	20	-	10-12
SET ZZ BASE MG-TS-CCTB 5.0 30° (Conical cemented titan base non hex and screw)	BSAD5135	20	-	10-12
SET ZZ BASE MG-TS-PCTBH 4.1 10° (Parallel cemented titan base hex and screw)	BSAB5182	20	-	10-12
SET ZZ BASE MG-TS-PCTBH 5.0 10° (Parallel cemented titan base hex and screw)	BSAB5184	20	-	10-12
SET ZZ BASE MG-TS-SM 3.4 10° (Scanmarker and screw)	ZBAD5131	-	-	10-12
SET ZZ BASE MG-TS-SM 4.1 10° (Scanmarker and screw)	ZBAD5132	-	-	10-12
SET ZZ BASE MG-TS-SM 4.1 30° (Scanmarker and screw)	ZBAD5133	-	-	10-12
SET ZZ BASE MG-TS-SM 5.0 10° (Scanmarker and screw)	ZBAD5134	-	-	10-12
SET ZZ BASE MG-TS-SM 5.0 30° (Scanmarker and screw)	ZBAD5135	-	-	10-12
ZZ BASE MG-TS-SC 3.4 10° (Screw)	BSAE5101	20	-	10-12
ZZ BASE MG-TS-SC 4.1 10° (Screw)	BSAE5102	20	-	10-12
ZZ BASE MG-TS-SC 4.1 30° (Screw)	BSAE5103	20	-	10-12
ZZ BASE MG-TS-SC 5.0 10° (Screw)	BSAE5104	20	-	10-12
ZZ BASE MG-TS-SC 5.0 30° (Screw)	BSAE5105	20	-	10-12
ZZ BASE MG-TS-ZSC 3.4 10° (Zirconia screw)	BSAE5131	20	-	10-12
ZZ BASE MG-TS-ZSC 4.1 10° (Zirconia screw)	BSAE5132	20	-	10-12
ZZ BASE MG-TS-ZSC 4.1 30° (Zirconia screw)	BSAE5133	20	-	10-12
ZZ BASE MG-TS-ZSC 5.0 10° (Zirconia screw)	BSAE5134	20	-	10-12
ZZ BASE MG-TS-ZSC 5.0 30° (Zirconia screw)	BSAE5135	20	-	10-12

NEO BIOTECH IS

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE NT-IS-ABCH 4.3 (Raw-Abutment® customizable hex and screw)	BSAF7381	30	-	10-12
SET ZZ BASE NT-IS-CCTB 4.3 (Conical cemented titan base non hex and screw)	BSAD7331	30	-	10-12
SET ZZ BASE NT-IS-PCTBH 4.3 (Parallel cemented titan base hex and screw)	BSAB7381	30	-	10-12
SET ZZ BASE NT-IS-SM 4.3 (Scanmarker and screw)	ZBAD7331	-	-	10-12
SET ZZ BASE NT-IS-WSM 4.3 (White Scanmarker and screw)	ZBCD7331	-	-	10-12
ZZ BASE NB-IS-ZSC 4.3 (Zirconia Screw)	BSAE7331	30	-	10-12
ZZ BASE NT-IS-SC 4.3 (Screw)	BSAE7301	30	-	10-12

NEOSS PROACTIVE®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE NE-PA-ABCH-ID 3.25 narrow (Raw-Abutment® customizable hex ID and screw)	BSAF3991	32	-	10-12
SET ZZ BASE NE-PA-ABCH-ID 3.5/4.0/4.5/5.0/5.5 tapered (Raw-Abutment® customizable hex ID and screw)	BSAF3992	32	-	10-12
SET ZZ BASE NE-PA-CCTB 3.25 narrow (Conical cemented titan base non hex and screw)	BSAD3931	32	-	10-12
SET ZZ BASE NE-PA-CCTB 3.5/4.0/4.5/5.0/5.5 tapered (Conical cemented titan base non hex and screw)	BSAD3932	32	-	10-12
SET ZZ BASE NE-PA-NTB 3.25 narrow (Narrow titan base non hex and screw)	BSAP3931	32	-	10-12
SET ZZ BASE NE-PA-NTBH 3.25 narrow (Narrow titan base hex and screw)	BSAP3981	32	-	10-12
SET ZZ BASE NE-PA-PCTBH 3.25 narrow (Parallel cemented titan base hex and screw)	BSAB3981	32	-	10-12
SET ZZ BASE NE-PA-PCTBH 3.5/4.0/4.5/5.0/5.5 tapered (Parallel cemented titan base hex and screw)	BSAB3982	32	-	10-12
SET ZZ BASE NE-PA-SM 3.25 narrow (Scanmarker and screw)	ZBAD3931	-	-	10-12
SET ZZ BASE NE-PA-SM 3.5/4.0/4.5/5.0/5.5 tapered (Scanmarker and screw)	ZBAD3932	-	-	10-12
ZZ BASE NE-PA-SC 3.25 narrow (Screw)	BSAE3901	32	-	10-12
ZZ BASE NE-PA-SC 3.5/4.0/4.5/5.0/5.5 tapered (Screw)	BSAE3902	32	-	10-12
ZZ BASE NE-PA-ZSC 3.25 narrow (Zirconia screw)	BSAE3931	32	-	10-12
ZZ BASE NE-PA-ZSC 3.5/4.0/4.5/5.0/5.5 tapered (Zirconia screw)	BSAE3932	32	-	10-12

NOBEL BIOCARE® BRÄNEMARK

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE NB-B-ABCH 3.5 (Raw-Abutment® customizable hex and screw)	BSAF0291	25	-	10-12
SET ZZ BASE NB-B-ABCH-ID 4.1 (Raw-Abutment® customizable hex ID and screw)	BSAF0292	35	-	10-12
SET ZZ BASE NB-B-ABCH-ID 5.1 (Raw-Abutment® customizable hex ID and screw)	BSAF0293	35	-	10-12
SET ZZ BASE NB-B-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD0231	25	-	10-12
SET ZZ BASE NB-B-CCTB 4.1 (Conical cemented titan base non hex and screw)	BSAD0232	35	-	10-12
SET ZZ BASE NB-B-CCTB 5.1 (Conical cemented titan base non hex and screw)	BSAD0233	35	-	10-12
SET ZZ BASE NB-B-CCTB-G 3.5 (Conical cemented titan base non hex and screw Gold)	BSCD0241	25	-	10-12
SET ZZ BASE NB-B-CCTB-G 4.1 (Conical cemented titan base non hex and screw Gold)	BSCD0242	35	-	10-12
SET ZZ BASE NB-B-CCTB-G 5.1 (Conical cemented titan base non hex and screw Gold)	BSCD0243	35	-	10-12
SET ZZ BASE NB-B-CCTB-L10 3.5 (Conical cemented titan base non hex L10 and screw)	BSAD0231L10	25	-	10-12
SET ZZ BASE NB-B-CCTB-L10 4.1 (Conical cemented titan base non hex L10 and screw)	BSAD0232L10	35	-	10-12
SET ZZ BASE NB-B-CCTB-L10 5.1 (Conical cemented titan base non hex L10 and screw)	BSAD0233L10	35	-	10-12
SET ZZ BASE NB-B-CCTB-L20 3.5 (Conical cemented titan base non hex L20 and screw)	BSAD0231L20	25	-	10-12
SET ZZ BASE NB-B-CCTB-L20 4.1 (Conical cemented titan base non hex L20 and screw)	BSAD0232L20	35	-	10-12
SET ZZ BASE NB-B-CCTB-L20 5.1 (Conical cemented titan base non hex L20 and screw)	BSAD0233L20	35	-	10-12
SET ZZ BASE NB-B-CCTB-L30 3.5 (Conical cemented titan base non hex L30 and screw)	BSAD0231L30	25	-	10-12
SET ZZ BASE NB-B-CCTB-L30 4.1 (Conical cemented titan base non hex L30 and screw)	BSAD0232L30	35	-	10-12
SET ZZ BASE NB-B-CCTB-L30 5.1 (Conical cemented titan base non hex L30 and screw)	BSAD0233L30	35	-	10-12
SET ZZ BASE NB-B-NTB 3.5 (Narrow titan base non hex and screw)	BSAP0231	25	-	10-12
SET ZZ BASE NB-B-NTB 4.1 (Narrow titan base non hex and screw)	BSAP0232	35	-	10-12
SET ZZ BASE NB-B-NTBH 3.5 (Narrow titan base hex and screw)	BSAP0281	25	-	10-12
SET ZZ BASE NB-B-NTBH 4.1 (Narrow titan base hex and screw)	BSAP0282	35	-	10-12
SET ZZ BASE NB-B-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB0281	25	-	10-12
SET ZZ BASE NB-B-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB0291	25	-	10-12
SET ZZ BASE NB-B-PCTBH 4.1 (Parallel cemented titan base hex and screw)	BSAB0282	35	-	10-12
SET ZZ BASE NB-B-PCTBH 4.1 (Parallel cemented titan base hex and screw)	BSAB0292	35	-	10-12
SET ZZ BASE NB-B-PCTBH 5.1 (Parallel cemented titan base hex and screw)	BSAB0283	35	-	10-12
SET ZZ BASE NB-B-PCTBH 5.1 (Parallel cemented titan base hex and screw)	BSAB0293	35	-	10-12
SET ZZ BASE NB-B-PCTBH-G 3.5 (Parallel cemented titan base hex and screw Gold)	BSCB0291	25	-	10-12
SET ZZ BASE NB-B-PCTBH-G 4.1 (Parallel cemented titan base hex and screw Gold)	BSCB0292	35	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

NOBEL BIOCARE® BRÄNEMARK

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE NB-B-PCTBH-G 5.1 (Parallel cemented titan base hex and screw Gold)	BSCB0293	35	-	10-12
SET ZZ BASE NB-B-PCTBH-L10 3.5 (Parallel cemented titan base hex L10 and screw)	BSAB0291L10	25	-	10-12
SET ZZ BASE NB-B-PCTBH-L10 4.1 (Parallel cemented titan base hex L10 and screw)	BSAB0292L10	35	-	10-12
SET ZZ BASE NB-B-PCTBH-L10 5.1 (Parallel cemented titan base hex L10 and screw)	BSAB0293L10	35	-	10-12
SET ZZ BASE NB-B-PCTBH-L20 3.5 (Parallel cemented titan base hex L20 and screw)	BSAB0291L20	25	-	10-12
SET ZZ BASE NB-B-PCTBH-L20 4.1 (Parallel cemented titan base hex L20 and screw)	BSAB0292L20	35	-	10-12
SET ZZ BASE NB-B-PCTBH-L20 5.1 (Parallel cemented titan base hex L20 and screw)	BSAB0293L20	35	-	10-12
SET ZZ BASE NB-B-PCTBH-L30 3.5 (Parallel cemented titan base hex L30 and screw)	BSAB0291L30	25	-	10-12
SET ZZ BASE NB-B-PCTBH-L30 4.1 (Parallel cemented titan base hex L30 and screw)	BSAB0292L30	35	-	10-12
SET ZZ BASE NB-B-PCTBH-L30 5.1 (Parallel cemented titan base hex L30 and screw)	BSAB0293L30	35	-	10-12
SET ZZ BASE NB-B-PCTBH-L40 3.5 (Parallel cemented titan base hex L40 and screw)	BSAB0291L40	25	-	10-12
SET ZZ BASE NB-B-PCTBH-L40 4.1 (Parallel cemented titan base hex L40 and screw)	BSAB0292L40	35	-	10-12
SET ZZ BASE NB-B-SA 3.5 (Scanalog and screw)	ZBDE0231	15	-	10-12
SET ZZ BASE NB-B-SA 4.1 (Scanalog and screw)	ZBDE0232	15	-	10-12
SET ZZ BASE NB-B-SA 5.1 (Scanalog and screw)	ZBDE0233	15	-	10-12
SET ZZ BASE NB-B-SM 3.5 (Scanmarker and screw)	ZBAD0231	-	-	10-12
SET ZZ BASE NB-B-SM 4.1 (Scanmarker and screw)	ZBAD0232	-	-	10-12
SET ZZ BASE NB-B-SM 5.1 (Scanmarker and screw)	ZBAD0233	-	-	10-12
SET ZZ BASE NB-B-WSM 3.5 (White Scanmarker and screw)	ZBCD0231	-	-	10-12
SET ZZ BASE NB-B-WSM 4.1 (White Scanmarker and screw)	ZBCD0232	12	-	10-12
SET ZZ BASE NB-B-WSM 5.1 (White Scanmarker and screw)	ZBCD0233	-	-	10-12
SET ZZ BASE NB-B-ZZMUA-L15 3.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0231L15	25	25	10-12
SET ZZ BASE NB-B-ZZMUA-L15 4.1 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0232L15	25	35	10-12
SET ZZ BASE NB-B-ZZMUA-L15 5.1 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0233L15	25	35	10-12
SET ZZ BASE NB-B-ZZMUA-L25 3.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0231L25	25	25	10-12
SET ZZ BASE NB-B-ZZMUA-L25 4.1 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0232L25	25	35	10-12
SET ZZ BASE NB-B-ZZMUA-L25 5.1 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0233L25	25	35	10-12
SET ZZ BASE NB-B-ZZMUA-L35 3.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0231L35	25	25	10-12
SET ZZ BASE NB-B-ZZMUA-L35 4.1 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0232L35	25	35	10-12
SET ZZ BASE NB-B-ZZMUA-L35 5.1 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0233L35	25	35	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

NOBEL BIOCARE® BRÄNEMARK

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE NB-B-ZZMUA-L45 3.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0231L45	25	25	10-12
SET ZZ BASE NB-B-ZZMUA-L45 4.1 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0232L45	25	35	10-12
SET ZZ BASE NB-B-ZZMUA-L45 5.1 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0233L45	25	35	10-12
ZZ BASE NB-B-MI 3.5 (Model implant)	MIAG0201	-	-	10-12
ZZ BASE NB-B-MI 4.1 (Model implant)	MIAG0202	-	-	10-12
ZZ BASE NB-B-MI 5.1 (Model implant)	MIAG0203	-	-	10-12
ZZ BASE NB-B-SC 3.5 (Screw)	BSAE0201	25	-	10-12
ZZ BASE NB-B-SC 4.1 (Screw)	BSAE0202	35	-	10-12
ZZ BASE NB-B-SC 5.1 (Screw)	BSAE0203	35	-	10-12
ZZ BASE NB-B-SC-G 3.5 (Screw Gold)	BSCE0201	25	-	10-12
ZZ BASE NB-B-SC-G 4.1 (Screw Gold)	BSCE0202	35	-	10-12
ZZ BASE NB-B-SC-G 5.1 (Screw Gold)	BSCE0203	35	-	10-12
ZZ BASE NB-B-ZSC 3.5 (Zirconia screw)	BSAE0231	25	-	10-12
ZZ BASE NB-B-ZSC 4.1 (Zirconia screw)	BSAE0232	35	-	10-12
ZZ BASE NB-B-ZSC 5.1 (Zirconia screw)	BSAE0233	35	-	10-12
ZZ BASE NB-B-ZSC-G 3.5 (Zirconia screw Gold)	BSCE0231	25	-	10-12
ZZ BASE NB-B-ZSC-G 4.1 (Zirconia screw Gold)	BSCE0232	35	-	10-12
ZZ BASE NB-B-ZSC-G 5.1 (Zirconia screw Gold)	BSCE0233	35	-	10-12

NOBEL BIO CARE® MULTI UNIT ABUTMENT

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE NB-MUA-CCTB NP/RP (Conical cemented titan base non hex and screw)	BSAD0441	15	-	10-12
SET ZZ BASE NB-MUA-CCTB NP/RP (Conical cemented titan base non hex and screw)	BSAD0431	15	-	10-12
SET ZZ BASE NB-MUA-CCTB WP (Conical cemented titan base non hex and screw)	BSAD0432	35	-	10-12
SET ZZ BASE NB-MUA-SA NP/RP (Scanalog and screw)	ZBDE0431	15	-	10-12
SET ZZ BASE NB-MUA-SA WP (Scanalog and screw)	ZBDE0432	15	-	10-12
SET ZZ BASE NB-MUA-SM NP/RP (Scanmarker and screw)	ZBAD0431	-	-	10-12
SET ZZ BASE NB-MUA-SM WP (Scanmarker and screw)	ZBAD0432	-	-	10-12
SET ZZ BASE NB-MUA-WSM NP/RP (White Scanmarker and screw)	ZBCD0431	-	-	10-12
SET ZZ BASE NB-MUA-WSM WP (White Scanmarker and screw)	ZBCD0432	-	-	10-12
ZZ BASE NB-MUA-MI NP/RP (Model implant)	MIAG0401	-	-	10-12
ZZ BASE NB-MUA-MI WP (Model implant)	MIAG0402	-	-	10-12
ZZ BASE NB-MUA-SC NP/RP (Screw)	BSAE0401	15	-	10-12
ZZ BASE NB-MUA-SC WP (Screw)	BSAE0402	35	-	10-12
ZZ BASE NB-MUA-SC-G NP/RP (Screw Gold)	BSCE0401	15	-	10-12
ZZ BASE NB-MUA-SC-G WP (Screw Gold)	BSCE0402	35	-	10-12
ZZ BASE NB-MUA-ZSC NP/RP (Zirconia screw)	BSAE0431	15	-	10-12
ZZ BASE NB-MUA-ZSC WP (Zirconia screw)	BSAE0432	35	-	10-12

NOBEL BIO CARE® NOBELACTIVE

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE NB-A-ABCH CC 3.0 (Raw-Abutment® customizable hex and screw)	BSAF2483	15	-	10-12
SET ZZ BASE NB-A-ABCH NP (Raw-Abutment® customizable hex and screw)	BSAF2481	25	-	10-12
SET ZZ BASE NB-A-ABCH RP (Raw-Abutment® customizable hex and screw)	BSAF2482	35	-	10-12
SET ZZ BASE NB-A-ABCH WP (Raw-Abutment® customizable hex and screw)	BSAF2484	35	-	10-12
SET ZZ BASE NB-A-ABCH-D14 RP (Raw-Abutment® customizable hex D14 and screw)	BSAF2482D14	35	-	10-12
SET ZZ BASE NB-A-ABCH-D14 WP (Raw-Abutment® customizable hex D14 and screw)	BSAF2484D14	35	-	10-12
SET ZZ BASE NB-A-ABH-L15 NP (Abutment hex L15 and screw)	BSAC2431L15	25	-	10-12
SET ZZ BASE NB-A-ABH-L15 RP (Abutment hex L15 and screw)	BSAC2432L15	35	-	10-12
SET ZZ BASE NB-A-ABH-L15 WP (Abutment hex L15 and screw)	BSAC2434L15	35	-	10-12
SET ZZ BASE NB-A-ABH-L30 NP (Abutment hex L30 and screw)	BSAC2431L30	25	-	10-12
SET ZZ BASE NB-A-ABH-L30 RP (Abutment hex L30 and screw)	BSAC2432L30	35	-	10-12
SET ZZ BASE NB-A-ABH-L30 WP (Abutment hex L30 and screw)	BSAC2434L30	35	-	10-12
SET ZZ BASE NB-A-CCTB CC 3.0 (Conical cemented titan base non hex and screw)	BSAD2433	15	-	10-12
SET ZZ BASE NB-A-CCTB NP (Conical cemented titan base non hex and screw)	BSAD2431	25	-	10-12
SET ZZ BASE NB-A-CCTB RP (Conical cemented titan base non hex and screw)	BSAD2432	35	-	10-12
SET ZZ BASE NB-A-CCTB WP (Conical cemented titan base non hex and screw)	BSAD2434	35	-	10-12
SET ZZ BASE NB-A-CCTB-L20 CC 3.0 (Conical cemented titan base non hex L20 and screw)	BSAD2433L20	15	-	10-12
SET ZZ BASE NB-A-CCTB-L20 NP (Conical cemented titan base non hex L20 and screw)	BSAD2431L20	25	-	10-12
SET ZZ BASE NB-A-CCTB-L20 RP (Conical cemented titan base non hex L20 and screw)	BSAD2432L20	35	-	10-12
SET ZZ BASE NB-A-CCTB-L30 CC 3.0 (Conical cemented titan base non hex L30 and screw)	BSAD2433L30	15	-	10-12
SET ZZ BASE NB-A-CCTB-L30 NP (Conical cemented titan base non hex L30 and screw)	BSAD2431L30	25	-	10-12
SET ZZ BASE NB-A-CCTB-L30 RP (Conical cemented titan base non hex L30 and screw)	BSAD2432L30	35	-	10-12
SET ZZ BASE NB-A-CCTB-L40 NP (Conical cemented titan base non hex L40 and screw)	BSAD2431L40	25	-	10-12
SET ZZ BASE NB-A-CCTB-L40 RP (Conical cemented titan base non hex L40 and screw)	BSAD2432L40	35	-	10-12
SET ZZ BASE NB-A-NTB CC 3.0 (Narrow titan base non hex and screw)	BSAP2433	15	-	10-12
SET ZZ BASE NB-A-NTB NP (Narrow titan base non hex and screw)	BSAP2431	25	-	10-12
SET ZZ BASE NB-A-NTB RP (Narrow titan base non hex and screw)	BSAP2432	35	-	10-12
SET ZZ BASE NB-A-NTBH CC 3.0 (Narrow titan base hex and screw)	BSAP2483	15	-	10-12
SET ZZ BASE NB-A-NTBH NP (Narrow titan base hex and screw)	BSAP2481	25	-	10-12
SET ZZ BASE NB-A-NTBH RP (Narrow titan base hex and screw)	BSAP2482	35	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

NOBEL BIO CARE® NOBELACTIVE

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE NB-A-PCTBH CC 3.0 (Parallel cemented titan base hex and screw)	BSAB2483	15	-	10-12
SET ZZ BASE NB-A-PCTBH NP (Parallel cemented titan base hex and screw)	BSAB2481	25	-	10-12
SET ZZ BASE NB-A-PCTBH RP (Parallel cemented titan base hex and screw)	BSAB2482	35	-	10-12
SET ZZ BASE NB-A-PCTBH WP (Parallel cemented titan base hex and screw)	BSAB2484	35	-	10-12
SET ZZ BASE NB-A-PCTBH-L20 CC 3.0 (Parallel cemented titan base hex L20 and screw)	BSAB2493L20	15	-	10-12
SET ZZ BASE NB-A-PCTBH-L20 NP (Parallel cemented titan base hex L20 and screw)	BSAB2491L20	25	-	10-12
SET ZZ BASE NB-A-PCTBH-L20 RP (Parallel cemented titan base hex L20 and screw)	BSAB2492L20	35	-	10-12
SET ZZ BASE NB-A-PCTBH-L30 CC 3.0 (Parallel cemented titan base hex L30 and screw)	BSAB2493L30	15	-	10-12
SET ZZ BASE NB-A-PCTBH-L30 NP (Parallel cemented titan base hex L30 and screw)	BSAB2491L30	25	-	10-12
SET ZZ BASE NB-A-PCTBH-L30 RP (Parallel cemented titan base hex L30 and screw)	BSAB2492L30	35	-	10-12
SET ZZ BASE NB-A-PCTBH-L40 NP (Parallel cemented titan base hex L40 and screw)	BSAB2491L40	25	-	10-12
SET ZZ BASE NB-A-PCTBH-L40 RP (Parallel cemented titan base hex L40 and screw)	BSAB2492L40	35	-	10-12
SET ZZ BASE NB-A-SA CC 3.0 (Scanalog and screw)	ZBDE2433	15	-	10-12
SET ZZ BASE NB-A-SA NP (Scanalog and screw)	ZBDE2431	15	-	10-12
SET ZZ BASE NB-A-SA RP (Scanalog and screw)	ZBDE2432	15	-	10-12
SET ZZ BASE NB-A-SA WP (Scanalog and screw)	ZBDE2434	15	-	10-12
SET ZZ BASE NB-A-SM CC 3.0 (Scanmarker and screw)	ZBAD2433	-	-	10-12
SET ZZ BASE NB-A-SM NP (Scanmarker and screw)	ZBAD2431	-	-	10-12
SET ZZ BASE NB-A-SM RP (Scanmarker and screw)	ZBAD2432	-	-	10-12
SET ZZ BASE NB-A-SM WP (Scanmarker and screw)	ZBAD2434	-	-	10-12
SET ZZ BASE NB-A-WSM CC 3.0 (White Scanmarker and screw)	ZBCD2433	-	-	10-12
SET ZZ BASE NB-A-WSM NP (White Scanmarker and screw)	ZBCD2431	-	-	10-12
SET ZZ BASE NB-A-WSM RP (White Scanmarker and screw)	ZBCD2432	-	-	10-12
SET ZZ BASE NB-A-WSM WP (White Scanmarker and screw)	ZBCD2434	-	-	10-12
SET ZZ BASE NB-A-ZZMUA-L15 CC 3.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH2433L15	25	15	10-12
SET ZZ BASE NB-A-ZZMUA-L15 NP (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH2431L15	25	25	10-12
SET ZZ BASE NB-A-ZZMUA-L15 RP (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH2432L15	25	35	10-12
SET ZZ BASE NB-A-ZZMUA-L15 WP (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH2434L15	25	35	10-12
SET ZZ BASE NB-A-ZZMUA-L25 CC 3.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH2433L25	25	15	10-12
SET ZZ BASE NB-A-ZZMUA-L25 NP (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH2431L25	25	25	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

NOBEL BIO CARE® NOBELACTIVE

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE NB-A-ZZMUA-L25 RP (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH2432L25	25	35	10-12
SET ZZ BASE NB-A-ZZMUA-L25 WP (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH2434L25	25	35	10-12
SET ZZ BASE NB-A-ZZMUA-L35 CC 3.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH2433L35	25	15	10-12
SET ZZ BASE NB-A-ZZMUA-L35 NP (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH2431L35	25	25	10-12
SET ZZ BASE NB-A-ZZMUA-L35 RP (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH2432L35	25	35	10-12
SET ZZ BASE NB-A-ZZMUA-L35 WP (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH2434L35	25	35	10-12
SET ZZ BASE NB-A-ZZMUA-L45 CC 3.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH2433L45	25	15	10-12
SET ZZ BASE NB-A-ZZMUA-L45 NP (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH2431L45	25	25	10-12
SET ZZ BASE NB-A-ZZMUA-L45 RP (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH2432L45	25	35	10-12
SET ZZ BASE NB-A-ZZMUA-L45 WP (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH2434L45	25	35	10-12
ZZ BASE NB-A-MI CC 3.0 (Model implant)	MIAG2403	-	-	10-12
ZZ BASE NB-A-MI NP (Model implant)	MIAG2401	-	-	10-12
ZZ BASE NB-A-MI RP (Model implant)	MIAG2402	-	-	10-12
ZZ BASE NB-A-MI WP (Model implant)	MIAG2404	-	-	10-12
ZZ BASE NB-A-SC CC 3.0 (Screw)	BSAE2403	15	-	10-12
ZZ BASE NB-A-SC NP (Screw)	BSAE2401	25	-	10-12
ZZ BASE NB-A-SC RP (Screw)	BSAE2402	35	-	10-12
ZZ BASE NB-A-SC WP (Screw)	BSAE2404	35	-	10-12
ZZ BASE NB-A-ZSC CC 3.0 (Zirconia screw)	BSAE2433	15	-	10-12
ZZ BASE NB-A-ZSC NP (Zirconia screw)	BSAE2431	25	-	10-12
ZZ BASE NB-A-ZSC RP (Zirconia screw)	BSAE2432	35	-	10-12
ZZ BASE NB-A-ZSC WP (Zirconia screw)	BSAE2434	35	-	10-12

NOBEL BIO CARE® NOBELREPLACE

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE NB-R-ABCH-D14-ID 3.5 (Raw-Abutment® customizable hex D14 ID and screw)	BSAF0391D14	35	-	10-12
SET ZZ BASE NB-R-ABCH-D14-ID 4.3 (Raw-Abutment® customizable hex D14 ID and screw)	BSAF0392D14	-	-	10-12
SET ZZ BASE NB-R-ABCH-D14-ID 5.0 (Raw-Abutment® customizable hex D14 ID and screw)	BSAF0393D14	-	-	10-12
SET ZZ BASE NB-R-ABCH-D14-ID 6.0 (Raw-Abutment® customizable hex D14 ID and screw)	BSAF0394D14	-	-	10-12
SET ZZ BASE NB-R-ABCH-ID 3.5 (Raw-Abutment® customizable hex ID and screw)	BSAF0391	35	-	10-12
SET ZZ BASE NB-R-ABCH-ID 4.3 (Raw-Abutment® customizable hex ID and screw)	BSAF0392	-	-	10-12
SET ZZ BASE NB-R-ABCH-ID 5.0 (Raw-Abutment® customizable hex ID and screw)	BSAF0393	-	-	10-12
SET ZZ BASE NB-R-ABCH-ID 6.0 (Raw-Abutment® customizable hex ID and screw)	BSAF0394	-	-	10-12
SET ZZ BASE NB-R-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD0331	35	-	10-12
SET ZZ BASE NB-R-CCTB 4.3 (Conical cemented titan base non hex and screw)	BSAD0332	-	-	10-12
SET ZZ BASE NB-R-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD0333	-	-	10-12
SET ZZ BASE NB-R-CCTB 6.0 (Conical cemented titan base non hex and screw)	BSAD0334	-	-	10-12
SET ZZ BASE NB-R-CCTB-L10 3.5 (Conical cemented titan base non hex L10 and screw)	BSAD0331L10	35	-	10-12
SET ZZ BASE NB-R-CCTB-L10 4.3 (Conical cemented titan base non hex L10 and screw)	BSAD0332L10	-	-	10-12
SET ZZ BASE NB-R-CCTB-L10 5.0 (Conical cemented titan base non hex L10 and screw)	BSAD0333L10	-	-	10-12
SET ZZ BASE NB-R-CCTB-L10 6.0 (Conical cemented titan base non hex L10 and screw)	BSAD0334L10	-	-	10-12
SET ZZ BASE NB-R-CCTB-L20 3.5 (Conical cemented titan base non hex L20 and screw)	BSAD0331L20	35	-	10-12
SET ZZ BASE NB-R-CCTB-L20 4.3 (Conical cemented titan base non hex L20 and screw)	BSAD0332L20	-	-	10-12
SET ZZ BASE NB-R-CCTB-L20 5.0 (Conical cemented titan base non hex L20 and screw)	BSAD0333L20	-	-	10-12
SET ZZ BASE NB-R-CCTB-L20 6.0 (Conical cemented titan base non hex L20 and screw)	BSAD0334L20	-	-	10-12
SET ZZ BASE NB-R-CCTB-L30 3.5 (Conical cemented titan base non hex L30 and screw)	BSAD0331L30	35	-	10-12
SET ZZ BASE NB-R-CCTB-L30 4.3 (Conical cemented titan base non hex L30 and screw)	BSAD0332L30	-	-	10-12
SET ZZ BASE NB-R-CCTB-L30 5.0 (Conical cemented titan base non hex L30 and screw)	BSAD0333L30	-	-	10-12
SET ZZ BASE NB-R-CCTB-L30 6.0 (Conical cemented titan base non hex L30 and screw)	BSAD0334L30	-	-	10-12
SET ZZ BASE NB-R-NTB 3.5 (Narrow titan base non hex and screw)	BSAP0331	35	-	10-12
SET ZZ BASE NB-R-NTB 4.3 (Narrow titan base non hex and screw)	BSAP0332	-	-	10-12
SET ZZ BASE NB-R-NTBH 3.5 (Narrow titan base hex and screw)	BSAP0381	35	-	10-12
SET ZZ BASE NB-R-NTBH 4.3 (Narrow titan base hex and screw)	BSAP0382	-	-	10-12
SET ZZ BASE NB-R-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB0381	35	-	10-12
SET ZZ BASE NB-R-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB0391	35	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

NOBEL BIO CARE® NOBELREPLACE

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE NB-R-PCTBH 4.3 (Parallel cemented titan base hex and screw)	BSAB0392	-	-	10-12
SET ZZ BASE NB-R-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB0383	-	-	10-12
SET ZZ BASE NB-R-PCTBH 6.0 (Parallel cemented titan base hex and screw)	BSAB0384	-	-	10-12
SET ZZ BASE NB-R-PCTBH-L10 3.5 (Parallel cemented titan base hex L10 and screw)	BSAB0391L10	35	-	10-12
SET ZZ BASE NB-R-PCTBH-L10 4.3 (Parallel cemented titan base hex L10 and screw)	BSAB0392L10	-	-	10-12
SET ZZ BASE NB-R-PCTBH-L10 5.0 (Parallel cemented titan base hex L10 and screw)	BSAB0393L10	-	-	10-12
SET ZZ BASE NB-R-PCTBH-L10 6.0 (Parallel cemented titan base hex L10 and screw)	BSAB0394L10	-	-	10-12
SET ZZ BASE NB-R-PCTBH-L20 3.5 (Parallel cemented titan base hex L20 and screw)	BSAB0391L20	35	-	10-12
SET ZZ BASE NB-R-PCTBH-L20 4.3 (Parallel cemented titan base hex L20 and screw)	BSAB0392L20	-	-	10-12
SET ZZ BASE NB-R-PCTBH-L20 5.0 (Parallel cemented titan base hex L20 and screw)	BSAB0393L20	-	-	10-12
SET ZZ BASE NB-R-PCTBH-L20 6.0 (Parallel cemented titan base hex L20 and screw)	BSAB0394L20	-	-	10-12
SET ZZ BASE NB-R-PCTBH-L30 3.5 (Parallel cemented titan base hex L30 and screw)	BSAB0391L30	35	-	10-12
SET ZZ BASE NB-R-PCTBH-L30 4.3 (Parallel cemented titan base hex L30 and screw)	BSAB0392L30	-	-	10-12
SET ZZ BASE NB-R-PCTBH-L30 5.0 (Parallel cemented titan base hex L30 and screw)	BSAB0393L30	-	-	10-12
SET ZZ BASE NB-R-PCTBH-L30 6.0 (Parallel cemented titan base hex L30 and screw)	BSAB0394L30	-	-	10-12
SET ZZ BASE NB-R-PCTBH-L40 3.5 (Parallel cemented titan base hex L40 and screw)	BSAB0391L40	35	-	10-12
SET ZZ BASE NB-R-PCTBH-L40 4.3 (Parallel cemented titan base hex L40 and screw)	BSAB0392L40	-	-	10-12
SET ZZ BASE NB-R-PCTBH-L40 5.0 (Parallel cemented titan base hex L40 and screw)	BSAB0393L40	-	-	10-12
SET ZZ BASE NB-R-SA 3.5 (Scanalog and screw)	ZBDE0331	15	-	10-12
SET ZZ BASE NB-R-SA 4.3 (Scanalog and screw)	ZBDE0332	15	-	10-12
SET ZZ BASE NB-R-SA 5.0 (Scanalog and screw)	ZBDE0333	15	-	10-12
SET ZZ BASE NB-R-SA 6.0 (Scanalog and screw)	ZBDE0334	15	-	10-12
SET ZZ BASE NB-R-SM 3.5 (Scanmarker and screw)	ZBAD0331	-	-	10-12
SET ZZ BASE NB-R-SM 4.3 (Scanmarker and screw)	ZBAD0332	-	-	10-12
SET ZZ BASE NB-R-SM 5.0 (Scanmarker and screw)	ZBAD0333	-	-	10-12
SET ZZ BASE NB-R-SM 6.0 (Scanmarker and screw)	ZBAD0334	-	-	10-12
SET ZZ BASE NB-R-WSM 3.5 (White Scanmarker and screw)	ZBCD0331	-	-	10-12
SET ZZ BASE NB-R-WSM 4.3 (White Scanmarker and screw)	ZBCD0332	-	-	10-12
SET ZZ BASE NB-R-WSM 5.0 (White Scanmarker and screw)	ZBCD0333	-	-	10-12
SET ZZ BASE NB-R-WSM 6.0 (White Scanmarker and screw)	ZBCD0334	-	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

NOBEL BIO CARE® NOBELREPLACE

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE NB-R-ZZMUA-L05 4.3 (Zirkonzahn Multi Unit Abutment non hex L05 and Screw)	BSAH0332L05	25	35	10-12
SET ZZ BASE NB-R-ZZMUA-L05 5.0 (Zirkonzahn Multi Unit Abutment non hex L05 and Screw)	BSAH0333L05	25	35	10-12
SET ZZ BASE NB-R-ZZMUA-L15 3.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0331L15	25	35	10-12
SET ZZ BASE NB-R-ZZMUA-L15 4.3 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0332L15	25	35	10-12
SET ZZ BASE NB-R-ZZMUA-L15 5.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0333L15	25	35	10-12
SET ZZ BASE NB-R-ZZMUA-L15 6.0 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0334L15	25	35	10-12
SET ZZ BASE NB-R-ZZMUA-L25 3.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0331L25	25	35	10-12
SET ZZ BASE NB-R-ZZMUA-L25 4.3 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0332L25	25	35	10-12
SET ZZ BASE NB-R-ZZMUA-L25 5.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0333L25	25	35	10-12
SET ZZ BASE NB-R-ZZMUA-L25 6.0 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0334L25	25	35	10-12
SET ZZ BASE NB-R-ZZMUA-L35 3.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0331L35	25	35	10-12
SET ZZ BASE NB-R-ZZMUA-L35 4.3 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0332L35	25	35	10-12
SET ZZ BASE NB-R-ZZMUA-L35 5.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0333L35	25	35	10-12
SET ZZ BASE NB-R-ZZMUA-L35 6.0 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0334L35	25	35	10-12
SET ZZ BASE NB-R-ZZMUA-L45 3.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0331L45	25	35	10-12
SET ZZ BASE NB-R-ZZMUA-L45 4.3 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0332L45	25	35	10-12
SET ZZ BASE NB-R-ZZMUA-L45 5.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0333L45	25	35	10-12
SET ZZ BASE NB-R-ZZMUA-L45 6.0 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0334L45	25	35	10-12
ZZ BASE NB-R-MI 3.5 (Model implant)	MIAG0301	-	-	10-12
ZZ BASE NB-R-MI 4.3 (Model implant)	MIAG0302	-	-	10-12
ZZ BASE NB-R-MI 5.0 (Model implant)	MIAG0303	-	-	10-12
ZZ BASE NB-R-MI 6.0 (Model implant)	MIAG0304	-	-	10-12
ZZ BASE NB-R-SC 3.5 (Screw)	BSAE0301	35	-	10-12
ZZ BASE NB-R-SC 4.3 (Screw)	BSAE0302	-	-	10-12
ZZ BASE NB-R-SC 5.0 (Screw)	BSAE0303	-	-	10-12
ZZ BASE NB-R-SC 6.0 (Screw)	BSAE0304	-	-	10-12
ZZ BASE NB-R-SC-G 3.5 (Screw Gold)	BSCE0301	35	-	10-12
ZZ BASE NB-R-SC-G 4.3 (Screw Gold)	BSCE0302	-	-	10-12
ZZ BASE NB-R-SC-G 5.0 (Screw Gold)	BSCE0303	-	-	10-12
ZZ BASE NB-R-SC-G 6.0 (Screw Gold)	BSCE0304	-	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

NOBEL BIO CARE® NOBELREPLACE

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
ZZ BASE NB-R-ZSC 3.5 (Zirconia screw)	BSAE0331	35	-	10-12
ZZ BASE NB-R-ZSC 4.3 (Zirconia screw)	BSAE0332	-	-	10-12
ZZ BASE NB-R-ZSC 5.0 (Zirconia screw)	BSAE0333	-	-	10-12
ZZ BASE NB-R-ZSC 6.0 (Zirconia screw)	BSAE0334	-	-	10-12

NOBEL BIO CARE® NOBELZYGOMA

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE NB-Z-ABCH-ID B4.1 (Raw-Abutment® customizable hex ID and screw)	BSAF9991	35	-	10-12
SET ZZ BASE NB-Z-CCTB B4.1 (Conical cemented titan base non hex and screw)	BSAD9931	35	-	10-12
SET ZZ BASE NB-Z-NTB B4.1 (Narrow titan base non hex and screw)	BSAP9931	35	-	10-12
SET ZZ BASE NB-Z-NTBH B4.1 (Narrow titan base hex and screw)	BSAP9981	35	-	10-12
SET ZZ BASE NB-Z-PCTBH B4.1 (Parallel cemented titan base hex and screw)	BSAB9981	35	-	10-12
SET ZZ BASE NB-Z-SM B4.1 (Scanmarker and screw)	ZBAD9931	-	-	10-12
SET ZZ BASE NB-Z-WSM B4.1 (White Scanmarker and screw)	ZBCD9931	-	-	10-12
ZZ BASE NB-Z-MI B4.1 (Model implant)	MIAG9901	-	-	10-12
ZZ BASE NB-Z-SC B4.1 (Screw)	BSAE9901	35	-	10-12
ZZ BASE NB-Z-ZSC B4.1 (Zirconia Screw)	BSAE9931	35	-	10-12

OSSTEM IMPLANT CONVERTIBLE ABUTMENT

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE OS-GS-CA-CCTB 4.0 (Conical cemented titan base non hex and screw)	BSAD5531	30	-	10-12
SET ZZ BASE OS-GS-CA-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD5532	30	-	10-12
SET ZZ BASE OS-GS-CA-CCTB 6.0 (Conical cemented titan base non hex and screw)	BSAD5533	30	-	10-12
SET ZZ BASE OS-GS-CA-PCTBH 4.0 (Parallel cemented titan base hex and screw)	BSAB5581	30	-	10-12
SET ZZ BASE OS-GS-CA-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB5582	30	-	10-12
SET ZZ BASE OS-GS-CA-PCTBH 6.0 (Parallel cemented titan base hex and screw)	BSAB5583	30	-	10-12
SET ZZ BASE OS-GS-CA-SM 4.0 (Scanmarker and screw)	ZBAD5531	-	-	10-12
SET ZZ BASE OS-GS-CA-SM 5.0 (Scanmarker and screw)	ZBAD5532	-	-	10-12
SET ZZ BASE OS-GS-CA-SM 6.0 (Scanmarker and screw)	ZBAD5533	-	-	10-12
ZZ BASE OS-GS-CA-SC 4.0 (Screw)	BSAE5501	30	-	10-12
ZZ BASE OS-GS-CA-SC 5.0 (Screw)	BSAE5502	30	-	10-12
ZZ BASE OS-GS-CA-SC 6.0 (Screw)	BSAE5503	30	-	10-12
ZZ BASE OS-GS-CA-ZSC 4.0 (Zirconia screw)	BSAE5531	30	-	10-12
ZZ BASE OS-GS-CA-ZSC 5.0 (Zirconia screw)	BSAE5532	30	-	10-12
ZZ BASE OS-GS-CA-ZSC 6.0 (Zirconia screw)	BSAE5533	30	-	10-12

OSSTEM IMPLANT GS/TS

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE OS-GS-ABCH 3.5 (Raw-Abutment® customizable hex and screw)	BSAF4282	20	-	10-12
SET ZZ BASE OS-GS-ABCH 4.0/4.5/5.0 (Raw-Abutment® customizable hex and screw)	BSAF4281	30	-	10-12
SET ZZ BASE OS-GS-ABCH-D14 3.5 (Raw-Abutment® customizable hex D14 and screw)	BSAF4282D14	20	-	10-12
SET ZZ BASE OS-GS-ABCH-D14 4.0/4.5/5.0 (Raw-Abutment® customizable hex D14 and screw)	BSAF4281D14	30	-	10-12
SET ZZ BASE OS-GS-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD4232	20	-	10-12
SET ZZ BASE OS-GS-CCTB 4.0/4.5/5.0 (Conical cemented titan base non hex and screw)	BSAD4231	30	-	10-12
SET ZZ BASE OS-GS-CCTB-L10 3.5 (Conical cemented titan base non hex L10 and screw)	BSAD4232L10	20	-	10-12
SET ZZ BASE OS-GS-CCTB-L10 4.0/4.5/5.0 (Conical cemented titan base non hex L10 and screw)	BSAD4231L10	30	-	10-12
SET ZZ BASE OS-GS-NTB 3.5 (Narrow titan base non hex and screw)	BSAP4232	20	-	10-12
SET ZZ BASE OS-GS-NTB 4.0/4.5/5.0 (Narrow titan base non hex and screw)	BSAP4231	30	-	10-12
SET ZZ BASE OS-GS-NTBH 3.5 (Narrow titan base hex and screw)	BSAP4282	20	-	10-12
SET ZZ BASE OS-GS-NTBH 4.0/4.5/5.0 (Narrow titan base hex and screw)	BSAP4281	30	-	10-12
SET ZZ BASE OS-GS-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB4282	20	-	10-12
SET ZZ BASE OS-GS-PCTBH 4.0/4.5/5.0 (Parallel cemented titan base hex and screw)	BSAB4281	30	-	10-12
SET ZZ BASE OS-GS-PCTBH-L10 3.5 (Parallel cemented titan base hex L10 and screw)	BSAB4292L10	20	-	10-12
SET ZZ BASE OS-GS-PCTBH-L10 4.0/4.5/5.0 (Parallel cemented titan base hex L10 and screw)	BSAB4291L10	30	-	10-12
SET ZZ BASE OS-GS-SM 3.5 (Scanmarker and screw)	ZBAD4232	-	-	10-12
SET ZZ BASE OS-GS-SM 4.0/4.5/5.0 (Scanmarker and screw)	ZBAD4231	-	-	10-12
SET ZZ BASE OS-GS-WSM 3.5 (White Scanmarker and screw)	ZBCD4232	-	-	10-12
SET ZZ BASE OS-GS-WSM 4.0/4.5/5.0 (White Scanmarker and screw)	ZBCD4231	-	-	10-12
ZZ BASE OS-GS-SC 3.5 (Screw)	BSAE4202	20	-	10-12
ZZ BASE OS-GS-SC 4.0/4.5/5.0 (Screw)	BSAE4201	30	-	10-12
ZZ BASE OS-GS-ZSC 3.5 (Zirconia screw)	BSAE4232	20	-	10-12
ZZ BASE OS-GS-ZSC 4.0/4.5/5.0 (Zirconia screw)	BSAE4231	30	-	10-12

OSSTEM IMPLANT GS/TS MULTI ABUTMENT

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE OS-GS-CCTB MA (Conical cemented titan base non hex and screw)	BSAD9131	20	-	10-12
SET ZZ BASE OS-GS-SM MA (Scanmarker and screw)	ZBAD9131	-	-	10-12
SET ZZ BASE OS-GS-WSM MA (White Scanmarker and screw)	ZBCD9131	-	-	10-12
ZZ BASE OS-GS-MA (Model implant)	MIAG9101	-	-	10-12
ZZ BASE OS-GS-SC MA (Screw)	BSAE9101	20	-	10-12
ZZ BASE OS-GS-ZSC MA (Zirconia Screw)	BSAE9131	20	-	10-12

OSSTEM IMPLANT US

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE OS-US-ABCH-D14-ID 3.5 (Raw-Abutment® customizable hex D14 ID and screw)	BSAF4391D14	30	-	10-12
SET ZZ BASE OS-US-ABCH-D14-ID 4.1 (Raw-Abutment® customizable hex D14 ID and screw)	BSAF4392D14	30	-	10-12
SET ZZ BASE OS-US-ABCH-D14-ID 5.1 (Raw-Abutment® customizable hex D14 ID and screw)	BSAF4393D14	30	-	10-12
SET ZZ BASE OS-US-ABCH-ID 3.5 (Raw-Abutment® customizable hex ID and screw)	BSAF4391	30	-	10-12
SET ZZ BASE OS-US-ABCH-ID 4.1 (Raw-Abutment® customizable hex ID and screw)	BSAF4392	30	-	10-12
SET ZZ BASE OS-US-ABCH-ID 5.1 (Raw-Abutment® customizable hex ID and screw)	BSAF4393	30	-	10-12
SET ZZ BASE OS-US-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD4331	30	-	10-12
SET ZZ BASE OS-US-CCTB 4.1 (Conical cemented titan base non hex and screw)	BSAD4332	30	-	10-12
SET ZZ BASE OS-US-CCTB 5.1 (Conical cemented titan base non hex and screw)	BSAD4333	30	-	10-12
SET ZZ BASE OS-US-CCTB-L10 3.5 (Conical cemented titan base non hex L10 and screw)	BSAD4331L10	30	-	10-12
SET ZZ BASE OS-US-CCTB-L10 4.1 (Conical cemented titan base non hex L10 and screw)	BSAD4332L10	30	-	10-12
SET ZZ BASE OS-US-NTB 3.5 (Narrow titan base non hex and screw)	BSAP4331	30	-	10-12
SET ZZ BASE OS-US-NTBH 3.5 (Narrow titan base hex and screw)	BSAP4381	30	-	10-12
SET ZZ BASE OS-US-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB4381	30	-	10-12
SET ZZ BASE OS-US-PCTBH 4.1 (Parallel cemented titan base hex and screw)	BSAB4382	30	-	10-12
SET ZZ BASE OS-US-PCTBH 5.1 (Parallel cemented titan base hex and screw)	BSAB4383	30	-	10-12
SET ZZ BASE OS-US-PCTBH-L10 3.5 (Parallel cemented titan base hex L10 and screw)	BSAB4391L10	30	-	10-12
SET ZZ BASE OS-US-PCTBH-L10 4.1 (Parallel cemented titan base hex L10 and screw)	BSAB4392L10	30	-	10-12
SET ZZ BASE OS-US-SM 3.5 (Scanmarker and screw)	ZBAD4331	-	-	10-12
SET ZZ BASE OS-US-SM 4.1 (Scanmarker and screw)	ZBAD4332	-	-	10-12
SET ZZ BASE OS-US-SM 5.1 (Scanmarker and screw)	ZBAD4333	-	-	10-12
SET ZZ BASE OS-US-WSM 3.5 (White Scanmarker and screw)	ZBCD4331	-	-	10-12
SET ZZ BASE OS-US-WSM 4.1 (White Scanmarker and screw)	ZBCD4332	-	-	10-12
SET ZZ BASE OS-US-WSM 5.1 (White Scanmarker and screw)	ZBCD4333	-	-	10-12
ZZ BASE OS-US-SC 3.5 (Screw)	BSAE4301	30	-	10-12
ZZ BASE OS-US-SC 4.1 (Screw)	BSAE4302	30	-	10-12
ZZ BASE OS-US-SC 5.1 (Screw)	BSAE4303	30	-	10-12
ZZ BASE OS-US-ZSC 3.5 (Zirconia screw)	BSAE4331	30	-	10-12
ZZ BASE OS-US-ZSC 4.1 (Zirconia screw)	BSAE4332	30	-	10-12
ZZ BASE OS-US-ZSC 5.1 (Zirconia screw)	BSAE4333	30	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

PHIBO® TSA®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE P-TS-CCTB S3 (Conical cemented titan base non hex and screw)	BSAD8431	35	-	10-12
SET ZZ BASE P-TS-CCTB S4 (Conical cemented titan base non hex and screw)	BSAD8432	35	-	10-12
SET ZZ BASE P-TS-CCTB S5 (Conical cemented titan base non hex and screw)	BSAD8433	35	-	10-12
SET ZZ BASE P-TS-PCTBH S3 (Parallel cemented titan base hex and screw)	BSAB8481	35	-	10-12
SET ZZ BASE P-TS-PCTBH S4 (Parallel cemented titan base hex and screw)	BSAB8482	35	-	10-12
SET ZZ BASE P-TS-PCTBH S5 (Parallel cemented titan base hex and screw)	BSAB8483	35	-	10-12
SET ZZ BASE P-TS-SM S3 (Scanmarker and screw)	ZBAD8431	-	-	10-12
SET ZZ BASE P-TS-SM S4 (Scanmarker and screw)	ZBAD8432	-	-	10-12
SET ZZ BASE P-TS-SM S5 (Scanmarker and screw)	ZBAD8433	-	-	10-12
SET ZZ BASE P-TS-WSM S3 (White Scanmarker and screw)	ZBCD8431	-	-	10-12
SET ZZ BASE P-TS-WSM S4 (White Scanmarker and screw)	ZBCD8432	-	-	10-12
SET ZZ BASE P-TS-WSM S5 (White Scanmarker and screw)	ZBCD8433	-	-	10-12
ZZ BASE P-TS-SC S3 (Screw)	BSAE8401	35	-	10-12
ZZ BASE P-TS-SC S4 (Screw)	BSAE8402	35	-	10-12
ZZ BASE P-TS-SC S5 (Screw)	BSAE8403	35	-	10-12
ZZ BASE P-TS-ZSC S3 (Zirconia Screw)	BSAE8431	35	-	10-12
ZZ BASE P-TS-ZSC S4 (Zirconia Screw)	BSAE8432	35	-	10-12
ZZ BASE P-TS-ZSC S5 (Zirconia Screw)	BSAE8433	35	-	10-12

PRE-MILLED ABUTMENT

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
Straumann® NC Pre-milled Abut. Blank, Meden. Ø11.5, TAN	ZZ-022-0028	35	-	10-12
Straumann® NC Pre-milled Abut. Blank, Meden. Ø15.8, TAN	ZZ-022-0029	35	-	10-12
Straumann® RC Pre-milled Abut. Blank, Meden. Ø11.5, TAN	ZZ-022-0030	35	-	10-12
Straumann® RC Pre-milled Abut. Blank, Meden. Ø15.8, TAN	ZZ-022-0031	35	-	10-12
Straumann® RN Pre-milled Abut. Blank, Meden. Ø11.5, TAN	ZZ-022-0034	35	-	10-12
Straumann® RN Pre-milled Abut. Blank, Meden. Ø15.8, TAN	ZZ-022-0035	35	-	10-12
Straumann® VN Pre-milled Abut. Blank, Meden. Ø11.5, TAN	ZZ-022-0036	35	-	10-12
Straumann® VN Pre-milled Abut. Blank, Meden. Ø15.8, TAN	ZZ-022-0037	35	-	10-12

SIC® SICACE

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE S-S-ABCH-ID 3.3 (Raw-Abutment® customizable hex ID and screw)	BSAF1291	20	-	10-12
SET ZZ BASE S-S-ABCH-ID 4.2 (Raw-Abutment® customizable hex ID and screw)	BSAF1292	20	-	10-12
SET ZZ BASE S-S-CCTB 3.3 (Conical cemented titan base non hex and screw)	BSAD1231	20	-	10-12
SET ZZ BASE S-S-CCTB 4.2 (Conical cemented titan base non hex and screw)	BSAD1232	20	-	10-12
SET ZZ BASE S-S-NTB 3.3 (Narrow titan base non hex and screw)	BSAP1231	20	-	10-12
SET ZZ BASE S-S-NTBH 3.3 (Narrow titan base hex and screw)	BSAP1281	20	-	10-12
SET ZZ BASE S-S-PCTBH 3.3 (Parallel cemented titan base hex and screw)	BSAB1281	20	-	10-12
SET ZZ BASE S-S-PCTBH 4.2 (Parallel cemented titan base hex and screw)	BSAB1282	20	-	10-12
SET ZZ BASE S-S-SM 3.3 (Scanmarker and screw)	ZBAD1231	-	-	10-12
SET ZZ BASE S-S-SM 4.2 (Scanmarker and screw)	ZBAD1232	-	-	10-12
ZZ BASE S-S-SC 3.3 (Screw)	BSAE1201	20	-	10-12
ZZ BASE S-S-SC 4.2 (Screw)	BSAE1202	20	-	10-12
ZZ BASE S-S-SC-G 3.3 (Screw Gold)	BSCE1201	20	-	10-12
ZZ BASE S-S-SC-G 4.2 (Screw Gold)	BSCE1202	20	-	10-12
ZZ BASE S-S-ZSC 3.3 (Zirconia screw)	BSAE1231	20	-	10-12
ZZ BASE S-S-ZSC 4.2 (Zirconia screw)	BSAE1232	20	-	10-12
ZZ BASE S-S-ZSC-G 3.3 (Zirconia Screw Gold)	BSCE1231	20	-	10-12
ZZ BASE S-S-ZSC-G 4.2 (Zirconia Screw Gold)	BSCE1232	20	-	10-12

SOUTHERN IMPLANTS DEEP CONICAL

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE SI-DC-ABCH 3.0 (Raw-Abutment® customizable hex and screw)	BSAF9781	15	-	10-12
SET ZZ BASE SI-DC-ABCH 3.5/4.0 (Raw-Abutment® customizable hex and screw)	BSAF9782	20	-	10-12
SET ZZ BASE SI-DC-CCTB 3.0 (Conical cemented titan base non hex and screw)	BSAD9731	15	-	10-12
SET ZZ BASE SI-DC-CCTB 3.5/4.0 (Conical cemented titan base non hex and screw)	BSAD9732	20	-	10-12
SET ZZ BASE SI-DC-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD9733	-	-	10-12
SET ZZ BASE SI-DC-PCTBH 3.0 (Parallel cemented titan base hex and screw)	BSAB9781	15	-	10-12
SET ZZ BASE SI-DC-PCTBH 3.5/4.0 (Parallel cemented titan base hex and screw)	BSAB9782	20	-	10-12
SET ZZ BASE SI-DC-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB9783	-	-	10-12
SET ZZ BASE SI-DC-SM 3.0 (Scanmarker and screw)	ZBAD9731	-	-	10-12
SET ZZ BASE SI-DC-SM 3.5/4.0 (Scanmarker and screw)	ZBAD9732	-	-	10-12
SET ZZ BASE SI-DC-SM 5.0 (Scanmarker and screw)	ZBAD9733	-	-	10-12
SET ZZ BASE SI-DC-WSM 3.0 (White Scanmarker and screw)	ZBCD9731	-	-	10-12
SET ZZ BASE SI-DC-WSM 3.5/4.0 (White Scanmarker and screw)	ZBCD9732	-	-	10-12
SET ZZ BASE SI-DC-WSM 5.0 (White Scanmarker and screw)	ZBCD9733	-	-	10-12
ZZ BASE SI-DC-SC 3.0 (Screw)	BSAE9701	15	-	10-12
ZZ BASE SI-DC-SC 3.5/4.0 (Screw)	BSAE9702	20	-	10-12
ZZ BASE SI-DC-SC 5.0 (Screw)	BSAE9703	-	-	10-12
ZZ BASE SI-DC-ZSC 3.0 (Zirconia Screw)	BSAE9731	15	-	10-12
ZZ BASE SI-DC-ZSC 3.5/4.0 (Zirconia Screw)	BSAE9732	20	-	10-12
ZZ BASE SI-DC-ZSC 5.0 (Zirconia Screw)	BSAE9733	-	-	10-12

SOUTHERN IMPLANTS INTERNAL HEX

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE SI-IH-ABCH 3.7/4.2/5.0 (Raw-Abutment® customizable hex and screw)	BSAF9881	32	-	10-12
SET ZZ BASE SI-IH-CCTB 3.7/4.2/5.0 (Conical cemented titan base non hex and screw)	BSAD9831	32	-	10-12
SET ZZ BASE SI-IH-PCTBH 3.7/4.2/5.0 (Parallel cemented titan base hex and screw)	BSAB9881	32	-	10-12
SET ZZ BASE SI-IH-SM 3.7/4.2/5.0 (Scanmarker and screw)	ZBAD9831	-	-	10-12
SET ZZ BASE SI-IH-WSM 3.7/4.2/5.0 (White Scanmarker and screw)	ZBCD9831	-	-	10-12
ZZ BASE SI-IH-SC 3.7/4.2/5.0 (Screw)	BSAE9801	32	-	10-12
ZZ BASE SI-IH-ZSC 3.7/4.2/5.0 (Zirconia Screw)	BSAE9831	32	-	10-12

SOUTHERN IMPLANTS IT CONNECTION

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE SI-IC-ABCH 4.8 (Raw-Abutment® customizable hex and screw)	BSAF9681	25	-	10-12
SET ZZ BASE SI-IC-ABCH 6.5 (Raw-Abutment® customizable hex and screw)	BSAF9682	25	-	10-12
SET ZZ BASE SI-IC-CCTB 4.8 (Conical cemented titan base non hex and screw)	BSAD9631	25	-	10-12
SET ZZ BASE SI-IC-CCTB 6.5 (Conical cemented titan base non hex and screw)	BSAD9632	25	-	10-12
SET ZZ BASE SI-IC-PCTBH 4.8 (Parallel cemented titan base hex and screw)	BSAB9681	25	-	10-12
SET ZZ BASE SI-IC-PCTBH 6.5 (Parallel cemented titan base hex and screw)	BSAB9682	25	-	10-12
SET ZZ BASE SI-IC-SM 4.8 (Scanmarker and screw)	ZBAD9631	-	-	10-12
SET ZZ BASE SI-IC-SM 6.5 (Scanmarker and screw)	ZBAD9632	-	-	10-12
SET ZZ BASE SI-IC-WSM 4.8 (White Scanmarker and screw)	ZBCD9631	-	-	10-12
SET ZZ BASE SI-IC-WSM 6.5 (White Scanmarker and screw)	ZBCD9632	-	-	10-12
ZZ BASE SI-IC-SC 4.8 (Screw)	BSAE9601	25	-	10-12
ZZ BASE SI-IC-SC 6.5 (Screw)	BSAE9602	25	-	10-12
ZZ BASE SI-IC-ZSC 4.8 (Zirconia Screw)	BSAE9631	25	-	10-12
ZZ BASE SI-IC-ZSC 6.5 (Zirconia Screw)	BSAE9632	25	-	10-12

STRAUMANN® BONE LEVEL

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE ST-BL-ABCH NC 3.3 (Raw-Abutment® customizable hex and screw)	BSAF0581	35	-	10-12
SET ZZ BASE ST-BL-ABCH RC 4.1 (Raw-Abutment® customizable hex and screw)	BSAF0582	35	-	10-12
SET ZZ BASE ST-BL-ABCH SC (Raw-Abutment® customizable hex and screw)	BSAF0583	35	-	10-12
SET ZZ BASE ST-BL-ABCH-D14 NC 3.3 (Raw-Abutment® customizable hex D14 and screw)	BSAF0581D14	35	-	10-12
SET ZZ BASE ST-BL-ABCH-D14 RC 4.1 (Raw-Abutment® customizable hex D14 and screw)	BSAF0582D14	35	-	10-12
SET ZZ BASE ST-BL-CCTB NC 3.3 (Conical cemented titan base non hex and screw)	BSAD0531	35	-	10-12
SET ZZ BASE ST-BL-CCTB RC 4.1 (Conical cemented titan base non hex and screw)	BSAD0532	35	-	10-12
SET ZZ BASE ST-BL-CCTB SC (Conical cemented titan base non hex and screw)	BSAD0533	35	-	10-12
SET ZZ BASE ST-BL-CCTB-L10 NC 3.3 (Conical cemented titan base non hex L10 and screw)	BSAD0531L10	35	-	10-12
SET ZZ BASE ST-BL-CCTB-L10 RC 4.1 (Conical cemented titan base non hex L10 and screw)	BSAD0532L10	35	-	10-12
SET ZZ BASE ST-BL-CCTB-L20 NC 3.3 (Conical cemented titan base non hex L20 and screw)	BSAD0531L20	35	-	10-12
SET ZZ BASE ST-BL-CCTB-L20 RC 4.1 (Conical cemented titan base non hex L20 and screw)	BSAD0532L20	35	-	10-12
SET ZZ BASE ST-BL-CCTB-L30 NC 3.3 (Conical cemented titan base non hex L30 and screw)	BSAD0531L30	35	-	10-12
SET ZZ BASE ST-BL-CCTB-L30 RC 4.1 (Conical cemented titan base non hex L30 and screw)	BSAD0532L30	35	-	10-12
SET ZZ BASE ST-BL-NTB NC 3.3 (Narrow titan base non hex and screw)	BSAP0531	35	-	10-12
SET ZZ BASE ST-BL-NTB RC 4.1 (Narrow titan base non hex and screw)	BSAP0532	35	-	10-12
SET ZZ BASE ST-BL-NTB SC (Narrow titan base non hex and screw)	BSAP0533	35	-	10-12
SET ZZ BASE ST-BL-NTBH NC 3.3 (Narrow titan base hex and screw)	BSAP0581	35	-	10-12
SET ZZ BASE ST-BL-NTBH RC 4.1 (Narrow titan base hex and screw)	BSAP0582	35	-	10-12
SET ZZ BASE ST-BL-NTBH SC (Narrow titan base hex and screw)	BSAP0583	35	-	10-12
SET ZZ BASE ST-BL-PCTBH NC 3.3 (Parallel cemented titan base hex and screw)	BSAB0591	35	-	10-12
SET ZZ BASE ST-BL-PCTBH RC 4.1 (Parallel cemented titan base hex and screw)	BSAB0592	35	-	10-12
SET ZZ BASE ST-BL-PCTBH SC (Parallel cemented titan base hex and screw)	BSAB0583	35	-	10-12
SET ZZ BASE ST-BL-PCTBH-L20 NC 3.3 (Parallel cemented titan base hex L20 and screw)	BSAB0591L20	35	-	10-12
SET ZZ BASE ST-BL-PCTBH-L20 RC 4.1 (Parallel cemented titan base hex L20 and screw)	BSAB0592L20	35	-	10-12
SET ZZ BASE ST-BL-PCTBH-L30 NC 3.3 (Parallel cemented titan base hex L30 and screw)	BSAB0591L30	35	-	10-12
SET ZZ BASE ST-BL-PCTBH-L30 RC 4.1 (Parallel cemented titan base hex L30 and screw)	BSAB0592L30	35	-	10-12
SET ZZ BASE ST-BL-PCTBH-L40 NC 3.3 (Parallel cemented titan base hex L40 and screw)	BSAB0591L40	35	-	10-12
SET ZZ BASE ST-BL-PCTBH-L40 RC 4.1 (Parallel cemented titan base hex L40 and screw)	BSAB0592L40	35	-	10-12
SET ZZ BASE ST-BL-SA NC 3.3 (Scananalog and screw)	ZBDE0531	15	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

STRAUMANN® BONE LEVEL

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE ST-BL-SA RC 4.1 (Scanalog and screw)	ZBDE0532	15	-	10-12
SET ZZ BASE ST-BL-SM NC 3.3 (Scanmarker and screw)	ZBAD0531	-	-	10-12
SET ZZ BASE ST-BL-SM RC 4.1 (Scanmarker and screw)	ZBAD0532	-	-	10-12
SET ZZ BASE ST-BL-SM SC (Scanmarker and screw)	ZBAD0533	-	-	10-12
SET ZZ BASE ST-BL-WSM NC 3.3 (White Scanmarker and screw)	ZBCD0531	-	-	10-12
SET ZZ BASE ST-BL-WSM RC 4.1 (White Scanmarker and screw)	ZBCD0532	-	-	10-12
SET ZZ BASE ST-BL-WSM SC (White Scanmarker and screw)	ZBCD0533	-	-	10-12
SET ZZ BASE ST-BL-ZZMUA-L15 NC 3.3 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0531L15	25	30	10-12
SET ZZ BASE ST-BL-ZZMUA-L15 RC 4.1 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0532L15	25	30	10-12
SET ZZ BASE ST-BL-ZZMUA-L25 NC 3.3 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0531L25	25	30	10-12
SET ZZ BASE ST-BL-ZZMUA-L25 RC 4.1 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0532L25	25	30	10-12
SET ZZ BASE ST-BL-ZZMUA-L35 NC 3.3 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0531L35	25	30	10-12
SET ZZ BASE ST-BL-ZZMUA-L35 RC 4.1 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0532L35	25	30	10-12
SET ZZ BASE ST-BL-ZZMUA-L45 NC 3.3 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0531L45	25	30	10-12
SET ZZ BASE ST-BL-ZZMUA-L45 RC 4.1 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0532L45	25	30	10-12
ZZ BASE ST-BL-MI NC 3.3 (Model implant)	MIAG0501	-	-	10-12
ZZ BASE ST-BL-MI RC 4.1 (Model implant)	MIAG0502	-	-	10-12
ZZ BASE ST-BL-MI SC (Model implant)	MIAG0503	-	-	10-12
ZZ BASE ST-BL-SC NC 3.3 (Screw)	BSAE0501	35	-	10-12
ZZ BASE ST-BL-SC RC 4.1 (Screw)	BSAE0502	35	-	10-12
ZZ BASE ST-BL-SC SC (Screw)	BSAE0503	35	-	10-12
ZZ BASE ST-BL-SC-G NC 3.3 (Screw Gold)	BSCE0501	35	-	10-12
ZZ BASE ST-BL-SC-G RC 4.1 (Screw Gold)	BSCE0502	35	-	10-12
ZZ BASE ST-BL-ZSC NC 3.3 (Zirconia screw)	BSAE0531	35	-	10-12
ZZ BASE ST-BL-ZSC RC 4.1 (Zirconia screw)	BSAE0532	35	-	10-12
ZZ BASE ST-BL-ZSC SC (Zirconia screw)	BSAE0533	35	-	10-12
ZZ BASE ST-BL-ZSC-G RC 4.1 (Zirconia Screw Gold)	BSCE0532	35	-	10-12

STRAUMANN® MULTI-BASE ABUTMENT

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE ST-MA-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD5331	15	-	10-12
SET ZZ BASE ST-MA-CCTB 4.5 (Conical cemented titan base non hex and screw)	BSAD5332	15	-	10-12
SET ZZ BASE ST-MA-SM 3.5 (Scanmarker and screw)	ZBAD5331	-	-	10-12
SET ZZ BASE ST-MA-SM 4.5 (Scanmarker and screw)	ZBAD5332	-	-	10-12
ZZ BASE ST-MA-SC 3.5 (Screw)	BSAE5301	15	-	10-12
ZZ BASE ST-MA-SC 4.5 (Screw)	BSAE5302	15	-	10-12
ZZ BASE ST-MA-ZSC 3.5 (Zirconia screw)	BSAE5331	15	-	10-12
ZZ BASE ST-MA-ZSC 4.5 (Zirconia screw)	BSAE5332	15	-	10-12

STRAUMANN® SCREW-RETAINED ABUTMENT

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE ST-SA-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD8031	15	-	10-12
SET ZZ BASE ST-SA-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD8041	15	-	10-12
SET ZZ BASE ST-SA-CCTB 4.6 (Conical cemented titan base non hex and screw)	BSAD8032	15	-	10-12
SET ZZ BASE ST-SA-CCTB 4.6 (Conical cemented titan base non hex and screw)	BSAD8042	15	-	10-12
SET ZZ BASE ST-SA-SM 3.5 (Scanmarker and screw)	ZBAD8031	-	-	10-12
SET ZZ BASE ST-SA-SM 4.6 (Scanmarker and screw)	ZBAD8032	-	-	10-12
SET ZZ BASE ST-SA-WSM 3.5 (White Scanmarker and screw)	ZBCD8031	-	-	10-12
SET ZZ BASE ST-SA-WSM 4.6 (White Scanmarker and screw)	ZBCD8032	-	-	10-12
ZZ BASE ST-SA-SC 3.5 (Screw)	BSAE8001	15	-	10-12
ZZ BASE ST-SA-SC 4.6 (Screw)	BSAE8002	15	-	10-12
ZZ BASE ST-SA-ZSC 3.5 (Zirconia Screw)	BSAE8031	15	-	10-12
ZZ BASE ST-SA-ZSC 4.6 (Zirconia Screw)	BSAE8032	15	-	10-12

STRAUMANN® TISSUE LEVEL SYNOCTA®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE ST-SO-ABCH NNC 3.5 (Raw-Abutment® customizable hex and screw)	BSAF0684	35	-	10-12
SET ZZ BASE ST-SO-ABCH RN 4.8 (Raw-Abutment® customizable hex and screw)	BSAF0682	35	-	10-12
SET ZZ BASE ST-SO-ABCH WN 6.5 (Raw-Abutment® customizable hex and screw)	BSAF0683	35	-	10-12
SET ZZ BASE ST-SO-ABCH-D14 NNC 3.5 (Raw-Abutment® customizable hex D14 and screw)	BSAF0684D14	35	-	10-12
SET ZZ BASE ST-SO-ABCH-D14 RN 4.8 (Raw-Abutment® customizable hex D14 and screw)	BSAF0682D14	35	-	10-12
SET ZZ BASE ST-SO-ABCH-D14 WN 6.5 (Raw-Abutment® customizable hex D14 and screw)	BSAF0683D14	35	-	10-12
SET ZZ BASE ST-SO-CCTB NN 3.5 (Conical cemented titan base non hex and screw)	BSAD0631	35	-	10-12
SET ZZ BASE ST-SO-CCTB NNC 3.5 (Conical cemented titan base non hex and screw)	BSAD0634	35	-	10-12
SET ZZ BASE ST-SO-CCTB RN 4.8 (Conical cemented titan base non hex and screw)	BSAD0642	35	-	10-12
SET ZZ BASE ST-SO-CCTB RN 4.8 (Conical cemented titan base non hex and screw)	BSAD0632	35	-	10-12
SET ZZ BASE ST-SO-CCTB WN 6.5 (Conical cemented titan base non hex and screw)	BSAD0643	35	-	10-12
SET ZZ BASE ST-SO-CCTB WN 6.5 (Conical cemented titan base non hex and screw)	BSAD0633	35	-	10-12
SET ZZ BASE ST-SO-CCTB-G NN 3.5 (Conical cemented titan base non hex and screw Gold)	BSCD0631	35	-	10-12
SET ZZ BASE ST-SO-CCTB-L15 RN 4.8 (Conical cemented titan base non hex L15 and screw)	BSAD0632L15	35	-	10-12
SET ZZ BASE ST-SO-CCTB-L15 WN 6.5 (Conical cemented titan base non hex L15 and screw)	BSAD0633L15	35	-	10-12
SET ZZ BASE ST-SO-CCTB-L20 NNC 3.5 (Conical cemented titan base non hex L20 and screw)	BSAD0634L20	35	-	10-12
SET ZZ BASE ST-SO-CCTB-L30 NNC 3.5 (Conical cemented titan base non hex L30 and screw)	BSAD0634L30	35	-	10-12
SET ZZ BASE ST-SO-NTB NNC 3.5 (Narrow titan base non hex and screw)	BSAP0634	35	-	10-12
SET ZZ BASE ST-SO-NTBH NNC 3.5 (Narrow titan base hex and screw)	BSAP0684	35	-	10-12
SET ZZ BASE ST-SO-PCTBH NN 3.5 (Parallel cemented titan base hex and screw)	BSAB0681	35	-	10-12
SET ZZ BASE ST-SO-PCTBH NNC 3.5 (Parallel cemented titan base hex and screw)	BSAB0684	35	-	10-12
SET ZZ BASE ST-SO-PCTBH RN 4.8 (Parallel cemented titan base hex and screw)	BSAB0682	35	-	10-12
SET ZZ BASE ST-SO-PCTBH RN 4.8 (Parallel cemented titan base hex and screw)	BSAB0672	35	-	10-12
SET ZZ BASE ST-SO-PCTBH RN 4.8 (Parallel cemented titan base hex and screw)	BSAB0692	35	-	10-12
SET ZZ BASE ST-SO-PCTBH WN 6.5 (Parallel cemented titan base hex and screw)	BSAB0683	35	-	10-12
SET ZZ BASE ST-SO-PCTBH WN 6.5 (Parallel cemented titan base hex and screw)	BSAB0693	35	-	10-12
SET ZZ BASE ST-SO-PCTBH-L10 NNC 3.5 (Parallel cemented titan base hex L10 and screw)	BSAB0694L10	35	-	10-12
SET ZZ BASE ST-SO-PCTBH-L15 RN 4.8 (Parallel cemented titan base hex L15 and screw)	BSAB0692L15	35	-	10-12
SET ZZ BASE ST-SO-PCTBH-L15 WN 6.5 (Parallel cemented titan base hex L15 and screw)	BSAB0693L15	35	-	10-12
SET ZZ BASE ST-SO-PCTBH-L20 NNC 3.5 (Parallel cemented titan base hex L20 and screw)	BSAB0694L20	35	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

STRAUMANN® TISSUE LEVEL SYNOCTA®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE ST-SO-PCTBH-L30 NNC 3.5 (Parallel cemented titan base hex L30 and screw)	BSAB0694L30	35	-	10-12
SET ZZ BASE ST-SO-SA NNC 3.5 (Scanalog and screw)	ZBDE0634	15	-	10-12
SET ZZ BASE ST-SO-SA RN 4.8 (Scanalog and screw)	ZBDE0632	15	-	10-12
SET ZZ BASE ST-SO-SA WN 6.5 (Scanalog and screw)	ZBDE0633	15	-	10-12
SET ZZ BASE ST-SO-SM NN 3.5 (Scanmarker and screw)	ZBAD0631	-	-	10-12
SET ZZ BASE ST-SO-SM NNC 3.5 (Scanmarker and screw)	ZBAD0634	-	-	10-12
SET ZZ BASE ST-SO-SM RN 4.8 (Scanmarker and screw)	ZBAD0632	-	-	10-12
SET ZZ BASE ST-SO-SM RN 4.8 (Scanmarker with Synocta and screw)	ZBAD0642	-	-	10-12
SET ZZ BASE ST-SO-SM WN 6.5 (Scanmarker and screw)	ZBAD0633	-	-	10-12
SET ZZ BASE ST-SO-SM WN 6.5 (Scanmarker with Synocta and screw)	ZBAD0643	-	-	10-12
SET ZZ BASE ST-SO-WSM NNC 3.5 (White Scanmarker and screw)	ZBCD0634	-	-	10-12
SET ZZ BASE ST-SO-WSM RN 4.8 (White Scanmarker and screw)	ZBCD0632	-	-	10-12
SET ZZ BASE ST-SO-WSM WN 6.5 (White Scanmarker and screw)	ZBCD0633	-	-	10-12
SET ZZ BASE ST-SO-ZZMUA-L05 WN 6.5 (Zirkonzahn Multi Unit Abutment non hex L05 and Screw)	BSAH0633L05	25	35	10-12
SET ZZ BASE ST-SO-ZZMUA-L15 NNC 3.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0634L15	25	35	10-12
SET ZZ BASE ST-SO-ZZMUA-L15 RN 4.8 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0632L15	25	35	10-12
SET ZZ BASE ST-SO-ZZMUA-L15 WN 6.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH0633L15	25	35	10-12
SET ZZ BASE ST-SO-ZZMUA-L25 NNC 3.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0634L25	25	35	10-12
SET ZZ BASE ST-SO-ZZMUA-L25 RN 4.8 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH0632L25	25	35	10-12
SET ZZ BASE ST-SO-ZZMUA-L35 NNC 3.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0634L35	25	35	10-12
SET ZZ BASE ST-SO-ZZMUA-L35 RN 4.8 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH0632L35	25	35	10-12
SET ZZ BASE ST-SO-ZZMUA-L45 NNC 3.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0634L45	25	35	10-12
SET ZZ BASE ST-SO-ZZMUA-L45 RN 4.8 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH0632L45	25	35	10-12
ZZ BASE ST-SO-BSC NN 3.5 (Black Screw)	BSIE0601	35	-	10-12
ZZ BASE ST-SO-BZSC NN 3.5 (Black Zirconia screw)	BSIE0631	35	-	10-12
ZZ BASE ST-SO-LBSC NN 3.5 (Lab Screw)	BSHE0601	35	-	10-12
ZZ BASE ST-SO-LBZSC NN 3.5 (Lab Zirconia screw)	BSHE0631	35	-	10-12
ZZ BASE ST-SO-LZSC RN 4.8 (Long zirconia screw)	BSAE0682	35	-	10-12
ZZ BASE ST-SO-LZSC WN 6.5 (Long zirconia screw)	BSAE0683	35	-	10-12
ZZ BASE ST-SO-MI NNC 3.5 (Model implant)	MIAG0604	-	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

STRAUMANN® TISSUE LEVEL SYNOCTA®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
ZZ BASE ST-SO-MI RN 4.8 (Model implant)	MIAG0602	-	-	10-12
ZZ BASE ST-SO-MI WN 6.5 (Model implant)	MIAG0603	-	-	10-12
ZZ BASE ST-SO-SC NN 3.5 (Screw)	BSAE0601	35	-	10-12
ZZ BASE ST-SO-SC NNC 3.5 (Screw)	BSAE0604	35	-	10-12
ZZ BASE ST-SO-SC RN 4.8 (Screw)	BSAE0602	35	-	10-12
ZZ BASE ST-SO-SC WN 6.5 (Screw)	BSAE0603	35	-	10-12
ZZ BASE ST-SO-SC-G 3.5 (Screw Gold)	BSCE0601	35	-	10-12
ZZ BASE ST-SO-SC-G 4.8 (Screw Gold)	BSCE0602	35	-	10-12
ZZ BASE ST-SO-SC-G 6.5 (Screw Gold)	BSCE0603	35	-	10-12
ZZ BASE ST-SO-ZSC NN 3.5 (Zirconia screw)	BSAE0631	35	-	10-12
ZZ BASE ST-SO-ZSC NNC 3.5 (Zirconia screw)	BSAE0634	35	-	10-12
ZZ BASE ST-SO-ZSC RN 4.8 (Zirconia screw)	BSAE0632	35	-	10-12
ZZ BASE ST-SO-ZSC WN 6.5 (Zirconia screw)	BSAE0633	35	-	10-12
ZZ BASE ST-SO-ZSC-G NN 3.5 (Zirconia Screw Gold)	BSCE0631	35	-	10-12

SWEDEN & MARTINA KOHNO/PREMIUM

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE S&M-K-SM 4.25 (Scanmarker and screw)	ZBAD1533	-	-	10-12
SET ZZ BASE S&M-K-ABCH-ID 3.30 (Raw-Abutment® customizable hex ID and screw)	BSAF1591	30	-	10-12
SET ZZ BASE S&M-K-ABCH-ID 3.80 (Raw-Abutment® customizable hex ID and screw)	BSAF1592	30	-	10-12
SET ZZ BASE S&M-K-ABCH-ID 4.25 (Raw-Abutment® customizable hex ID and screw)	BSAF1593	35	-	10-12
SET ZZ BASE S&M-K-ABCH-ID 5.0 (Raw-Abutment® customizable hex ID and screw)	BSAF1594	35	-	10-12
SET ZZ BASE S&M-K-ABCH-ID 6.0 (Raw-Abutment® customizable hex ID and screw)	BSAF1595	35	-	10-12
SET ZZ BASE S&M-K-CCTB 3.30 (Conical cemented titan base non hex and screw)	BSAD1531	30	-	10-12
SET ZZ BASE S&M-K-CCTB 3.80 (Conical cemented titan base non hex and screw)	BSAD1532	30	-	10-12
SET ZZ BASE S&M-K-CCTB 4.25 (Conical cemented titan base non hex and screw)	BSAD1533	35	-	10-12
SET ZZ BASE S&M-K-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD1534	35	-	10-12
SET ZZ BASE S&M-K-CCTB 6.0 (Conical cemented titan base non hex and screw)	BSAD1535	35	-	10-12
SET ZZ BASE S&M-K-NTB 3.3 (Narrow titan base non hex and screw)	BSAP1531	30	-	10-12
SET ZZ BASE S&M-K-NTBH 3.3 (Narrow titan base hex and screw)	BSAP1581	30	-	10-12
SET ZZ BASE S&M-K-PCTBH 3.30 (Parallel cemented titan base hex and screw)	BSAB1581	30	-	10-12
SET ZZ BASE S&M-K-PCTBH 3.80 (Parallel cemented titan base hex and screw)	BSAB1582	30	-	10-12
SET ZZ BASE S&M-K-PCTBH 4.25 (Parallel cemented titan base hex and screw)	BSAB1583	35	-	10-12
SET ZZ BASE S&M-K-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB1584	35	-	10-12
SET ZZ BASE S&M-K-PCTBH 6.0 (Parallel cemented titan base hex and screw)	BSAB1585	35	-	10-12
SET ZZ BASE S&M-K-SM 3.3 (Scanmarker and screw)	ZBAD1531	-	-	10-12
SET ZZ BASE S&M-K-SM 3.8 (Scanmarker and screw)	ZBAD1532	-	-	10-12
SET ZZ BASE S&M-K-SM 5.0 (Scanmarker and screw)	ZBAD1534	-	-	10-12
SET ZZ BASE S&M-K-SM 6.0 (Scanmarker and screw)	ZBAD1535	-	-	10-12
SET ZZ BASE S&M-K-WSM 3.30 (White Scanmarker and screw)	ZBCD1531	-	-	10-12
SET ZZ BASE S&M-K-WSM 3.80 (White Scanmarker and screw)	ZBCD1532	-	-	10-12
SET ZZ BASE S&M-K-WSM 4.25 (White Scanmarker and screw)	ZBCD1533	-	-	10-12
SET ZZ BASE S&M-K-WSM 5.0 (White Scanmarker and screw)	ZBCD1534	-	-	10-12
SET ZZ BASE S&M-K-WSM 6.0 (White Scanmarker and screw)	ZBCD1535	-	-	10-12
ZZ BASE S&M-K-SC 3.30 (Screw)	BSAE1501	30	-	10-12
ZZ BASE S&M-K-SC 3.80 (Screw)	BSAE1502	30	-	10-12
ZZ BASE S&M-K-SC 4.25 (Screw)	BSAE1503	35	-	10-12

SWEDEN & MARTINA KOHNO/PREMIUM

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
ZZ BASE S&M-K-ZSC 3.30 (Zirconia screw)	BSAE1531	30	-	10-12
ZZ BASE S&M-K-ZSC 3.80 (Zirconia screw)	BSAE1532	30	-	10-12
ZZ BASE S&M-K-ZSC 4.25 (Zirconia screw)	BSAE1533	35	-	10-12
ZZ BASE S&M-K-SC 5.0 (Screw)	BSAE1504	35	-	10-12
ZZ BASE S&M-K-SC 6.0 (Screw)	BSAE1505	35	-	10-12
ZZ BASE S&M-K-SC-G 3.30 (Screw Gold)	BSCE1501	30	-	10-12
ZZ BASE S&M-K-SC-G 3.80 (Screw Gold)	BSCE1502	30	-	10-12
ZZ BASE S&M-K-SC-G 4.25 (Screw Gold)	BSCE1503	35	-	10-12
ZZ BASE S&M-K-SC-G 5.0 (Screw Gold)	BSCE1504	35	-	10-12
ZZ BASE S&M-K-SC-G 6.0 (Screw Gold)	BSCE1505	35	-	10-12
ZZ BASE S&M-K-ZSC 5.0 (Zirconia screw)	BSAE1534	35	-	10-12
ZZ BASE S&M-K-ZSC 6.0 (Zirconia screw)	BSAE1535	35	-	10-12
ZZ BASE S&M-K-ZSC-G 3.3 (Zirconia Screw Gold)	BSCE1531	30	-	10-12
ZZ BASE S&M-K-ZSC-G 3.8 (Zirconia Screw Gold)	BSCE1532	30	-	10-12
ZZ BASE S&M-K-ZSC-G 4.25 (Zirconia Screw Gold)	BSCE1533	35	-	10-12
ZZ BASE S&M-K-ZSC-G 5.00 (Zirconia Screw Gold)	BSCE1534	35	-	10-12
ZZ BASE S&M-K-ZSC-G 6.00 (Zirconia Screw Gold)	BSCE1535	35	-	10-12

SWEDEN & MARTINA OUT-LINK

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE S&M-O-SM 3.30 (Scanmarker and screw)	ZBAD1331	-	-	10-12
SET ZZ BASE S&M-O-SM 4.10 (Scanmarker and screw)	ZBAD1332	-	-	10-12
SET ZZ BASE S&M-O-CCTB 3.30 (Conical cemented titan base non hex and screw)	BSAD1331	30	-	10-12
SET ZZ BASE S&M-O-CCTB 4.10 (Conical cemented titan base non hex and screw)	BSAD1332	35	-	10-12
SET ZZ BASE S&M-O-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD1333	35	-	10-12
SET ZZ BASE S&M-O-NTB 3.3 (Narrow titan base non hex and screw)	BSAP1331	30	-	10-12
SET ZZ BASE S&M-O-NTB 4.1 (Narrow titan base non hex and screw)	BSAP1332	35	-	10-12
SET ZZ BASE S&M-O-NTBH 3.3 (Narrow titan base hex and screw)	BSAP1381	30	-	10-12
SET ZZ BASE S&M-O-NTBH 4.1 (Narrow titan base hex and screw)	BSAP1382	35	-	10-12
SET ZZ BASE S&M-O-PCTBH 3.3 (Parallel cemented titan base hex and screw)	BSAB1381	30	-	10-12
SET ZZ BASE S&M-O-PCTBH 4.1 (Parallel cemented titan base hex and screw)	BSAB1382	35	-	10-12
SET ZZ BASE S&M-O-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB1383	35	-	10-12
SET ZZ BASE S&M-O-SM 5.0 (Scanmarker and screw)	ZBAD1333	-	-	10-12
ZZ BASE S&M-O-SC 3.30 (Screw)	BSAE1301	30	-	10-12
ZZ BASE S&M-O-SC 4.10 (Screw)	BSAE1302	35	-	10-12
ZZ BASE S&M-O-ZSC 3.30 (Zirconia screw)	BSAE1331	30	-	10-12
ZZ BASE S&M-O-ZSC 4.10 (Zirconia screw)	BSAE1332	35	-	10-12
ZZ BASE S&M-O-SC 5.0 (Screw)	BSAE1303	35	-	10-12
ZZ BASE S&M-O-SC-G 3.3 (Screw Gold)	BSCE1301	30	-	10-12
ZZ BASE S&M-O-SC-G 4.1 (Screw Gold)	BSCE1302	35	-	10-12
ZZ BASE S&M-O-SC-G 5.0 (Screw Gold)	BSCE1303	35	-	10-12
ZZ BASE S&M-O-ZSC 5.0 (Zirconia screw)	BSAE1333	35	-	10-12
ZZ BASE S&M-O-ZSC-G 3.3 (Zirconia Screw Gold)	BSCE1331	30	-	10-12

SWEDEN & MARTINA P.A.D® MULTI UNIT ABUTMENT

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE S&M-MUA-SM (Scanmarker and screw)	ZBAD1431	-	-	10-12
SET ZZ BASE S&M-MUA-CCTB (Conical cemented titan base non hex and screw)	BSAD1431	15	-	10-12
ZZ BASE S&M-MUA-SC (Screw)	BSAE1401	15	-	10-12
ZZ BASE S&M-MUA-SC-G (Screw Gold)	BSCE1401	15	-	10-12
ZZ BASE S&M-MUA-ZSC (Zirconia screw)	BSAE1431	15	-	10-12

TEKKA® IN-KONE®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE T-IK-ABCH 3.5/4.0/4.5/5.5 (Raw-Abutment® customizable hex and screw)	BSAF9081	15	-	10-12
SET ZZ BASE T-IK-CCTB 0° (Conical cemented titan base non hex and screw)	BSAD9032	10	-	10-12
SET ZZ BASE T-IK-CCTB 17°/30° (Conical cemented titan base non hex and screw)	BSAD9033	10	-	10-12
SET ZZ BASE T-IK-CCTB 3.5/4.0/4.5/5.5 (Conical cemented titan base non hex and screw)	BSAD9031	15	-	10-12
SET ZZ BASE T-IK-PCTBH 3.5/4.0/4.5/5.5 (Parallel cemented titan base hex and screw)	BSAB9081	15	-	10-12
SET ZZ BASE T-IK-SM 0° (Scanmarker and screw)	ZBAD9032	-	-	10-12
SET ZZ BASE T-IK-SM 17°/30° (Scanmarker and screw)	ZBAD9033	-	-	10-12
SET ZZ BASE T-IK-SM 3.5/4.0/4.5/5.5 (Scanmarker and screw)	ZBAD9031	-	-	10-12
SET ZZ BASE T-IK-WSM 0° (White Scanmarker and screw)	ZBCD9032	-	-	10-12
SET ZZ BASE T-IK-WSM 17°/30° (White Scanmarker and screw)	ZBCD9033	-	-	10-12
SET ZZ BASE T-IK-WSM 3.5/4.0/4.5/5.5 (White Scanmarker and screw)	ZBCD9031	-	-	10-12
ZZ BASE T-IK-MI 0° (Model implant)	MIAG9002	-	-	10-12
ZZ BASE T-IK-MI 17°/30° (Model implant)	MIAG9003	-	-	10-12
ZZ BASE T-IK-MI 3.5/4.0/4.5/5.5 (Model implant)	MIAG9001	-	-	10-12
ZZ BASE T-IK-SC 0° (Screw)	BSAE9002	10	-	10-12
ZZ BASE T-IK-SC 17°/30° (Screw)	BSAE9003	10	-	10-12
ZZ BASE T-IK-SC 3.5/4.0/4.5/5.5 (Screw)	BSAE9001	15	-	10-12
ZZ BASE T-IK-ZSC 0° (Zirconia Screw)	BSAE9032	10	-	10-12
ZZ BASE T-IK-ZSC 17°/30° (Zirconia Screw)	BSAE9033	10	-	10-12
ZZ BASE T-IK-ZSC 3.5/4.0/4.5/5.5 (Zirconia Screw)	BSAE9031	15	-	10-12

THOMMEN SPI®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE T-SPI-ABCH-ID 3.5 (Raw-Abutment® customizable hex ID and screw)	BSAF3491	15	-	10-12
SET ZZ BASE T-SPI-ABCH-ID 4.0 (Raw-Abutment® customizable hex ID and screw)	BSAF3492	25	-	10-12
SET ZZ BASE T-SPI-ABCH-ID 4.5 (Raw-Abutment® customizable hex ID and screw)	BSAF3493	25	-	10-12
SET ZZ BASE T-SPI-ABCH-ID 5.0 (Raw-Abutment® customizable hex ID and screw)	BSAF3494	25	-	10-12
SET ZZ BASE T-SPI-ABCH-ID 6.0 (Raw-Abutment® customizable hex ID and screw)	BSAF3495	25	-	10-12
SET ZZ BASE T-SPI-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD3431	15	-	10-12
SET ZZ BASE T-SPI-CCTB 4.0 (Conical cemented titan base non hex and screw)	BSAD3432	25	-	10-12
SET ZZ BASE T-SPI-CCTB 4.5 (Conical cemented titan base non hex and screw)	BSAD3433	25	-	10-12
SET ZZ BASE T-SPI-CCTB 5.0 (Conical cemented titan base non hex and screw)	BSAD3434	25	-	10-12
SET ZZ BASE T-SPI-CCTB 6.0 (Conical cemented titan base non hex and screw)	BSAD3435	25	-	10-12
SET ZZ BASE T-SPI-CCTB-L20 4.0 (Conical cemented titan base non hex L20 and screw)	BSAD3432L20	25	-	10-12
SET ZZ BASE T-SPI-CCTB-L20 4.5 (Conical cemented titan base non hex L20 and screw)	BSAD3433L20	25	-	10-12
SET ZZ BASE T-SPI-CCTB-L20 5.0 (Conical cemented titan base non hex L20 and screw)	BSAD3434L20	25	-	10-12
SET ZZ BASE T-SPI-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB3481	15	-	10-12
SET ZZ BASE T-SPI-PCTBH 4.0 (Parallel cemented titan base hex and screw)	BSAB3482	25	-	10-12
SET ZZ BASE T-SPI-PCTBH 4.5 (Parallel cemented titan base hex and screw)	BSAB3483	25	-	10-12
SET ZZ BASE T-SPI-PCTBH 5.0 (Parallel cemented titan base hex and screw)	BSAB3484	25	-	10-12
SET ZZ BASE T-SPI-PCTBH 6.0 (Parallel cemented titan base hex and screw)	BSAB3485	25	-	10-12
SET ZZ BASE T-SPI-SM 3.5 (Scanmarker and screw)	ZBAD3431	-	-	10-12
SET ZZ BASE T-SPI-SM 4.0 (Scanmarker and screw)	ZBAD3432	-	-	10-12
SET ZZ BASE T-SPI-SM 4.5 (Scanmarker and screw)	ZBAD3433	-	-	10-12
SET ZZ BASE T-SPI-SM 5.0 (Scanmarker and screw)	ZBAD3434	-	-	10-12
SET ZZ BASE T-SPI-SM 6.0 (Scanmarker and screw)	ZBAD3435	-	-	10-12
SET ZZ BASE T-SPI-WSM 3.5 (White Scanmarker and screw)	ZBCD3431	-	-	10-12
SET ZZ BASE T-SPI-WSM 4.0 (White Scanmarker and screw)	ZBCD3432	-	-	10-12
SET ZZ BASE T-SPI-WSM 4.5 (White Scanmarker and screw)	ZBCD3433	-	-	10-12
SET ZZ BASE T-SPI-WSM 5.0 (White Scanmarker and screw)	ZBCD3434	-	-	10-12
SET ZZ BASE T-SPI-WSM 6.0 (White Scanmarker and screw)	ZBCD3435	-	-	10-12
ZZ BASE T-SPI-MI 3.5 (Model implant)	MIAG3401	-	-	10-12
ZZ BASE T-SPI-MI 4.0 (Model implant)	MIAG3402	-	-	10-12

THOMMEN SPI®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
ZZ BASE T-SPI-MI 4.5 (Model implant)	MIAG3403	-	-	10-12
ZZ BASE T-SPI-MI 5.0 (Model implant)	MIAG3404	-	-	10-12
ZZ BASE T-SPI-MI 6.0 (Model implant)	MIAG3405	-	-	10-12
ZZ BASE T-SPI-SC 3.5 (Screw)	BSAE3401	15	-	10-12
ZZ BASE T-SPI-SC 4.0 (Screw)	BSAE3402	25	-	10-12
ZZ BASE T-SPI-SC 4.5 (Screw)	BSAE3403	25	-	10-12
ZZ BASE T-SPI-SC 5.0 (Screw)	BSAE3404	25	-	10-12
ZZ BASE T-SPI-SC 6.0 (Screw)	BSAE3405	25	-	10-12
ZZ BASE T-SPI-ZSC 3.5 (Zirconia screw)	BSAE3431	25	-	10-12
ZZ BASE T-SPI-ZSC 4.0 (Zirconia screw)	BSAE3432	25	-	10-12
ZZ BASE T-SPI-ZSC 4.5 (Zirconia screw)	BSAE3433	25	-	10-12
ZZ BASE T-SPI-ZSC 5.0 (Zirconia screw)	BSAE3434	25	-	10-12
ZZ BASE T-SPI-ZSC 6.0 (Zirconia screw)	BSAE3435	25	-	10-12

THOMMEN VARIOMULTI

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE T-VM-CCTB 4.8 (Conical cemented titan base non hex and screw)	BSAD6131	25	-	10-12
SET ZZ BASE T-VM-SM 4.8 (Scanmarker and screw)	ZBAD6131	-	-	10-12
SET ZZ BASE T-VM-WSM 4.8 (White Scanmarker and screw)	ZBCD6131	-	-	10-12
ZZ BASE T-VM-SC 4.8 (Screw)	BSAE6101	25	-	10-12
ZZ BASE T-VM-ZSC 4.8 (Zirconia screw)	BSAE6131	25	-	10-12

VARIOBASE® FOR CROWN

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
Straumann® Basal Screw NC	ZZ-025-2900	35	-	10-12
Straumann® Basal Screw NNC	ZZ-048-313	35	-	10-12
Straumann® Basal Screw RC	ZZ-025-4900	35	-	10-12
Straumann® Basal Screw Synocta®	ZZ-048-356	35	-	10-12
Straumann® NC Variobase® Abutment for crown AH 3.5 with screw	ZZ-025-2921	35	-	10-12
Straumann® NC Variobase® for crown AH 5.5 with screw	ZZ-022-0027	35	-	10-12
Straumann® NC Variobase® for Crown incl. screw, Ø 3.9mm, AH 3.5mm, GH 2mm, TAN	ZZ-022-0102	35	-	10-12
Straumann® NC Variobase® for Crown incl. screw, Ø 3.9mm, AH 3.5mm, GH 3mm, TAN	ZZ-022-0104	35	-	10-12
Straumann® NC Variobase® for Crown incl. screw, Ø 3.9mm, AH 5.5mm, GH 2mm, TAN	ZZ-022-0106	35	-	10-12
Straumann® NC Variobase® for Crown incl. screw, Ø 3.9mm, AH 5.5mm, GH 3mm, TAN	ZZ-022-0108	35	-	10-12
Straumann® NNC Variobase® Abutment for crown AH 3.75 with screw	ZZ-048-709	35	-	10-12
Straumann® NNC Variobase® for crown AH 5.75 with screw	ZZ-022-0021	35	-	10-12
Straumann® RC Variobase® Abutment for crown AH 3.5 with screw	ZZ-025-4921	35	-	10-12
Straumann® RC Variobase® for crown AH 5.5 with screw	ZZ-022-0026	35	-	10-12
Straumann® RC Variobase® for Crown incl. screw, Ø 4.5mm, AH 3.5mm, GH 2mm, TAN	ZZ-022-0103	35	-	10-12
Straumann® RC Variobase® for Crown incl. screw, Ø 4.5mm, AH 3.5mm, GH 3mm, TAN	ZZ-022-0105	35	-	10-12
Straumann® RC Variobase® for Crown incl. screw, Ø 4.5mm, AH 5.5mm, GH 2mm, TAN	ZZ-022-0107	35	-	10-12
Straumann® RC Variobase® for Crown incl. screw, Ø 4.5mm, AH 5.5mm, GH 3mm, TAN	ZZ-022-0109	35	-	10-12
Straumann® RN Variobase® Abutment for crown AH 4.0 with screw	ZZ-048-710	35	-	10-12
Straumann® RN Variobase® for crown AH 6 with screw	ZZ-022-0022	35	-	10-12
Straumann® WN Variobase® Abutment for crown AH 4.5 with screw	ZZ-048-711	35	-	10-12
Straumann® WN Variobase® for crown AH 6.5 with screw	ZZ-022-0023	35	-	10-12

WARANTEC ONEPLANT

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE WT-OI-ABCH 3.6-5.3 (Raw-Abutment® customizable hex and screw)	BSAF7481	30	-	10-12

ZERAMEX® P6

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
ZZ BASE ZX-P6-SM RN (Scanmarker)	ZBAD10232	-	-	10-12
ZZ BASE ZX-P6-SM SN (Scanmarker)	ZBAD10231	-	-	10-12
ZZ BASE ZX-P6-WSM RN (White Scanmarker)	ZBCD10232	-	-	10-12
ZZ BASE ZX-P6-WSM SN (White Scanmarker)	ZBCD10231	-	-	10-12

ZIMMER TAPERED SCREW-VENT®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE Z-TS-ABCH 3.5 (Raw-Abutment® customizable hex and screw)	BSAF1181	30	-	10-12
SET ZZ BASE Z-TS-ABCH 4.5 (Raw-Abutment® customizable hex and screw)	BSAF1182	30	-	10-12
SET ZZ BASE Z-TS-ABCH 5.7 (Raw-Abutment® customizable hex and screw)	BSAF1183	30	-	10-12
SET ZZ BASE Z-TS-CCTB 3.5 (Conical cemented titan base non hex and screw)	BSAD1131	30	-	10-12
SET ZZ BASE Z-TS-CCTB 4.5 (Conical cemented titan base non hex and screw)	BSAD1132	30	-	10-12
SET ZZ BASE Z-TS-CCTB 5.7 (Conical cemented titan base non hex and screw)	BSAD1133	30	-	10-12
SET ZZ BASE Z-TS-CCTB-L10 3.5 (Conical cemented titan base non hex L10 and screw)	BSAD1131L10	30	-	10-12
SET ZZ BASE Z-TS-CCTB-L10 4.5 (Conical cemented titan base non hex L10 and screw)	BSAD1132L10	30	-	10-12
SET ZZ BASE Z-TS-CCTB-L10 5.7 (Conical cemented titan base non hex L10 and screw)	BSAD1133L10	30	-	10-12
SET ZZ BASE Z-TS-CCTB-L20 3.5 (Conical cemented titan base non hex L20 and screw)	BSAD1131L20	30	-	10-12
SET ZZ BASE Z-TS-CCTB-L20 4.5 (Conical cemented titan base non hex L20 and screw)	BSAD1132L20	30	-	10-12
SET ZZ BASE Z-TS-CCTB-L20 5.7 (Conical cemented titan base non hex L20 and screw)	BSAD1133L20	30	-	10-12
SET ZZ BASE Z-TS-CCTB-L30 3.5 (Conical cemented titan base non hex L30 and screw)	BSAD1131L30	30	-	10-12
SET ZZ BASE Z-TS-CCTB-L30 4.5 (Conical cemented titan base non hex L30 and screw)	BSAD1132L30	30	-	10-12
SET ZZ BASE Z-TS-CCTB-L30 5.7 (Conical cemented titan base non hex L30 and screw)	BSAD1133L30	30	-	10-12
SET ZZ BASE Z-TS-CCTB-L40 3.5 (Conical cemented titan base non hex L40 and screw)	BSAD1131L40	30	-	10-12
SET ZZ BASE Z-TS-CCTB-L40 4.5 (Conical cemented titan base non hex L40 and screw)	BSAD1132L40	30	-	10-12
SET ZZ BASE Z-TS-CCTB-L40 5.7 (Conical cemented titan base non hex L40 and screw)	BSAD1133L40	30	-	10-12
SET ZZ BASE Z-TS-NTB 3.5 (Narrow titan base non hex and screw)	BSAP1131	30	-	10-12
SET ZZ BASE Z-TS-NTBH 3.5 (Narrow titan base hex and screw)	BSAP1181	30	-	10-12
SET ZZ BASE Z-TS-PCTBH 3.5 (Parallel cemented titan base hex and screw)	BSAB1181	30	-	10-12
SET ZZ BASE Z-TS-PCTBH 4.5 (Parallel cemented titan base hex and screw)	BSAB1182	30	-	10-12
SET ZZ BASE Z-TS-PCTBH 5.7 (Parallel cemented titan base hex and screw)	BSAB1183	30	-	10-12
SET ZZ BASE Z-TS-PCTBH-L10 3.5 (Parallel cemented titan base hex L10 and screw)	BSAB1191L10	30	-	10-12
SET ZZ BASE Z-TS-PCTBH-L10 4.5 (Parallel cemented titan base hex L10 and screw)	BSAB1192L10	30	-	10-12
SET ZZ BASE Z-TS-PCTBH-L10 5.7 (Parallel cemented titan base hex L10 and screw)	BSAB1193L10	30	-	10-12
SET ZZ BASE Z-TS-PCTBH-L20 3.5 (Parallel cemented titan base hex L20 and screw)	BSAB1191L20	30	-	10-12
SET ZZ BASE Z-TS-PCTBH-L20 4.5 (Parallel cemented titan base hex L20 and screw)	BSAB1192L20	30	-	10-12
SET ZZ BASE Z-TS-PCTBH-L20 5.7 (Parallel cemented titan base hex L20 and screw)	BSAB1193L20	30	-	10-12
SET ZZ BASE Z-TS-PCTBH-L30 3.5 (Parallel cemented titan base hex L30 and screw)	BSAB1191L30	30	-	10-12

All information is subject to change. Errors and omissions excepted. Version: 27.04.2018

ZIMMER TAPERED SCREW-VENT®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE Z-TS-PCTBH-L30 4.5 (Parallel cemented titan base hex L30 and screw)	BSAB1192L30	30	-	10-12
SET ZZ BASE Z-TS-PCTBH-L30 5.7 (Parallel cemented titan base hex L30 and screw)	BSAB1193L30	30	-	10-12
SET ZZ BASE Z-TS-PCTBH-L40 3.5 (Parallel cemented titan base hex L40 and screw)	BSAB1191L40	30	-	10-12
SET ZZ BASE Z-TS-PCTBH-L40 4.5 (Parallel cemented titan base hex L40 and screw)	BSAB1192L40	30	-	10-12
SET ZZ BASE Z-TS-PCTBH-L40 5.7 (Parallel cemented titan base hex L40 and screw)	BSAB1193L40	30	-	10-12
SET ZZ BASE Z-TS-SM 3.5 (Scanmarker and screw)	ZBAD1131	-	-	10-12
SET ZZ BASE Z-TS-SM 4.5 (Scanmarker and screw)	ZBAD1132	-	-	10-12
SET ZZ BASE Z-TS-SM 5.7 (Scanmarker and screw)	ZBAD1133	-	-	10-12
SET ZZ BASE Z-TS-WSM 3.5 (White Scanmarker and screw)	ZBCD1131	-	-	10-12
SET ZZ BASE Z-TS-WSM 4.5 (White Scanmarker and screw)	ZBCD1132	-	-	10-12
SET ZZ BASE Z-TS-WSM 5.7 (White Scanmarker and screw)	ZBCD1133	-	-	10-12
SET ZZ BASE Z-TS-ZZMUA-L05 3.5 (Zirkonzahn Multi Unit Abutment non hex L05 and Screw)	BSAH1131L05	25	30	10-12
SET ZZ BASE Z-TS-ZZMUA-L15 3.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH1131L15	25	30	10-12
SET ZZ BASE Z-TS-ZZMUA-L15 4.5 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH1132L15	25	30	10-12
SET ZZ BASE Z-TS-ZZMUA-L15 5.7 (Zirkonzahn Multi Unit Abutment non hex L15 and Screw)	BSAH1133L15	25	30	10-12
SET ZZ BASE Z-TS-ZZMUA-L25 3.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH1131L25	25	30	10-12
SET ZZ BASE Z-TS-ZZMUA-L25 4.5 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH1132L25	25	30	10-12
SET ZZ BASE Z-TS-ZZMUA-L25 5.7 (Zirkonzahn Multi Unit Abutment non hex L25 and Screw)	BSAH1133L25	25	30	10-12
SET ZZ BASE Z-TS-ZZMUA-L35 3.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH1131L35	25	30	10-12
SET ZZ BASE Z-TS-ZZMUA-L35 4.5 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH1132L35	25	30	10-12
SET ZZ BASE Z-TS-ZZMUA-L35 5.7 (Zirkonzahn Multi Unit Abutment non hex L35 and Screw)	BSAH1133L35	25	30	10-12
SET ZZ BASE Z-TS-ZZMUA-L45 3.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH1131L45	25	30	10-12
SET ZZ BASE Z-TS-ZZMUA-L45 4.5 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH1132L45	25	30	10-12
SET ZZ BASE Z-TS-ZZMUA-L45 5.7 (Zirkonzahn Multi Unit Abutment non hex L45 and Screw)	BSAH1133L45	25	30	10-12
ZZ BASE Z-TS-MI 3.5 (Model implant)	MIAG1101	-	-	10-12
ZZ BASE Z-TS-MI 4.5 (Model implant)	MIAG1102	-	-	10-12
ZZ BASE Z-TS-MI 5.7 (Model implant)	MIAG1103	-	-	10-12
ZZ BASE Z-TS-SC 3.5 (Screw)	BSAE1101	30	-	10-12
ZZ BASE Z-TS-SC 4.5 (Screw)	BSAE1102	30	-	10-12
ZZ BASE Z-TS-SC 5.7 (Screw)	BSAE1103	30	-	10-12

ZIMMER TAPERED SCREW-VENT®

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
ZZ BASE Z-TS-SC-G 3.5 (Screw Gold)	BSCE1101	30	-	10-12
ZZ BASE Z-TS-SC-G 4.5 (Screw Gold)	BSCE1102	30	-	10-12
ZZ BASE Z-TS-SC-G 5.7 (Screw Gold)	BSCE1103	30	-	10-12
ZZ BASE Z-TS-ZSC 3.5 (Zirconia screw)	BSAE1131	30	-	10-12
ZZ BASE Z-TS-ZSC 4.5 (Zirconia screw)	BSAE1132	30	-	10-12
ZZ BASE Z-TS-ZSC 5.7 (Zirconia screw)	BSAE1133	30	-	10-12
ZZ BASE Z-TS-ZSC-G 3.5 (Zirconia Screw Gold)	BSCE1131	20	-	10-12
ZZ BASE Z-TS-ZSC-G 4.5 (Zirconia Screw Gold)	BSCE1132	20	-	10-12
ZZ BASE Z-TS-ZSC-G 5.7 (Zirconia Screw Gold)	BSCE1133	20	-	10-12

ZIMMER TAPERED SCREW-VENT® MULTI UNIT ABUTI

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE Z-TS-MUA-CCTB (Conical cemented titan base non hex and screw)	BSAD6031	30	-	10-12
SET ZZ BASE Z-TS-MUA-SM (Scanmarker and screw)	ZBAD6031	-	-	10-12
SET ZZ BASE Z-TS-MUA-WSM (White Scanmarker and screw)	ZBCD6031	-	-	10-12
ZZ BASE Z-TS-MUA-SC (Screw)	BSAE6001	30	-	10-12
ZZ BASE Z-TS-MUA-ZSC (Zirconia screw)	BSAE6031	30	-	10-12

ZIRKONZAHN MULTI UNIT ABUTMENT

Item description	Item number	Recommended torque in Ncm SC	Recommended torque in Ncm MUA	Composite resin temporaries in Ncm
SET ZZ BASE ZZ-MUA-CCTB (Conical cemented titan base non hex and screw)	BSAD9231	30	-	10-12
SET ZZ BASE ZZ-MUA-CCTB-L10 (Conical cemented titan base non hex L10 and screw)	BSAD9231L10	30	-	10-12
SET ZZ BASE ZZ-MUA-IC (Impression coping and impression coping screw)	BSAQ9231	10	-	10-12
SET ZZ BASE ZZ-MUA-IC-OP (Impression coping and impression coping screw)	BSAQ9233	10	-	10-12
SET ZZ BASE ZZ-MUA-IC-SP (Impression coping and impression coping screw)	BSAQ9232	10	-	10-12
SET ZZ BASE ZZ-MUA-NTBH-OP (Narrow titan base hex and screw)	BSAP9283	30	-	10-12
SET ZZ BASE ZZ-MUA-NTBH-SP (Narrow titan base hex and screw)	BSAP9282	30	-	10-12
SET ZZ BASE ZZ-MUA-SA (Scanalog and screw)	ZBDE9231	15	-	10-12
SET ZZ BASE ZZ-MUA-SM (Scanmarker and screw)	ZBAD9231	-	-	10-12
SET ZZ BASE ZZ-MUA-SM OP (Scanmarker and screw)	ZBAD9233	-	-	10-12
SET ZZ BASE ZZ-MUA-SM SP (Scanmarker and screw)	ZBAD9232	-	-	10-12
SET ZZ BASE ZZ-MUA-WSM (White scanmarker and screw)	ZBCD9231	-	-	10-12
SET ZZ BASE ZZ-MUA-WSM OP (White Scanmarker and screw)	ZBCD9233	-	-	10-12
SET ZZ BASE ZZ-MUA-WSM SP (White Scanmarker and screw)	ZBCD9232	-	-	10-12
ZZ BASE ZZ-MUA-MI (Model implant)	MIAG9201	-	-	10-12
ZZ BASE ZZ-MUA-SC (Screw)	BSAE9201	30	-	10-12
ZZ BASE ZZ-MUA-ZSC (Zirconia Screw)	BSAE9231	30	-	10-12