



Penn

CAD/CAM-Ceramic-Center

Director: Michael Bergler, MDT
School of Dental Medicine,
Robert Schattner Building
240 S. 40th Street,
Philadelphia, PA 19104
T. 215-746-2075 F. 215-898-9981
e-mail. mbergler@dental.upenn.edu

Flexural Strength Test Result of Aged Anterior Zirconia, Zirkozahn GmbH, Gais, Italy

Material	Flexural Strength (MPa)	
	As received	Thermocycled (4 year Aging)
Anterior T 0	670.0	666.2

Dr. Francis K. Mante
University of Pennsylvania
School of Dental Medicine

Michael Bergler, CDT, MDT
University of Pennsylvania
School of Dental Medicine