



100% PRETTAU[®] ZIRCONIA - MADE BY CAD/CAM

"Perfection through perseverance" Enrico Steger



Zirkonzahn World Wide



Zirkonzahn Germany



Zirkonzahn USA



Zirkonzahn Mexiko


Zirkon zahn

Prettau[®] Zirconia

In order to demonstrate the functionality of these types of construction, we collaborate with Dr. Fernando Rojas-Vizcaya from Spain and manufacture several works each month. The advantages of occlusally screwed Prettau[®] bridges are their easy fabrication and their stability.

Using wax-up and double-scan, the construction is realized in a short period of time. The system creates the occlual screw accesses automatically. The buccal part of the front teeth has been veneered, resulting in a very high esthetic in the visible part and high stability in the lateral tooth area.







The Prettau[®] Bridge

Dr. Fernando Rojas-Vizcaya, Spain Georg Walcher – Zirkonzahn Education Center







Zirkon zahn[®]



HUMAN ZIRCONIUM TECHNOLOGY

ZIRKONZAHN World Wide - Tel +39 0474 066 680 - info@zirkonzahn.com - www.zirkonzahn.com

	CAD CAM MIS	SCANNER BSCO		
. 6		a.	÷	4
		¢.	2	-
	Zirkonzahn	5	Zirkonzahn	
		Milled with	Zirkonzahn CAD/CAM System	