


Zirkonzahn®



AESTHETICS REQUIRES SKILL AND PERFECTION


In addition to a drive for perfection and a high level of dedication, a good dental technician needs to have, above all, a precise aesthetic goal in mind.



GINGIVA-COMPOSITES

The Gingiva-Composites veneering composites have been used to veneer the gingival aspect of a denture base made of Denture Gingiva Basic, a PMMA resin. Depending on the technician's liking, either a viscous or paste-like composite can be used for this purpose. To achieve particularly high aesthetics, the structure can further be high-gloss polished with Composite Polishing Paste.

DT Alexander Lichtmannegger – Zirkonzahn Education Center Brunico, South Tyrol, Italy



Zirkonzahn®



MORE
INFORMATION
ON THE
GINGIVA-
COMPOSITES



GINGIVA-COMPOSITES

- For veneering resin structures
- Enable the individual design of the gingival areas of resin temporaries as well as final restorations
- Colour spectrum based on the ICE Ceramics Tissue colours (light to dark). Dentists and patients get an immediate aesthetic impression of the final restoration
- Especially suitable for long-term temporaries

HUMAN ZIRCONIUM TECHNOLOGY

Zirkonzahn Worldwide – Tel +39 0474 066 680 – info@zirkonzahn.com – www.zirkonzahn.com



Ø125
mm



Ø95
mm



TELESKOPER
ORBIT WITH
Ø 125 mm
IN THE M2
MILLING UNIT
COMFORT LINE



WEAB3545EN=

This is an international advertisement. Not all of the products mentioned are available in all countries. The product's field of application may vary depending on the country. Please contact your sales team for more information. All information is subject to change. Errors and omissions excepted. Version: 28/04/2022