THE SOUTH TYROLEAN DOLOMITES SOURCE OF INSPIRATION & STRENGTH

As unique as this natural monolith, is zirconia the dental material of choice. To imbue it with unique properties. To improve it. To make it resistant beyond our lifetime. To infuse luminosity and all its extraordinary attributes, gives us, as producers of this noble material, inspired by our South Tyrolean homeland and scientific spirit, the unbelievably immense feeling of doing the right thing.

PRETTAU®

ZIRCONIA FROM THE DOLOMITES

THE MOST EXPENSIVE - FOR THE PATIENTS YOU VALUE

Zirkonzahn

PRETTAU[®] ZIRCONIA MANUFACTURED WITH PRIDE

Without zirconia there would be no Zirkonzahn. Strongly convinced of its base properties, I decided to make the best of it for us dental technicians. To achieve that, much more was necessary than initially thought – premium raw materials, expensive manufacturing and cleaning processes, colouring concepts, processing tools, milling strategies both manual and digital, sintering furnaces and consistently new solutions. But we never gave up.


We studied everything to the last detail and assembled the necessary puzzle that is Prettau[®] zirconia.

Prettau[®] – a name deeply connected to my South Tyrolean homeland. A commitment and a promise at the same time. Now, we are pushing the path of our Prettau[®] zirconia toward monolithic design.

Zirconia is and remains our passions – and our pride.

You can count on that!

finico Stejen



****Gradual-Triplex-Technology:** The zirconia is already provided with colour, tanslucency and bending strength gradients during the production process. This ensures a natural colouring, a highly translucent cutting edge and an extremely high bending strength at the tooth neck.

*Average value of the biaxial bending strength from several test series



PRETTAU® 4 ANTTURN'S DIS.

MORE INFORMATION

WORLD NOVELTY

ZIRCONIA SHADE GUIDES FOR UPPER AND LOWER **ANTERIORS, POSTERIORS AND WITH MINIMAL CUTBACK**

NEW! ZIRKONZAHN SHADE GUIDES

PRETTAU® 3

High-quality monolithic zirconia shade guides for Prettau[®] 2 Dispersive[®], Prettau[®] 3 Dispersive[®] and Prettau[®] 4 Anterior[®] Dispersive[®]

All information is subject to change. Errors and omissions excepted. Version: 19/01/2021

THE MOST EXPERSIVE - DESCENSE PROM THE PRODUCTS



Zirkonzahn Worldwide – Tel +39 0474 066 680 – info@zirkonzahn.com – www.zirkonzahn.com

LC ADJS AN RT

tion .

113

140

(C)

6

itter -At

DISPUSA

PRETTAU® 2

PRETTAU® 2