



**Zirkonzahn®**

*Human Zirconium Technology*



## **BITE SPLINTS**

*Software-module CAD/CAM*











## SOFTWARE-MODULE CAD/CAM BITE SPLINTS

*Module for the manufacture of bite splints*

- *precise modelation of the necessary contact points in conjunction with the virtual articulator*
- *individual shaping*
- *adjustable parameters: occlusal thickness, wall thickness, smoothness (blocking-out)*
- *millable in various plastic materials*

*Item number: SY0530*





**Zirkonzahn®**



## HUMAN ZIRCONIUM TECHNOLOGY

*Zirkonzahn Worldwide – Tel +39 0474 066 680 – info@zirkonzahn.com – www.zirkonzahn.com*

*Milled with the M1 Wet Heavy Metal CAD/CAM system*



WEAA5493=

*All information is subject to change. Errors and omissions excepted. Version: 20/03/2018*