

# Zirkonzahn®


*Human Zirconium Technology*




**NEW**



## EGGSHELL TEMPORARY

*Software-module CAD/CAM*










## SOFTWARE-MODULE CAD/CAM EGGSHELL TEMPORARY

*Module for the realization of individual eggshell temporaries as esthetically appealing immediate restoration for dentist and patient.*

- *Module for the creation of individual eggshell temporaries*
- *Import of intraoral scan data possible*
- *Adjustable parametres: depth and type of preparation, wall thickness*
- *Can be modeled very thinly (0,3 mm) and is easily adaptable*
- *Esthetically appealing and individually modeled immediate restoration for dentist and patient*

*Item number: SY0540*





## HUMAN ZIRCONIUM TECHNOLOGY

ZIRKONZAHN Worldwide - Tel +39 0474 066 680 - [info@zirkonzahn.com](mailto:info@zirkonzahn.com) - [www.zirkonzahn.com](http://www.zirkonzahn.com)



WEAA5494=