


# Zirkonzahn<sup>®</sup>

*Human Zirconium Technology*



**ALL ZIRCONIA**

## *Combined technology made of all-zirconia*

Upper jaw – Prettau bridge on 4 ZrO<sub>2</sub>- implants and locators (Camlog)

Lower jaw – Prettau bridge on 4 hybrid telescopes (2 thereof on implants) and activatable telescope attachments (TK1)

Both ensure good wear comfort and excellent cleaning properties for the patient.


The layering of ceramic is limited to the gingiva which was also milled out of one piece of ZrO<sub>2</sub>, together with the dental structure. The photographs show the basal view of the supports for the telescopes or locator matrices and the homogeneous all-zirconia occlusal support on the gingiva.



## *Final result of the removable Prettau bridges*

Master Technician Alexander Kanotscher, Austria

See [www.zahnprofi.at](http://www.zahnprofi.at) for the complete work sequence.



# Zirkonzahn®



## HUMAN ZIRCONIUM TECHNOLOGY

ZIRKONZAHN USA Inc. - Tel 1 800 989 8931 - Tel 678 441 9419 - [infousa@zirkonzahn.com](mailto:infousa@zirkonzahn.com)

ZIRKONZAHN World Wide - Tel +39 0474 066 680 - [info@zirkonzahn.com](mailto:info@zirkonzahn.com) - [www.zirkonzahn.com](http://www.zirkonzahn.com)



*Milled with the Zirkonzahn Milling System*