

Zirkonzahn[®]

Human Zirconium Technology



ALL ZIRCONIA

Prettau Zirconia

This Prettau Bridge was 2 toned coloured in the pre-sintered stage by using the special Prettau Colour Liquids applied with a metal-free brush.

The premolar abutments diverge severely calling for a primary substructure. The secondary part is fixed to the primary with a screw. The entire construction is made from 100% zirconia. Only the pink tissue part was completed with ICE Ceramic veneering porcelain.



Prettau Zirconia

Georg Walcher, South Tyrol and Dr. Fernando Rojas-Vizcaya, Spain



Zirkonzahn®



HUMAN ZIRCONIUM TECHNOLOGY

ZIRKONZAHN USA Inc. - Tel 1 800 989 8931 - Tel 678 441 9419 - infousa@zirkonzahn.com

ZIRKONZAHN World Wide - Tel +39 0474 066 680 - info@zirkonzahn.com - www.zirkonzahn.com



WEAA5408=

Milled with Zirkonzahn Milling System