

Zirkonzahn®

Human Zirconium Technology



PRETTAU® ANTERIOR®



INCREASE YOUR PRICE WITH PRETTAU® ANTERIOR®

Transparency combined with high flexural strength! The new zirconia material has been specifically developed for anterior tooth restorations, but is just as well suited for posterior areas. Prettau® Anterior® can be milled dry with regular zirconia burs and therefore offers an ideal alternative to lithium disilicate!

- Double the flexural strength compared to lithium disilicate: 670 MPa (School of Dental Medicine, University of Pennsylvania)
- Excellent esthetic properties
- Higher transparency allows perfect light transmission and good background masking
- Extra resistant to ageing
- Suitable for fully anatomical contours or framework for ceramic veneering
- Shorter sinter times due to sintering at 1,500 °C (2,732 °F)
- Available in various block sizes
- Indications: metal-free partial and individual crowns, inlays, onlays, veneers, 3-unit bridges (max).





Zirkonzahn®



Perfect esthetics also in the posterior area with Prettau® Anterior®, colored with ICE Zirkon Malfarben 3D by Enrico Steger surface stains.

HUMAN ZIRCONIUM TECHNOLOGY

ZIRKONZAHN Worldwide – Tel +39 0474 066 680 – info@zirkonzahn.com – www.zirkonzahn.com

ZIRKONZAHN USA Inc. – Phone +1 800 989 8931 – Phone +1 678 441 9419 – infousa@zirkonzahn.com



Prettau® Anterior®



Prettau® Zirconia